



## **AIR FRAME WHEELCHAIR**

**Owner's Manual**  
**MOB1062**



To see all FAQs in one place  
visit **[vhealth.link/3ew](https://vhealth.link/3ew)**

# INSTRUCTIONS



For a video demonstration  
check out **[vhealth.link/zpf](https://vhealth.link/zpf)**

1. Unfold the chair by pushing outward against the armrests. (to refold the chair it is easiest to do by pulling up on the seat of the chair)



2. Pull up on the backrest handles and make sure they click lock into place on both sides so the backrest is upright. (To refold this back down there are black release tabs on the back side of the folding joint.)



3. \*If Necessary\* You can add the included foot securing strap. This can be used as needed along with the seatbelt. The seatbelt is strongly recommended while the chair is in use. You can also adjust the footrest with the included tools (also only if necessary).



## USING THE BRAKES

If braking while in the chair a hand brake is provided for each wheel. Rear brakes are also provided for use if being assisted by another person. Note that the handbrakes can lock the wheels til no longer needed but the rear brakes cannot lock they can only brake while the brake handles are being squeezed.



### **WARNING** ⚠

It is best to always use the brakes while not moving. **ALWAYS** lock the brakes when getting out of the chair or while stopping on any type of incline or decline.

### **CAUTION** ⚠

Please note that this chair is designed to use for transporting another individual and is not meant for self transportation. This Transport Chair is designed to safely be a stable center of gravity but it is important to try and keep a stable center of gravity of all times while performing activities that may change your center of gravity or a passenger's center of gravity. Never reach out so far that it would require you to shift in the seat or bend over out of the chair. When you have to reach for objects that require you to lean, be certain your front wheels are pointing forward as this will stabilize the chair.

### **RAMPS**

Be sure of your own capabilities, limitations and strength before attempting to push another individual up a ramp. Correct techniques on ramp negotiation should be shown to you by a qualified professional. If you must stop on an incline, avoid any and all sudden or abrupt movement when you resume movement, because of backward angle of the chair, sudden abrupt movement could cause instability. Do not attempt to change direction when moving in incline or decline, this could result in tipping.

### **TIRE MAINTENANCE**

Tires should be cleaned occasionally with a damp cloth, Replace if tires become severely worn or cracking appears. Front casters should be checked periodically for cracking and lubricated occasionally rear wheels should be inspected at least once a month to ensure proper tightness.

### **GENERAL MAINTENANCE**

Wipe metal parts with a soft cloth at least once a week. All aluminum parts can be cleaned with auto wax or similar products. If exposed to moisture, dry immediately. Do not use an abrasive cleanser that will scratch the finish.

## WARNINGS

- Wheelchair will support up to a 240lbs individual. Do not exceed chair weight limit.
- Do not attempt to use wheelchair without reading carefully all instructions in this manual.
- This wheelchair has not been approved as a seat that can be used within vehicles of any kind.
- Always transfer wheelchair user to a vehicle manufacturers approved seat and use safety restraints available.
- Do not operate this wheelchair on roadways, streets or any surfaces with vehicular traffic.
- Do not attempt to operate chair on inclines without proper assistance.
- Wheelchair should be on a smooth, stable and level area with the hand brakes engaged by the transporter before transferring the passenger in or out of the wheelchair.
- When transferring in or out of wheelchair, do not stand on or apply weight to foot rests.
- Do not lean over or apply weight to the top of wheelchair back. This can cause wheelchair to tip over.
- Do not stand on wheelchair.
- Abruptly changing directions while going down an incline could cause instability and/or increase the possibility of tipping over.

## HAVE MORE QUESTIONS?

Check out our list of Frequently Asked Questions at **vhealth.link/3ew** for helpful answers.



And if that doesn't answer your question, our customer service team would love to help! Feel free to connect with them by phone, e-mail, or chat on our website.



**sales@aracent.com**



**1-888-550-7750**



**aracent.com**

Distributed by

**vive  
health**

8955 Fontana Del Sol Way  
Naples, FL 34109  
Made in China