

INVACARE®
9000 SL
WHEELCHAIR



The **Invacare® 9000 SL Wheelchair** is the perfect high-strength lightweight rental product for those who want quality and durability at a moderate price. This low-maintenance wheelchair also has a full complement of options and accessories. Its contemporary design and high-performance frame make the 9000 SL wheelchair a superior value.

Features

Stylish, contemporary carbon steel frame is long-lasting
 Multiple frame color variations for customizing
 Dual axle positions allow variations of seat-to-floor heights
 Urethane rear tires mounted on "no flex" wheels offer superior performance
 Urethane casters for better ride and lighter weight
 Low-maintenance, dual-sealed precision bearings
 Heavy-duty inner liners keep seat and back from stretching
 Space-saver arms reduce overall frame width without limiting seat width
 Meets RESNA standards for manual wheelchairs in product performance

HCPCS Code: K0004



1. Dual axle positions
2. Urethane casters
3. Space saver arms
4. Seating options available
5. Builder options available



Aracent Healthcare
www.aracent.com

E: sales@aracent.com
 P: 888-550-7750
 F: 888-550-7752

©2014
 Invacare Corporation.
 All rights reserved.
 Trademarks are identified by the symbols™ and®.
 All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.
 Form No. 11-226
 Rev. 04/14
 131067

Model #	Overall Width	Overall Height	Overall Depth w/o riggings	Folded Width	Seat Dimensions	Seat-to-Floor Height	Back Height
9SL	23", 25"	36"	30"	15"	16"-20"x16" or 18"	15.5"-19.5"	15"-17", 17"-19"

Caster Height	Wheel Height	Weight w/o Rigging	Weight Capacity	Shipping Weight	Carton Dimensions	Limited Warranty
6", 8"	20", 22", 24"	33 lb.	250 lb.	43 lb.	32"x13"x 37"	5 Years

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

