

INVACARE®
9000 TOPAZ™
WHEELCHAIR



The **Invacare® 9000 Topaz™ Wheelchair** is designed specifically to suit the unique needs of bariatric patients. The Topaz wheelchair combines the strength and durability of a carbon steel frame with the comfort of heavy-duty front riggings and arms for proper positioning. With a variety of seat widths and depths, frame reinforcements that include gusseted sides, caster journals and dual heavy-duty crossbraces, this chair provides not only ease in maneuverability but proven stability.

Features

Available in vinyl or nylon upholstery

Durable low-maintenance, carbon steel frame is long-lasting

Reinforced frame includes gusseted sides, caster journals and dual heavy-duty crossbraces

Multiple frame color variations available for customizing

Heavy-duty front riggings provide proper positioning for feet with footrest and elevating legrest

Heavy-duty arms designed to accommodate weight capacity and provide proper positioning

Dual axle positions allow for quick and simple seat-to-floor height variations, 17.5" and 19.5"

Heavy-duty package option for 1,000 lb. weight capacity

Heavy-duty triple inner liners keep seat and back from stretching

HCPCS Code: K0007



1. Dual axle positions
2. Heavy-duty front riggings
3. Heavy-duty arms
4. Bariatric sizes
5. Seating available
6. Builder options



Aracent Healthcare
 www.aracent.com
 E: sales@aracent.com
 P: 888-550-7750
 F: 888-550-7752

Invacare Corporation
 www.invacare.com

USA
 One Invacare Way
 Elyria, Ohio
 44035-4190
 (800) 333-6900

Canada
 570 Matheson Blvd. E.,
 Unit 8
 Mississauga, Ontario
 L4Z 4G4 Canada
 (800) 668-5324

©2014
 Invacare Corporation.
 All rights reserved.
 Trademarks are
 identified by the
 symbols ™ and ®.
 All trademarks are
 owned by or licensed
 to Invacare Corporation
 unless otherwise
 noted. Specifications
 are subject to change
 without notification.
 Form No. 11-230
 Rev. 04/14
 131072



Model #	Overall Width	Overall Height	Overall Depth w/o riggings	Folded Width	Seat Dimensions	Seat-to-Floor Height	Back Height
9TPZ	39"	32"	30"	16"	20"-30" (2" inc) x 18", 20"	15.5"-19.5"	16"-18"
Caster Height	Wheel Height	Weight w/o Rigging	Weight Capacity	Shipping Weight	Carton Dimensions	Limited Warranty	
9TPZ	6", 8"	20", 24"	55 lb.	700 lb., 1,000 lb	97 lb.	35"x15"x 46.5"	3 Years

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.