

## Features



Padded armrest



Modern aesthetics with backrest detail



Easy to attach and remove legrests



Elevating Legrest Option

# Essentials™ SX3

## Technical data

For more information about this product, including the user manual, please visit our website: [www.invacareamerica.com](http://www.invacareamerica.com)

	Seat width	Overall width	Folded width	Seat depth	Overall depth (without riggings)	Overall height
Essentials™ SX3	20 inches (50.8 cm) 18 inches (45.7 cm) 16 inches (40.6 cm)	Seat width + 7.1 inches (18 cm)	11.8 inches (30 cm)	16 inches (40.6 cm) plus a 2-inch (5.1 cm) extended	31 inches (78.7 cm)	35.6 inches (90.5 cm)
	Seat to Floor Height	Back height	Rear wheel size	Weight	Shipping weight	Weight limit
Essentials™ SX3	19.5 inches (49.5 cm)	17 inches (43 cm)	24 inches (61 cm)	Chair only: 35 lb (16 kg) Chair with front rigging: 39 lb (17.7 kg)	43 lb (19.5 kg)	300 lb (136 kg)

Item No.	Description
ISX366DASA	Essential SX3 16x16 Desk Arms Swing-Away Footrests
ISX386DASA	Essential SX3 18x16 Desk Arms Swing-Away Footrests
ISX306DASA	Essential SX3 20x16 Desk Arms Swing-Away Footrests
ISX366DAEL	Essential SX3 16x16 Desk Arms Swing-Away Legrests
ISX386DAEL	Essential SX3 18x16 Desk Arms Swing-Away Legrests
ISX306DAEL	Essential SX3 20x16 Desk Arms Swing-Away Legrests

Silver vein  
Frame color



Invacare America  
One Invacare Way  
Elyria, Ohio 44035-4190  
(800) 333-6900

©2025 Invacare America. All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare® America or its subsidiaries unless otherwise noted. Specifications are subject to change without notification. Form No. 24-037rev, 250121



Lightweight and compact

for effortless everyday mobility

### ONLINE ORDERING:

Aracent Healthcare  
[www.aracent.com](http://www.aracent.com)  
E: [sales@aracent.com](mailto:sales@aracent.com)  
P: 888-550-7750  
F: 888-550-7752



# Introducing the new Essentials™ SX3

A reliable and robust manual wheelchair that ensures optimal comfort wherever you go.

Designed for ease of use and portability, the Essentials SX3 is ideal for everyday use. Swing-away footrests and flip-back desk arms facilitate easy transfers and flexibility for daily activities.

The durable build ensures you feel confident with a reliable drive every time. Customizable adjustments allow the Essentials SX3 to give a more personalized fit for optimal comfort.



## Personalized Fit

Adjustable in seat depth by two inches and extendable footrests help to achieve a more personalized fit and comfortable seating position.



## Maximum Comfort

Benefiting from a subtly padded backrest and seat, the Essentials SX3 allows you to immerse into the chair so you feel comfortable and well supported. Also easy to maintain in a wipeable vinyl finish.



## Foldable & Compact

The lightweight frame and cross folding design helps to ensure folding is simple and effortless. Compact when folded also makes storage and transportation easy.



**Black vinyl seat & back upholstery** for easy cleaning

**Flip-back armrests** for easy transfer

**Desk length armpad** gives easier access up to and under tables

**Seat depth** adjusts by 2 inches for a more precise fit

**Swing-away or elevating legrests**



\*Product specifications may vary from those displayed.