


INVACARE® MATRX®

ELAN

HEADREST



The **Invacare® Matrx® Elan® Headrest** — Takes head support to the next level. A new elegance in head support and comfort, the Elan Headrest is compact, lightweight and highly adjustable, and compliments a full range of pediatric through adult applications.

Call to order

1.800.333.6900
invacare.com

1.888.433.6818
motionconcepts.com



matrx®
SEATING SERIES

Elan Headrest Hardware Features and Functional Benefits


Elegant – Clean, compact, streamlined aesthetics.

Adjustable – Six points of adjustment – infinite positioning opportunities to meet user needs.

Easy – Single tool set up. Optional tool-less adjustment levers provided for ease of position changes.

Reliable – Multi-disk locking provides superior holding strength.

Sizes – Standard adult (15" post) or mini (11" post) available.



Elan Headrest Pads Features and Functional Benefits

Comfortable – Dual density foam provides immersion, prevents bottoming out.

Durable – Co-molded, dual firmness, HR foam surrounds backing plates, can be easily removed for modification, and are made with Ultra-fresh antimicrobial additive.

Lightweight – 45% Lighter than competitive pads.

Versatile – Standard, occipital and multi-adjustable 4-point support models available.
– Aluminum plates can be formed to desired shape.

Covers – Zippered, moisture-resistant, breathable Startex fabric
– Standard or Infection control cover (reverse Startex) available.



Model Description	Product	Dimensions Width X Height	Startex Cover Model #	Reversed Startex Cover (Infection Control) Model #	Limited Warranty	HCPCS Code
ESP-6	6" Standard Pad	6" x 3.5"	ESP6	ESP6-IC	2 Years	E0955
ESP-10	10" Standard Pad	10" x 5"	ESP10	ESP10-IC	2 Years	E0955
ESP-14	14" Standard Pad	14" x 5"	ESP14	ESP14-IC	2 Years	E0955
EOP-9	9" Occipital Pad	9" x 7"	EOP9	EOP9-IC	2 Years	E0955
EOP-12	12" Occipital Pad	12" x 8"	EOP12	EOP12-IC	2 Years	E0955
E4POINT	4 Point Pad	11" x 10"	E4POINT	E4POINT-IC	2 Years	E0955
MEHW	Elan Standard Hardware	15" Vertical Post	N/A	N/A	2 Years	E1028
MEMH	Elan Mini Hardware	11" Vertical Post	N/A	N/A	2 Years	E1028

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program and should consult an attorney or other advisor to discuss specific situations in further detail.

Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd. E.
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Motion Concepts
www.motionconcepts.com
info@motionconcepts.com

USA
700 Ensminger Rd., Suite 112
Tonawanda, NY 14150
(888) 433-6818
Fax: (888) 433-6834

Canada
84 Citation Dr., Unit 1-7
Concord, ON L4K 3C1
(866) 748-7943
(905)-695-0134
Fax: 905-695-0138

© 2016
Invacare Corporation.
All rights reserved.
Trademarks are identified by the symbols™ and ©. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Matrix is a registered trademark of Motion Concepts. Specifications are subject to change without notification.
Form No. 16-083
Rev. 09/16
160956

