

INVACARE®
TRACER® EX2
WHEELCHAIR



The **Invacare Tracer EX2 Wheelchair** sets the benchmark in manual wheelchairs. A true dual axle position allows for repositioning the 24" rear wheels and 8" casters for adult and hemi seat-to-floor heights. Offset permanent arms offer full seat width up to 20" with a permanent arm, hemi-spaced footrest and full range of accessories. Incorporating the design and technology of the Invacare 9000 series wheelchair, the Tracer EX2 wheelchair has interchangeable components with the 9000 series chairs.

Features

- Durable, low-maintenance, carbon steel frame is long-lasting
- Dual axle positions allow variation of seat-to-floor heights
- Offset permanent arms available in 16", 18" and 20" widths
- 14-gauge crossbraces add strength and durability
- Urethane rear tires mounted on "no flex" wheels offer superior performance
- Low-maintenance, dual-sealed precision bearings
- Heavy-duty inner liners keep seat and back from stretching
- Upholstery is dual-embossed

HCPCS Code: K0001/K0002



1. Dual axle positions
2. Dual-embossed upholstery
3. Offset permanent arms
4. Seating options available



Aracent Healthcare
 www.aracent.com
 E: sales@aracent.com
 P: 888-550-7750
 F: 888-550-7752

Invacare Corporation
 www.invacare.com

USA
 One Invacare Way
 Elyria, Ohio
 44035-4190
 (800) 333-6900

Canada
 570 Matheson Blvd. E.,
 Unit 8
 Mississauga, Ontario
 L4Z 4G4 Canada
 (800) 668-5324

©2017
 Invacare Corporation.
 All rights reserved.
 Trademarks are
 identified by the
 symbols™ and®.
 All trademarks are
 owned by or licensed
 to Invacare Corporation
 unless otherwise
 noted. Specifications
 are subject to change
 without notification.
 Form No. 11-224
 Rev. 06/17
 170613



Model #	Overall Width	Overall Height	Overall Depth w/o riggings	Folded Width	Seat Dimensions	Seat-to-Floor Height	Back Height
TREX26	23.25"-28.25"	34"	31.88"	13"	16"x16"	17.5", 19.5"	16"
TREX28	23.25"-28.25"	34"	31.88"	13"	18"x16"	17.5", 19.5"	16"
TREX20	23.25"-28.25"	34"	31.88"	13"	20"x16"	17.5", 19.5"	16"



	Caster Height	Wheel Height	Weight w/o Rigging	Weight Capacity	Shipping Weight	Carton Dimensions	Limited Warranty
TREX26	8"	24"	36 lb.	250 lb.	43 lb.	32"x13"x 38"	5 Years
TREX28	8"	24"	36 lb.	250 lb.	43 lb.	32"x13"x 38"	5 Years
TREX20	8"	24"	36 lb.	250 lb.	43 lb.	32"x13"x 38"	5 Years

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

