

INVACARE®
TRACER® IV
WHEELCHAIR



The **Invacare® Tracer® IV Wheelchair** features superior durability, rollability and streamlined looks. It is designed for the active individual who needs the strength of a reinforced heavy-duty frame. Its many options include components interchangeable within the Tracer series. The Tracer IV wheelchair is long-lasting and low-maintenance.

Features

- Dual axle position allows variation of seat-to-floor heights
- Durable, low-maintenance, triple chrome-plated carbon steel frame is long-lasting
- Reinforced frame includes gusseted sides, caster journals and dual crossbraces
- Corded urethane rear tires (350 lb.) and pneumatic flat free tires (450 lb.) mounted on "no flex" wheels offer superior performance
- Low-maintenance, dual-sealed precision bearings
- Heavy-duty inner liners keep seat and back from stretching
- Upholstery is dual-embossed

HCPCS Code: K0007



- Dual axle positions
- Corded urethane tires
- Dual-embossed upholstery
- Bariatric sizes available
- Seating available
- Builder options



Aracent Healthcare
 www.aracent.com
 E: sales@aracent.com
 P: 888-550-7750
 F: 888-550-7752

Invacare Corporation
 www.invacare.com

USA
 One Invacare Way
 Elyria, Ohio
 44035-4190
 (800) 333-6900

Canada
 570 Matheson Blvd. E.,
 Unit 8
 Mississauga, Ontario
 L4Z 4G4 Canada
 (800) 668-5324

©2014
 Invacare Corporation.
 All rights reserved.
 Trademarks are
 identified by the
 symbols™ and®.
 All trademarks are
 owned by or licensed
 to Invacare Corporation
 unless otherwise
 noted. Specifications
 are subject to change
 without notification.
 Form No. 11-228
 Rev. 04/14
 131070



| Model # | Overall Width | Overall Height | Overall Depth w/o riggings | Folded Width | Seat Dimensions | Seat-to-Floor Height | Back Height |
|--------------------------------|---------------|----------------|----------------------------|--------------|------------------|----------------------|-------------|
| T420, T422, T424, T4X22, T4X24 | 29" | 36" | 30" | 13" | 18"-24"x18", 20" | 17.5", 19.5" | 16" |



| | Caster Height | Wheel Height | Weight w/o Rigging | Weight Capacity | Shipping Weight | Carton Dimensions | Limited Warranty |
|--------------------------------|---------------|--------------|--------------------|------------------|-----------------|-------------------|------------------|
| T420, T422, T424, T4X22, T4X24 | 8" | 24" | 42 lb. | 350 lb., 450 lb. | 53 lb. | 32"x13"x 37" | 5 Years |

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.