

Astra Aluminum Transport Chair

The **Astra Aluminum Transport Chair** is a compact, extra lightweight design with an easy-to-maintain finish. The back-release hinge allows the chair to fold down flat and compact for easy transport and storage. The Astra comes complete with seat belt, tool-free swing-away composite footrests and heel loops.



Seat belt



Back-release hinge



Folds down flat



Push-to-lock rear wheel locks

FEATURES

- Lightweight aluminum construction
- Flared handles and back-release hinge allow the back to fold down flat and compact for easy storage and transport
- Stylish aluminum frame provides extra strength and durability
- Padded full-length armrests
- Padded nylon upholstery is comfortable, lightweight and easy to clean
- Seat belt provides added safety
- Nylon seat straps allow for ease of folding and lifting
- Push-to-lock rear wheel locks
- Composite 8" wheels with solid rubber tires are durable and maintenance free
- Tool-free, swing-away composite footrests with heel loops included
- HCPCS Code: E1038

COLORS

Black Finish



Blue Finish



SPECIFICATIONS

Seat Width:	19"
Seat Depth:	16"
Seat-to-Floor Height:	19"
Back Height:	19"
Back Style:	Folding
Overall Width w/Arms:	21"
Weight Capacity:	250 lbs.
Product Weight:	19 lbs.
Shipping Weight:	23 lbs.
Warranty:	Limited 3 Year

Item #	Description	Color
TC19A-BLK	Astra Transport Chair 19" Aluminum w/Nylon Seat - 250 lb. Wt-Cap.	Black
TC19A-BLU	Astra Transport Chair 19" Aluminum w/Nylon Seat - 250 lb. Wt-Cap.	Blue