

Chariot II Wheelchair



The **Chariot II Standard Hemi Wheelchair** is designed to meet a variety of patient needs. With a weight capacity of 300-350 lbs., a dual axle for hemi-height adjustment and multi-position caster forks, it is the most versatile wheelchair on the market. Standard features include removable padded armrests, easy-to-clean vinyl upholstery, chart pocket and heel loops on swing-away footrests. Available with swing-away footrests or elevating legrests.



Removable armrests



Easy-to-operate wheel locks



Swing-away footrests with heel loops

FEATURES

- Dual axle provides easy transition of seat height to hemi-level
- Powder-coated silver vein steel frame
- Includes either tool-free swing-away footrests with heel loops or padded elevating legrests with composite footplates
- Elevating legrests are equipped with a pneumatic piston for infinite elevating positions
- Removable, padded, upholstered armrests
- Easy-to-clean, comfortable embossed vinyl upholstery resists mildew and bacteria
- Easy push-to-lock wheel locks
- Precision sealed wheel bearings in front and rear provides long lasting performance and reliability
- Composite handrims with 24" solid urethane tires
- 8" x 1" solid rubber front casters adjustable in 3 positions
- Chart pocket on back for additional convenience
- HCPCS Code: K0002, K0195 OR E0990 (elevating legrests only)

SPECIFICATIONS

Seat Width:	16"/18"/20"
Seat Depth:	16"
Seat-to-Floor Height:	19.5"/Hemi 17.5"
Back Height:	18"
Back Style:	Fixed
Overall Width w/Arms:	24"/26"/28"
Weight Capacity:	300/300/350 lbs.
Product Weight:	42/43/44 lbs.
Shipping Weight:	46/47/48 lbs.
Warranty:	Limited 3 Year

Item #	Description
WCK216DAELR	Chariot II 16" K2 Wheelchair w/Removable Desk Arms, Elevating Leg Rests, Vinyl Seat - 300 lb. Wt-Cap.
WCK216DASF	Chariot II 16" K2 Wheelchair w/Removable Desk Arms, Swing Away Foot Rests, Vinyl Seat - 300 lb. Wt-Cap.
WCK216FAELR	Chariot II 16" K2 Wheelchair w/Removable Full Arms, Elevating Leg Rests, Vinyl Seat - 300 lb. Wt-Cap.
WCK216FASF	Chariot II 16" K2 Wheelchair w/Removable Full Arms, Swing Away Foot Rests, Vinyl Seat - 300 lb. Wt-Cap.
WCK218DAELR	Chariot II 18" K2 Wheelchair w/Removable Desk Arms, Elevating Leg Rests, Vinyl Seat - 300 lb. Wt-Cap.
WCK218DASF	Chariot II 18" K2 Wheelchair w/Removable Desk Arms, Swing Away Foot Rests, Vinyl Seat - 300 lb. Wt-Cap.
WCK218FAELR	Chariot II 18" K2 Wheelchair w/Removable Full Arms, Elevating Leg Rests, Vinyl Seat - 300 lb. Wt-Cap.
WCK218FASF	Chariot II 18" K2 Wheelchair w/Removable Full Arms, Swing Away Foot Rests, Vinyl Seat - 300 lb. Wt-Cap.
WCK220DAELR	Chariot II 20" K2 Wheelchair w/Removable Desk Arms, Elevating Leg Rests, Vinyl Seat - 350 lb. Wt-Cap.
WCK220DASF	Chariot II 20" K2 Wheelchair w/Removable Desk Arms, Swing Away Foot Rests, Vinyl Seat - 350 lb. Wt-Cap.
WCK220FAELR	Chariot II 20" K2 Wheelchair w/Removable Full Arms, Elevating Leg Rests, Vinyl Seat - 350 lb. Wt-Cap.
WCK220FASF	Chariot II 20" K2 Wheelchair w/Removable Full Arms, Swing Away Foot Rests, Vinyl Seat - 350 lb. Wt-Cap.

270 Washington Street, Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) * Fax: 845-205-6005

www.proactivemedical.com * sales@proactivemedical.com