

Chariot IV XTC Pediatric Wheelchair 12" / 14"



FEATURES

- Chariot IV XTC Pediatric weighs 32 lbs. (excluding front rigging)
- Easier to propel and transport
- Frame style eliminates seat guides and allows for custom backs and accessories
- Adjustable height back adjusts from 16" to 18"
- Carbon steel black powder-coated frame double coated for an attractive, chip-proof, maintainable finish
- Durable, flame retardant nylon upholstery resists mildew and bacteria
- Composite 20" mag-style wheels with aluminum hand rims are lightweight and maintenance free
- Urethane rear tires and casters offer superior performance, smoother ride, and are lightweight
- 6" front casters have three height adjustments
- Nylon, padded armrests provide added patient comfort
- Flip-back arms allow for easy transfer
- Composite steel lined front receiver is angled so the arm will lock into place every time

The **Chariot IV XTC Pediatric Wheelchair** is a versatile wheelchair right out of the box and is easier to propel and transport than an average weight wheelchair. The Chariot IV XTC Pediatric comes standard with padded arms, multiple seat to floor height positions, anti-tippers with flip-up wheels which never have to be removed, and push-to-lock wheel locks. The back is height adjustable to provide added patient comfort. The carbon steel frame is high strength while remaining lightweight. The wheelchair upholstery is made of durable nylon for long lasting comfort and cleanliness.

SPECIFICATIONS

Seat Width:	12"/14"
Seat Depth:	14"
Back Height:	16"-18"
Height 12"/14":	32"
Length w/Riggings 12"/14":	36"
Weight Capacity:	250/250 lbs.
Product Weight:	32 lbs.
Shipping Weight:	37 lbs.
Warranty:	Limited Lifetime

- Front mounted release lever can be flipped up or down to unlock flip-up arm
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- Dual axle provides easy transition of seat height to hemi-level
- Elevating leg rests with tool-free adjustable length riggings
- Composite foot plates with heel loops are durable and lightweight
- Comes standard with flip-up anti-tippers
- Comes with push-to-lock wheel locks
- HCPCS Code: K0004+E1236+K0195 or E0990

Item #	Description
WCK4X12DAELR	Chariot IV XTC 12" K4 Wheelchair w/Flip-Back Desk Arms, Elevating Leg Rests, Nylon Seat, Anti-Tippers - 250 lb. Wt-Cap.
WCK4X14DAELR	Chariot IV XTC 14" K4 Wheelchair w/Flip-Back Desk Arms, Elevating Leg Rests, Nylon Seat, Anti-Tippers - 250 lb. Wt-Cap.

270 Washington Street, Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) * Fax: 845-205-6005

www.proactivemedical.com * sales@proactivemedical.com