



FEATURES

- State-of-the art hydraulic reclining mechanism
- Backrest reclines from 90 to 180 degrees
- Double-coated chrome carbon steel frame for maximum durability
- Padded flip-back desk arms
- Composite lightweight 20" mag-style rear wheels with urethane tires and 6" front rubber tires
- Offset rear wheel position provides additional stability and evenly distributes users weight in all reclining positions
- Precision sealed front and rear wheel bearings
- Front caster forks are adjustable in three positions
- Standard dual axle provides easy transition of seat to hemi-level
- Push-to-lock wheel locks
- Spreader bar functions as a back stabilizer and a push handle for caregiver
- Standard with height adjustable swing-away elevating leg rests
- Durable and lightweight composite footplates
- Flip-up anti-tippers with wheels
- HCPCS Code: E1236+E1226+K0195 or E0990

Chariot IV XTC-RC Pediatric 12"/14" Reclining Wheelchair

The **Chariot IV XTC-RC Pediatric Reclining Wheelchair** has a headrest that extends past the backrest and is detachable. This reclining wheelchair also features a spreader bar that acts as a back stabilizer and a push handle for an attendant.

The Chariot IV XTC-RC Pediatric Wheelchair has offset rear wheels to add additional stability to evenly distribute the patient's weight. The rear wheels are supported by composite mag-style wheels to shave off added weight, while in front are two 6 inch wheels. The rear wheels include locks conveniently located next to the chair.

This chair is equipped with a hydraulic reclining mechanism to allow for adjustments up to 180 degrees. The seat is upholstered with flame-retardant nylon which is resistant to mildew and bacteria. The frame is constructed of carbon steel and has double-coated chrome and a chip-proof finish. The arm rests are padded with durable stitched nylon to provide an additional level of comfort and will flip back to aid in transferring in and out of the chair. It also comes with composite foot plates that are durable and lightweight.

SPECIFICATIONS

Seat Width:	12"/14"
Seat Depth:	12"/14"
Seat-to-Floor Height:	15"/16"
Back Height:	26"
Back Style:	Reclining
Weight Capacity:	250/250 lbs.
Product Weight:	41 lbs.
Shipping Weight:	46 lbs.
Warranty:	Limited Lifetime

Item #	Description
WCK4XRC12DAELR	Chariot IV XTC-RC 12" K4 Reclining Wheelchair, Flip-Back Desk Arms, Elevating Leg Rests, Nylon Seat, Anti-Tippers - 250 lb. Wt-Cap.
WCK4XRC14DAELR	Chariot IV XTC-RC 14" K4 Reclining Wheelchair, Flip-Back Desk Arms, Elevating Leg Rests, Nylon Seat, Anti-Tippers - 250 lb. Wt-Cap.

270 Washington Street, Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) * Fax: 845-205-6005

www.proactivemedical.com * sales@proactivemedical.com