

# Chariot IV XTC-RC Reclining Wheelchair



The **Chariot IV XTC-RC Reclining Wheelchair** has a padded headrest that extends past the backrest and is detachable. The seat back extends to cover the full back and comes in optional widths, from 16 to 22 inches. This reclining wheelchair also features a spreader bar that acts as a back stabilizer and a push handle for an attendant.

The Chariot IV XTC-RC has offset rear wheels to add additional stability to evenly distribute the patient's weight. The rear wheels are supported by composite mag-style wheels to shave off added weight, while in front are two 8 inch wheels. The rear wheels include locks conveniently located next to the chair. This chair is equipped with a hydraulic reclining mechanism to allow for infinite adjustments up to 180 degrees. The frame is constructed of carbon steel and has double-coated chrome and a chip proof finish.

## FEATURES

- Built-in seat rail extension and extendable upholstery easily adjusts seat depth
- Attendant pull to lock wheel locks are located in a convenient and easily accessible location
- Spreader bar acts as a back stabilizer and a convenient push handle for the attendant
- Offset rear wheel position adds additional stability and equally distributes the patients weight while in the various reclining positions
- Hydraulic reclining mechanism allows for infinite adjustments up to 180 degrees
- Durable, flame retardant nylon upholstery resists mildew and bacteria
- Carbon steel frame with double-coated chrome for an attractive, chip-proof, maintainable finish
- Composite 24" mag-style wheels are lightweight and maintenance free, featuring 8" front casters
- Nylon padded armrests provide added patient comfort, and feature the flip-back style for easy transfers
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- Standard dual axle provides easy transition of seat-to-floor height from adult to hemi-level

## SPECIFICATIONS

|                       |                      |
|-----------------------|----------------------|
| Seat Width:           | 16"/18"/20"/22"      |
| Seat Depth:           | 16"/18"              |
| Seat-to-Floor Height: | 19.5"/Hemi 17.5"     |
| Back Height:          | 23"/33"              |
| Back Style:           | Reclining            |
| Overall Width w/Arms: | 24"/26"/28"/29"      |
| Weight Capacity:      | 300/300/300/300 lbs. |
| Product Weight:       | 47/49/51/53 lbs.     |
| Shipping Weight:      | 52/53/56/58 lbs.     |
| Warranty:             | Limited Lifetime     |

- Front caster forks are adjustable in three positions
- Comes standard with swing-away elevating leg rests
- Composite foot plates are durable and lightweight
- Urethane rear tires and casters offer superior performance, smoother ride, and are lightweight
- HCPCS Code: (16"-20" Chairs) K0004+E1226+K0195 or E0990  
(22" Chairs) K0007+E1226+K0195 or E0990

| Item #         | Description  |
|----------------|--|
| WCK4XRC16DAELR | Chariot IV XTC-RC 16" K4 Reclining Wheelchair w/Removable Flip-Back Desk Arms, ELR, Nylon Seat, Seat Extension - 300 lb. Wt-Cap. |
| WCK4XRC16FAELR | Chariot IV XTC-RC 16" K4 Reclining Wheelchair w/Removable Flip-Back Full Arms, ELR, Nylon Seat, Seat Extension - 300 lb. Wt-Cap. |
| WCK4XRC18DAELR | Chariot IV XTC-RC 18" K4 Reclining Wheelchair w/Removable Flip-Back Desk Arms, ELR, Nylon Seat, Seat Extension - 300 lb. Wt-Cap. |
| WCK4XRC18FAELR | Chariot IV XTC-RC 18" K4 Reclining Wheelchair w/Removable Flip-Back Full Arms, ELR, Nylon Seat, Seat Extension - 300 lb. Wt-Cap. |
| WCK4XRC20DAELR | Chariot IV XTC-RC 20" K4 Reclining Wheelchair w/Removable Flip-Back Desk Arms, ELR, Nylon Seat, Seat Extension - 300 lb. Wt-Cap. |
| WCK4XRC20FAELR | Chariot IV XTC-RC 20" K4 Reclining Wheelchair w/Removable Flip-Back Full Arms, ELR, Nylon Seat, Seat Extension - 300 lb. Wt-Cap. |
| WCK4XRC22DAELR | Chariot IV XTC-RC 22" K4 Reclining Wheelchair w/Removable Flip-Back Desk Arms, ELR, Nylon Seat, Seat Extension - 300 lb. Wt-Cap. |

270 Washington Street, Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) \* Fax: 845-205-6005

www.proactivemedical.com \* sales@proactivemedical.com