


Chariot-RC Wheelchair

The **Chariot-RC Reclining Wheelchair** provides full back and headrest support as patient reclines from 90 to 180 degrees. The headrest and padded armrests are removable. Black, padded nylon upholstery resists mildew and bacteria. The frame has a chip-resistant silver vein finish. Mag-style wheels are lightweight and set back on frame to prevent user from tipping over when in the reclined position.



Reclines from
90° - 180°



Recliner adjustment
mechanism



Headrest supports
patient when in the
reclined position

FEATURES

- Slide tube steel frame with a silver vein finish that is chip resistant and easy to maintain
- Padded, detachable armrests provide added patient comfort
- Easy-to-clean, comfortable nylon upholstery resists mildew and bacteria
- Tool-free elevating padded legrests are equipped with a pneumatic piston for infinite positions
- Composite footplates on 16"/18"/20" & aluminum footplates on 22"/24" chairs
- Composite handrims with mag-style wheels are lightweight and maintenance free
- Wheels are set back on frame to prevent tipping
- Rear anti-tippers come standard
- Headrest extension comes standard with removable head immobilizer
- 8" x 1" solid rubber front casters
- Easy push-to-lock wheel locks
- HCPCS Code: K0001, E1226 (reclining back), E0990 (elevating legrests only)

SPECIFICATIONS

Seat Width:	16"-18"-20" / 22"-24"
Seat Depth:	16" / 18"
Seat-to-Floor Height:	20"
Back Height:	32"
Back Style:	Reclining
Overall Width w/Arms:	24"/26"/28"/30"/32"
Weight Capacity:	300/300/350/500/500 lbs.
Product Weight:	58/59/60/61/62 lbs.
Shipping Weight:	64/65/66/67/68 lbs.
Warranty:	Limited 3 Year

Item #	Description
WCRC16DAELR	16" Reclining Wheelchair w/Removable Desk Arms, Rear Anti-Tippers, Elevating Leg Rests, Nylon Seat - 300 lb. Wt-Cap.
WCRC16FAELR	16" Reclining Wheelchair w/Removable Full Arms, Rear Anti-Tippers, Elevating Leg Rests, Nylon Seat - 300 lb. Wt-Cap.
WCRC18DAELR	18" Reclining Wheelchair w/Removable Desk Arms, Rear Anti-Tippers, Elevating Leg Rests, Nylon Seat - 300 lb. Wt-Cap.
WCRC18FAELR	18" Reclining Wheelchair w/Removable Full Arms, Rear Anti-Tippers, Elevating Leg Rests, Nylon Seat - 300 lb. Wt-Cap.
WCRC20DAELR	20" Reclining Wheelchair w/Removable Desk Arms, Rear Anti-Tippers, Elevating Leg Rests, Nylon Seat - 350 lb. Wt-Cap.
WCRC20FAELR	20" Reclining Wheelchair w/Removable Full Arms, Rear Anti-Tippers, Elevating Leg Rests, Nylon Seat - 350 lb. Wt-Cap.
WCRC22DAELR	22" Reclining Wheelchair w/Removable Desk Arms, Rear Anti-Tippers, Elevating Leg Rests, Nylon Seat - 500 lb. Wt-Cap.
WCRC22FAELR	22" Reclining Wheelchair w/Removable Full Arms, Rear Anti-Tippers, Elevating Leg Rests, Nylon Seat - 500 lb. Wt-Cap.
WCRC24DAELR	24" Reclining Wheelchair w/Removable Desk Arms, Rear Anti-Tippers, Elevating Leg Rests, Nylon Seat - 500 lb. Wt-Cap.

270 Washington Street, Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) * Fax: 845-205-6005

www.proactivemedical.com * sales@proactivemedical.com