


## Transport Chair/Wheelchair Combo

The Protekt Cambria™ High Strength Lightweight Convertible Transport Chair/Wheelchair can be used as a standard, self-propelled wheelchair or as a transport chair, all-in-one for your convenience. The quick release 24" wheels can be easily removed to transition from a wheelchair to a transport chair. Also comes with two sets of aluminum wheel locks, one for use with wheelchair and one for use with transport chair, for added safety. The deluxe back release folds down the backrest for easy and convenient storage and transportation. Dual brakes, desk length removable armrests, a fold-down backrest and swing-away removable footrests all combine to make for a portable and high-quality chair.



### FEATURES

- **Dual functionality:** Easily converts between self-propelled wheelchair and transport chair
- **Lightweight and sturdy:** High-strength, all-aluminum frame with quick-release wheels
- **Versatile braking system:** Two sets of brakes for different chair configurations
- Quick release 24" wheels can be removed to transition self-propelled chair to a transport chair
- Includes both 24" and 8" rear wheels
- Comes standard with swing away footrests with composite footplates and heel loops
- Swing-away footrests and anti-tippers included
- Fold down back for easy storage
- Nylon easy to clean upholstery
- 8" casters provide a smooth, sturdy, and secure ride
- HCPCS Code: K0004/E0951/E0971

SPECIFICATIONS	TRANSPORT CHAIR	WHEELCHAIR
Seat Width:	16"/18"/20"	16"/18"/20"
Seat Depth:	16"	16"
Seat-to-Floor Height:	19.5"	19.5"
Back Height:	16"	16"
Rear Wheel:	8"	24"
Overall Width:	21"-25.5"	24.5"-28.5"
Back Style:	Fold-down	Fold-down
Weight Capacity:	275 lbs.	275 lbs.
Product Weight:	23/24/25 lbs.	31.5/32.5/33.5 lbs.
Warranty:	Limited 3 Year	

Item #	Description
WCTC16SF	Protekt Cambria™ 16" Wheelchair/Transport Chair Combo, Desk Arms, Swing Away Foot Rests - 275 lb. Wt-Cap.
WCTC18SF	Protekt Cambria™ 18" Wheelchair/Transport Chair Combo, Desk Arms, Swing Away Foot Rests - 275 lb. Wt-Cap.
WCTC20SF	Protekt Cambria™ 20" Wheelchair/Transport Chair Combo, Desk Arms, Swing Away Foot Rests - 275 lb. Wt-Cap.