

Titus Bariatric Wheelchair

The **Titus Heavy-Duty Wheelchair** features an ultra sturdy design to accommodate patients with wider seat and higher weight capacity requirements. With seat widths ranging from 22 to 30 inches, a weight capacity of 500-700 pounds and a heavy-duty slide tube frame, the Titus offers a substantial yet attractive option for the bariatric patient population.



Removable, padded arms



Heavy-duty upholstery



Aluminum footplates

FEATURES

- Heavy-duty, slide tube silver vein steel frame
- 8" x 2" semi-pneumatic casters are adjustable in 3 positions
- Standard dual axle (22"/24") provides easy transition of seat to hemi-level
- Optional elevating legrests or swing-away footrests
- 24" rear urethane wheels with durable hand rims
- Removable arms with padded armrests
- Heavy-duty padded upholstery
- Tool-free elevating padded legrests are equipped with a pneumatic piston for infinite elevating positions
- Extra wide aluminum footplates on 28" & 30" chairs
- Easy push-to-lock wheel locks
- HCPCS Code: K0007 or K0195

SPECIFICATIONS

Seat Width:	22"-24" / 26"-28"-30"
Seat Depth:	18" / 20"
Seat-to-Floor Height:	20.5"/Hemi 18.5"
Back Height:	16"/18"
Back Style:	Fixed
Overall Width w/Arms:	30"/32"/36"/38"/40"
Weight Capacity:	500/500/700/700/700 lbs.
Product Weight:	56/57/59/61/64 lbs.
Shipping Weight:	62/63/66/68/70 lbs.
Warranty:	Limited 3 Year

Item #	Description
WCK722DAELR	Titus 22" Heavy Duty Wheelchair w/Removable Desk Arms, Elevating Leg Rests, Vinyl Seat - 500 lb. Wt-Cap.
WCK722DASF	Titus 22" Heavy Duty Wheelchair w/Removable Desk Arms, Swing Away Foot Rests, Vinyl Seat - 500 lb. Wt-Cap.
WCK722FAELR	Titus 22" Heavy Duty Wheelchair w/Removable Full Arms, Elevating Leg Rests, Vinyl Seat - 500 lb. Wt-Cap.
WCK722FASF	Titus 22" Heavy Duty Wheelchair w/Removable Full Arms, Swing Away Foot Rests, Vinyl Seat - 500 lb. Wt-Cap.
WCK724DAELR	Titus 24" Heavy Duty Wheelchair w/Removable Desk Arms, Elevating Leg Rests, Vinyl Seat - 500 lb. Wt-Cap.
WCK724DASF	Titus 24" Heavy Duty Wheelchair w/Removable Desk Arms, Swing Away Foot Rests, Vinyl Seat - 500 lb. Wt-Cap.
WCK724FAELR	Titus 24" Heavy Duty Wheelchair w/Removable Full Arms, Elevating Leg Rests, Vinyl Seat - 500 lb. Wt-Cap.
WCK724FASF	Titus 24" Heavy Duty Wheelchair w/Removable Full Arms, Swing Away Foot Rests, Vinyl Seat - 500 lb. Wt-Cap.
WCK726DAELR	Titus 26" Heavy Duty Wheelchair w/Removable Desk Arms, Elevating Leg Rests, Nylon Seat - 700 lb. Wt-Cap.
WCK726DASF	Titus 26" Heavy Duty Wheelchair w/Removable Desk Arms, Swing Away Foot Rests, Nylon Seat - 700 lb. Wt-Cap.
WCK726FAELR	Titus 26" Heavy Duty Wheelchair w/Removable Full Arms, Elevating Leg Rests, Nylon Seat - 700 lb. Wt-Cap.
WCK726FASF	Titus 26" Heavy Duty Wheelchair w/Removable Full Arms, Swing Away Foot Rests, Nylon Seat - 700 lb. Wt-Cap.
WCK728DAELR	Titus 28" Heavy Duty Wheelchair w/Removable Desk Arms, Elevating Leg Rests, Nylon Seat - 700 lb. Wt-Cap.
WCK728DASF	Titus 28" Heavy Duty Wheelchair w/Removable Desk Arms, Swing Away Foot Rests, Nylon Seat - 700 lb. Wt-Cap.
WCK728FAELR	Titus 28" Heavy Duty Wheelchair w/Removable Full Arms, Elevating Leg Rests, Nylon Seat - 700 lb. Wt-Cap.
WCK728FASF	Titus 28" Heavy Duty Wheelchair w/Removable Full Arms, Swing Away Foot Rests, Nylon Seat - 700 lb. Wt-Cap.
WCK730DAELR	Titus 30" Heavy Duty Wheelchair w/Removable Desk Arms, Elevating Leg Rests, Nylon Seat - 700 lb. Wt-Cap.
WCK730DASF	Titus 30" Heavy Duty Wheelchair w/Removable Desk Arms, Swing Away Foot Rests, Nylon Seat - 700 lb. Wt-Cap.
WCK730FAELR	Titus 30" Heavy Duty Wheelchair w/Removable Full Arms, Elevating Leg Rests, Nylon Seat - 700 lb. Wt-Cap.
WCK730FASF	Titus 30" Heavy Duty Wheelchair w/Removable Full Arms, Swing Away Foot Rests, Nylon Seat - 700 lb. Wt-Cap.

270 Washington Street, Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) * Fax: 845-205-6005

www.proactivemedical.com * sales@proactivemedical.com