

# Protekt® Two-Button Folding Lightweight Aluminum Walker

The **Protekt® Lightweight Aluminum Walkers** are an ideal and dependable solution for those who require a supportive and durable walking aid for their daily activities. The easy-to-use push button mechanism allows for opening and closing with the user's fingers or palms. The walker's height can be adjusted with a tool-free push button, providing additional stability. Available with or without 5" wheels and in junior or adult sizes. Supports patient weights up to 350 pounds.



Rear glide caps standard on walker with wheels



Push button provides audible click when locked



Height adjustable

## FEATURES

- Attractive, lightweight Aluminum frame construction provides superior strength and durability
- Easy-to-use push button folding mechanism provides audible click when locked
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing
- 5" wheels provide added stability & support
- Front cross brace allows for full user stride
- Side bar adjusts with height to provide additional support
- Adjustable height settings
- Vinyl contoured hand grip for added comfort and long lasting wear
- Rear glide cap allows for walker to slide easily and smoothly over most surfaces (come standard on models with 5" wheels)
- Folds flat for transport and storage

## SPECIFICATIONS

Width Inside Hand Grips:	17"	
Width Inside Base Legs:	Adult 19.5" / Junior 19"	
Depth Opened at Base:	Adult 17" / Junior 16"	
Depth Folded:	4.5"	
Wheel Size:	5"	
Patient Height Range:	Junior:	4'4" - 5'7"
	Adult:	5'2" - 6'4"
Weight Capacity:	350 lbs.	
Product Weight:	PM1051A	6.5 lbs.
	PM1052A	7.5 lbs.
	PM1052AJ	5 lbs.
Carton Dimensions:	PM1051A	33" x 21.5" x 9.5"
	PM1052A	30.5" x 22" x 9.5"
	PM1052AJ	24" x 21" x 10.5"
Shipping Weight:	PM1051A	34.5 lbs.
	PM1052A	36.5 lbs.
	PM1052AJ	27 lbs.
Warranty:	Limited Lifetime	

Item #	Description	HCPSC Code
PM1051A	Protekt Walker Adult Without Wheels, Aluminum - 350 lb. Wt-Cap.	E0135 (without wheels)
PM1052A	Protekt Walker Adult With Wheels, Aluminum - 350 lb. Wt-Cap.	E0143 (with wheels)
PM1052AJ	Protekt Walker Junior With Wheels, Aluminum - 350 lb. Wt-Cap.	E0143 (with wheels)