

# PROTEKT® 350 HYDRAULIC STAND

**HCPCS  
EO630/  
EO635**

**Hydraulic Sit-To-Stand  
Patient Lift \*Optional  
Electric Upgrade**

Model: 34350H

**Electric  
Upgrade Kit**

Model:  
LP-34350HE-K

**Hydraulic Sit-To-Stand  
Patient Lift with Electric  
Upgrade**

Model: 34350HE

The Protekt® 350 Hydraulic Stand is Proactive's answer to at home safe patient handling. It is easy to maneuver, transport, and maintain while providing safe and secure patient transfers. Designed for patients who are partially weight bearing, cognitively motivated and need rehabilitation support. The hydraulic pump style provides a low-cost option for caregivers to easily care for patients. An optional electric kit upgrades the hydraulic system to a full powered electric lift.

- **Durable powder coated steel tube construction.**
- **Adjustable cushioned leg supports with lateral contours provide security and comfort.**
- **Low friction casters and wide push handle for easy movement.**
- **Manual open/close base legs.**
- **Retractable, non-slip footplate provides a stable base for transfers.**
- **Multi-function slings for transfers and toileting.**
- **Under-bed clearance of 4.5" allows use with most home care beds.**
- **Pump handle can rotate to either side for caregiver convenience.**
- **3" front casters/5" rear casters with locking brakes.**
- **Optional upgrade to battery power.**
- **350 lb. weight capacity.**

**WARRANTY: Electronics 2 years  
Frame 3 years  
Battery 1 year**

## SPECIFICATION FOR ITEM # 34350H

Weight Capacity	350 lbs.
Unit Weight	84 lbs.
External Base Width - Closed	24"
Internal Base Width - Opened	36.6"
Upright Mast Height	42.3"
Overall Base Length	37.4"
Base Clearance	4.5"
Minimum Lifting Height with Spreader Bar	42.7"
Maximum Lifting Height with Spreader Bar	62.8"
Maximum Boom Width	19.7"
Knee Pad Height Extension Range	24"
Shipping Weight	93.5 lbs.

Battery	LP-34350-BPAK
Charging Power Adapter	LP-34350-ADP
Hand Pendant	LP-34350-PEN



**34350H**



**ELECTRIC  
UPGRADE KIT**



**34350HE**