

Protekt® Aluminum Rollator

The **Protekt® Aluminum Rollator** is designed to fit the majority of your patients' needs with a focus on style, affordability, and independence. Constructed of lightweight aluminum, it offers great maneuverability and stability. Equipped with easy-to-operate ergonomic hand brakes, a removable padded backrest and a water-resistant nylon pouch. Easily folds for storage and transport with built-in fold strap. Supports patient weights up to 300 pounds.



Easy-to-operate
Ergonomic
Hand Brakes



6" Wheels



Water-Resistant
Nylon Pouch

FEATURES

- Lightweight aluminum frame available in Blue or Black
- Two-button removable curved padded backrest
- Comfortable, flip-up padded seat
- Easy-to-operate ergonomic hand brakes
- Height-adjustable for a custom fit
- Nylon pouch attaches under seat for added privacy and security of personal items
- Equipped with a folding mechanism for easy storage and transport
- 6" black casters with soft grip tires, ideal for both indoor and outdoor use

SPECIFICATIONS

Overall Width:	24"
Width b/w Hand Grips:	20"
Seat Dimensions:	13.5"(W) X 12"(D)
Seat-to-Floor Height:	20"
Weight Capacity:	300 lbs.
Adjustable Handle Height:	31" - 36.5"
Wheel Size:	6"
Product Weight:	13 lbs.
Carton Dimensions:	31.5" x 24.5" x 7.5"
Shipping Weight:	16 lbs.
Warranty:	Limited Lifetime

COLORS

Blue
Finish



Black
Finish



Item #	Description	Color	HCPCS Code
RLADAL6-BLU	Protekt Rollator, Aluminum Standard	Blue	E0143 + E0156
RLADAL6-BLK	Protekt Rollator, Aluminum Standard	Black	E0143 + E0156