

# Protekt® Bariatric Two-Button Folding Walker



The **Protekt® Bariatric Two-Button Folding Walkers** provide maximum strength and durability for the larger individual. This walker is an ideal and dependable solution for those who require a supportive yet durable aid for their daily activities. The easy-to-use push button mechanism allows for opening and closing with the user's fingers or palms. The walker's height can be adjusted with a tool-free push button, providing additional stability. Available with or without 5" wheels. Supports patient weights up to 500 pounds.



Rear glide caps standard on walker with wheels



Push button provides audible click when locked



Height adjustable

## FEATURES

- Durable steel frame
- Adjustable height settings
- Easy-to-use push button folding mechanism provides audible click when locked
- Front cross brace allows for full user stride
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing
- Vinyl contoured hand grip for added comfort and long lasting wear
- Rear glide cap allows for walker to slide easily and smoothly over most surfaces (come standard on models with 5" wheels)
- Folds flat for transport and storage

## SPECIFICATIONS

Width Inside Hand Grips:	20"
Width Inside Base Legs:	23.5"
Depth Opened at Base:	17"
Depth Folded:	4.5"
Wheel Size:	5"
Patient Height Range:	Adult: 5'2"-6'4"
Weight Capacity:	500 lbs.
Product Weight:	W/out Wheels: 12 lbs. With Wheels: 13 lbs.
Carton Dimensions:	W/out Wheels: 34" x 26.5" x 6.5" With Wheels: 29" x 27" x 6.5"
Shipping Weight:	W/out Wheels: 26.5 lbs. With Wheels: 28 lbs.
Warranty:	Limited Lifetime

Item #	Description	HCPCS Code
PM1051B	Protekt Bariatric Walker Without 5" Wheels, Steel - 500 lb. Wt-Cap.	E0148 (without wheels)
PM1052B	Protekt Bariatric Walker With 5" Wheels, Steel - 500 lb. Wt-Cap.	E0149 (with wheels)