

PROTEKT® STS COMPACT 500

34500-COM | SIT-TO-STAND PATIENT LIFT

The Proakt® STS Compact 500 Patient Lift is a user friendly, durable, and affordable sit-to-stand lift. It is extremely compact and easy to maneuver and transport anyone 500 lbs. or under. It has a low base height of 4.9" which is ideal for sliding under most furniture or chairs. The foot pedal easily opens the base to allow wheelchair or toilet entrance.

The 34500-COM is intended for patients who are lacking leg strength in their lower extremities but are still able to bear weight. It is designed to assist with standing from a seated position while using a toilet/commode, changing clothes and/or to relieve pressure. The 34500-COM is intended to have one caregiver operate the lift to perform the transfer.

The 34500-COM allows the patient to be lifted and transferred safely with minimal physical effort provided by the operator. The 34500-COM uses an electric actuator to lift patients from a sitting to a standing position from either a bed or chair for transferring. The lift has a standing platform to be used while transferring the patient from one surface to another, also included is a cushioned knee pad for additional comfort and support.

The performance class actuator includes an emergency quick release in case of a power outage.

**HCPCS
E0637**

**WARRANTY: Electronics 2 years
Frame 3 years
Battery 1 year**

SPECIFICATION FOR ITEM # 34500-COM

Weight Capacity	500 lbs.
Unit Weight	111 lbs.
Base Ext. Width - Closed	24.9"
Mast Height	40"
Base Length	33.5"
Front Caster Height	4.7"
Knee Pad Depth Range	4"
Knee Pad Height Range	21", 20" or 19"
Footplate Width	13.9"
Footplate Length	14.4"
Base Opened Turning Diameter	45"
Base Closed Turning Diameter	41.9"

Battery	LP-TC12-BPAK
Charger (Plug)	LP-TP7C-ADP
Charging Cradle (No plug)	LP-TBC2-CRDL
Hand Pendant	LP-TH10-PEN



270 Washington Street, Mount Vernon, NY 10553
Tel: 855-BE-PROACTIVE (855-237-7622)
Fax: 845-205-6005
www.proactivemedical.com
sales@proactivemedical.com