

PROTEKT® TAKE-A-LONG

33400P | FOLDING ELECTRIC PATIENT LIFT

The Proakt® Take-A-Long Folding Electric Patient Lift represents a new standard in full body lifts for individuals who travel or have storage issues in their home. The small foot print is ideal for home use. It easily folds up to fit in most automobile trunks for transport or travel.

HCPCS
E0635



(Optional Scale
Item # 30300)

WARRANTY:
Electronics 2 years
Frame 3 years
Battery 1 year

SPECIFICATION FOR ITEM # 33400P

Weight Capacity	400 lbs.
Unit Weight	87.3 lbs.
Base Length	41.0"
Base Ext. Width - Closed	25.6"
Base Int. Width - Opened	32.0"
Mast Height	48.3"
Front Caster Height	4.7"
Lifting Height - Lowest	18.5"
Lifting Height - Highest	67.2"
Base Opened Turning Diameter	55.0"
Base Closed Turning Diameter	51.0"

Battery	LP-TC12-BPAK
Charger (Plug)	LP-TP7C-ADP
Charging Cradle (No plug)	LP-TBC2-CRDL
Hand Pendant	LP-TH10-PEN

- 2 point spreader bar.
- Portable and foldable for easy storage and transportation.
- Lightweight and easy to move.
- 4" rear casters with brake & 3" dual front casters.
- Tool-Free assembly and disassembly.
- Battery powered with a high-capacity battery pack.
- Low battery and overload warning alarm.
- Built in battery gauge.
- Emergency stop button.
- Optional spare battery and wall charger available.
- Optional digital scale.
- 400 lb. weight capacity.



270 Washington Street, Mount Vernon, NY 10553
Tel: 855-237-7622 • Fax: 845-205-6005
www.proactivemedical.com