

# PROTEKT® TRANSFORMER

## Hydraulic Patient Lift \*Optional Electric Upgrade

Model: 33400H

## Electric Upgrade Kit

Model: LP-33400HE-K

## Hydraulic Patient Lift with Electric Upgrade

Model: 33400HE

The Protekt® Transformer is Proactive's answer to at home safe patient handling. It is easy to maneuver, transport, and maintain while providing safe and secure patient transfers. The hydraulic pump style provides a low-cost option for caregivers to easily care for patients. An optional upgrade kit "transforms" the hydraulic system to a full powered electric lift. Given the demand for patient lifts in a variety of situations, the Protekt® Transformer is a low cost yet highly effective solution for patients in a homecare setting.

- Durable powder coated steel tube construction.
- Small & quick disassembly for transport.
- Six-point spreader bar.
- Manual open/close base legs.
- 4" rear casters with locking brakes.
- 3" front casters.
- Tool free assembly.
- 400 lb. weight capacity.
- Optional upgrade to battery power.
- Optional digital scale.

**WARRANTY:** Electronics 2 years  
Frame 3 years  
Battery 1 year

### SPECIFICATION FOR ITEM # 33400H/33400HE

	33400H	33400HE
Weight Capacity	400 lbs.	400 lbs.
Unit Weight	82 lbs.	82 lbs.
Base Length	40.3"	40.3"
Base Ext. Width - Closed	25.3"	25.3"
Base Int. Width - Opened	32.5"	32.5"
Mast Height	47.0"	47.0"
Front Caster Height	4.7"	4.7"
Lifting Height - Lowest	20.9"	25.3"
Lifting Height - Highest	68.2"	65.7"
Base Opened Turning Diameter	52.8"	54.7"
Base Closed Turning Diameter	49.5"	52.6"

Battery	LP-TC12-BPAK
Charger (Plug)	LP-TP7C-ADP
Charging Cradle (No plug)	LP-TBC2-CRDL
Hand Pendant	LP-TH10-PEN



33400H

33400HE