

vive  
health



## VIVE 4 Wheel Mobility Scooter

MOB1027

The Vive 4 Wheel mobility scooter provides safe and stable mobility assistance. With a turning radius of only 4.3' the scooter easily maneuvers around corners. 8" flat-free tires traverse over smooth indoor and outdoor surfaces such as concrete, pavement, tile, and carpeting. The padded seat features flip-back armrests and the ability to swivel and lock in place. The front mounted storage basket disconnects for optional use. The scooter has a speed control knob that adjusts up to 3.7mph and has a range of up to 12.4 miles per charge. Quickly disassemble into four pieces for easy transportation. The heaviest piece is only 29lbs.

Body Colors



Seat Color



## Best Economy in its class



### WHAT'S INCLUDE

- A. Front Section
- B. Rear Section
- C. Seat Unit
- D. Battery
- E. Charging Cable
- F. 4x Screw Knobs

- G. Basket
- H. Battery Securing Plate
- I. Basket Holder
- J. 2x Bolts
- K. Seat Pin
- L. Arm Rest

(888) 550-7750

aracent.com

sales@aracent.com



## Specifications

Overall Dimensions (L x W x H)	37.7" x 19.6" x 33"
Total Weight	88.2 lbs.
Weight of Heaviest Piece	29lbs. (rear base)
Seat Height (from ground)	18.5" - 21.5"
Seat Width	17"
Seat Depth	15"
Armrest Height (above seat)	6.75"
Seat Back Height	13"
Net Weight	88.2 lbs.
Maximum Speed	3.7 mph
Maximum Braking Distance	≤5"
Minimum Turning Radius	≤4.3ft.
Weight Limit	265 lbs.
Maximum Travel Distance (in Theory)	12.4 mi
Static Stability	≥ 9°
Dynamic Stability	≥ 6°
Climbing Ability	≤ 9°
Battery	Lead-acid battery (24V/20AH)
Controller Maximum Output Current	45A
Controller Minimum Output Current	2A
Front Wheel Diameter	7.5" (Solid)

### SAFELY NAVIGATE SMOOTH SURFACES:

Easy maneuver over smooth, indoor and outdoor surfaces with comfort and ease with the VIVE 4 Wheel mobility scooter. Designed to safely support up to 265 pounds, the mobility scooter features 8" flat-free, non-marking tires, a maximum speed of 3.7mph and a driving distance up to 12.4 miles per charge.

### COMFORTABLE SEATING:

Comfortably padded, the wide, 15.75" seat and backrest provide gentle cushioning and support. The seat rotates and is also height adjustable for a customizable fit. Additionally, the scooter seat is surrounded by two padded armrests that adjust and flip upward for greater convenience when exiting or boarding.

### EASY TO TRANSPORT:

Quickly disassemble the VIVE 4 Wheel mobility scooter into four separate pieces for easy portability. Featuring a folding tiller and removable seat, each section easily fits in the trunk of standard size vehicles. The rear section weighs 29 pounds and is the heaviest of the four pieces.

### INTUITIVE CONTROL PANEL:

Located on the tiller console, the control panel allows you to set the maximum speed with the variable speed dial, change from a forward to a reverse direction, turn on the headlight as well as produce a warning tone as needed.

The display also includes a color-coded battery life indicator for checking the status of the battery at a glance.