

vive health



VIVE

AIR FRAME ALUMINUM POWER WHEELCHAIR MOB1107AH12

The Air Frame Aluminum Power Wheelchair is a lightweight (44 pounds with battery) and portable mobility solution. Its foldable design allows for easy transport and storage in most vehicle trunks.

With a 12Ah battery, it offers a range of up to 9 miles on a single charge, ideal for daily use. The adjustable-sensitivity joystick and tight turning radius provide precise maneuverability, while multiple speed settings enhance performance indoors and outdoors. Its reinforced aluminum frame supports up to 225 pounds for durability and reliability.

WHAT'S INCLUDED

- 1x Airframe Power Wheelchair
- 1x 24V12AH battery

**COMPACT, FOLDABLE,
AND LIGHTWEIGHT
ENGINEERED FOR EVERYDAY
MOBILITY.**



aracent.com



(888) 550-7750



sales@aracent.com

SPECIFICATIONS

Max. Weight Capacity	225lbs
Seat Size	18.1" x 17.1" x 18.5"
Obstacle Height Climbing	≤1.6Inch
Maximum Speed	3.73 mph
Estimated Range	10 miles
Braking Distance at Flat Surface	≤59.1 Inch
Turning Radius	≤32"
Motor Type	24V 180W
Battery	DC 25.2V6AH Lithium
Brake	Electronic Brakes
Front Tire	7"
Rear Tire	12"
Wheelchair Size (Unfold)	39" x 24.4" x 36.6"
Folding Size	31.5" x 11" x 29"
Weight	39 lbs
Battery Type	Li-ion battery
Voltage	25.2V
Capacity	6Ah
Cut-off voltage	21V

EASY TO STORE AND TRANSPORT

At just 44 lbs with the battery, this electric wheelchair is light enough to lift and fold with ease. Its compact build fits into most car trunks, making it ideal for everyday mobility and travel.

FOLDABLE AND SPACE-SAVING

Designed for convenience on the go, this portable wheelchair folds down in seconds. Perfect for storing in tight spaces at home or while traveling.

EXTENDED BATTERY LIFE

Powered by a 12AH battery, the Airframe delivers up to 9 miles on a single charge. Enjoy reliable, long-range mobility without frequent charging stops.

RESPONSIVE JOYSTICK STEERING

Control your ride with an intuitive joystick featuring adjustable sensitivity. Effortlessly navigate tight spaces and uneven surfaces with a tight turning radius and multiple speed settings.

SUPPORTIVE AND LASTING COMFORT

Stay comfortable with a padded seat and durable frame built to hold up to 225 lbs. This wheelchair is engineered for smooth, stable use both indoors and outdoors.

