

vive health



VIVE

AIR FRAME WHEELCHAIR MOB1062BLKBLU

Discover the super lightweight folding manual wheelchair, a game-changing mobility solution that weighs only 28lbs, ensuring effortless travel and convenience. The adjustable footrest on the Vive Air Frame Wheelchair caters to various adult heights, ensuring a personalized fit. Built with durable aluminum and 24" wheels, it offers a stable ride, while the comfortable 18" W x 15.5" deep seat enhances comfort. Safety features like a foot securing strap, seatbelt, hand brakes, and rear brakes guarantee security while supporting up to 240 pounds, making it an ideal choice for a wide range of users.

WHAT'S INCLUDED

- Chair with seatbelt, leg straps, footrests, handbrakes.

COLORS

- Black with black and blue upholstery.

DURABLE CONSTRUCTION FOR PREMIUM COMFORT



aracent.com



(888) 550-7750



sales@aracent.com



Specifications

Product Dimensions:	40" x 38"x 26"
Seat width:	18"
Seat height:	19.5"
Seat deep:	15.5"
Wheels	24"
Product Weight:	28lbs
Weight capacity:	240 pounds
Weight limit	240 pounds



LIGHTWEIGHT FOLDING FRAME

The Vive Air Frame Wheelchair is a remarkable mobility solution designed to provide a portable way to regain your independence and safely navigate any environment. Weighing only 28 pounds, this durable yet lightweight chair revolutionizes your ability to travel, easily folding down for greater convenience when traveling or storing the device when not in use.

VERSATILE DESIGN WITH AN ADJUSTABLE FOOTREST

Enjoy personalized comfort with our adjustable footrest feature, built to accommodate various heights. This adaptable design ensures optimum comfort and support for a variety of users.

DURABLE CONSTRUCTION FOR PREMIUM COMFORT

Crafted from high-quality aluminum and equipped with 22" wheels, the Air Frame wheelchair guarantees a smooth and stable ride on diverse terrains. The 18" W by 15.5" deep seat adds to the comfort factor, making extended usage comfortable and enjoyable.

COMPREHENSIVE SAFETY FEATURES

The Vive Air Frame wheelchair comes with an array of safety features. From the foot securing strap and seat belt to the hand brakes and rear brakes, every detail is meticulously designed for maximum security and control. With a remarkable weight capacity of up to 240 pounds, this wheelchair assures both comfort and safety for a wide range of users.



8955 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808

SR23-4R2



aracent.com



(888) 550-7750



sales@aracent.com