

# vive health



## VIVE LIGHTWEIGHT ROLLATOR MOB1079

Crafted with precision and weighing only 14 pounds, the lightweight aluminum frame ensures effortless transportation and storage when not in use while safely accommodating up to 300 pounds. The ergonomic dual-braking system provides enhanced safety with secure wheel locking when not in use. Enjoy superior comfort and support with the padded seat and backrest, complete with a handy storage compartment for your essentials. Whether you are navigating indoor spaces or exploring the outdoors, the 6-inch wheels guarantee smooth movement, making it an ideal companion for travel or everyday use.

### WHAT'S INCLUDED

- 1 Lightweight Rollator

COLORS  

## 14 lbs LIGHTWEIGHT



[www.aracent.com](http://www.aracent.com)



(888) 550-7750



[sales@aracent.com](mailto:sales@aracent.com)



#### Specifications

|                            |                  |
|----------------------------|------------------|
| Width                      | 22" x 24"        |
| Adjustable Height          | 31" to 36"       |
| Seat                       | 14 x 12"         |
| Seat height from the floor | 21.5"            |
| Wheels                     | 6"               |
| Materials                  | Aluminum         |
| Product Weight:            | 14lbs            |
| Weight Capacity            | up to 300 pounds |

#### LIGHTWEIGHT ALUMINUM FRAME

The carefully crafted, lightweight aluminum frame of the Lightweight Rollator weighs a mere 14 pounds, allowing you to easily transport and store this rollator when it is not in use. The strength of the aluminum frame safely supports up to 300 pounds, for greater peace of mind.

#### ERGONOMIC DUAL-BRAKING SYSTEM

With comfortable handles that allow a natural grip and relieve tension on your wrists and hands, this rollator provides an ergonomic braking system that is easy to use. With dual braking, the wheels lock securely in place when not in use, for enhanced safety.

#### SUPERIOR COMFORT AND SUPPORT

Sit with greater comfort thanks to the padded seat and backrest. The integrated storage compartment allows you to keep your essentials close at hand, for greater convenience while on the go.

#### SMOOTH MOVEMENT INDOORS AND OUTDOORS

The 6" wheels are crafted to seamlessly glide over any surface, whether you are exploring new outdoor spaces or moving about a familiar indoor area. Take this rollator walker with you while traveling or use it in the comfort of your own home.



SR23-12R1



[www.aracent.com](http://www.aracent.com)



(888) 550-7750



[sales@aracent.com](mailto:sales@aracent.com)