

vive health



VIVE COMPACT POWER WHEELCHAIR MOB1029S

The Vive Compact Electric Wheelchair provides comfortable transportation at home and on the go. The power wheelchair features a dense, padded seat and a universal joystick controller that can be mounted on either the left or right armrest for added convenience. Equipped with dual motors, the electric wheelchair is powered with a rechargeable lithium battery that allows the wheelchair to navigate mild inclines and travel up to twelve miles per charge. Supporting up to 220 pounds, the durable motorized wheelchair frame is constructed with lightweight aluminum and is designed to easily fold for storage and travel.

WHAT'S INCLUDED

- Vive Compact Power Wheelchair
- Rechargeable battery
- 3-Pin charging cable

SAFELY AND SMOOTHLY NAVIGATING OVER ANY SURFACE



aracent..com



(888) 550-7750



sales@aracent.com

Specifications

Key Specs:

Top Speed	4 mph
Maximum Drive Range Per Charge	12 mi
Weight Capacity	220 lbs.
Total Weight	95 lbs.

General Specs:

Front Tires	7.5"
Rear Tires	15.7"
Maximum Speed	Up to 4 mph
Turning Radius	47"
Length	39.5"
Width	26"
Height	38"
Dimensions when folded	31" x 28.5" x 16.5"
Seat-to-Ground Height Range	20.5"
Seat Dimensions	16" x 18"
Seat back height	16.14"
Seat Material	Dense foam base with a breathable cover
Ground Clearance	4"
Total Weight	95 lbs.
Battery Requirements	24V 10AH
Maximum Range Per Charge	Up to 12miles
Battery Charger	On-board charger, 2A
Motor	24-volt/250W
Maximum Weight Capacity	220 lbs.



COMPACT FOR EASY TRAVEL

Perfect for traveling, the Vive compact electric wheelchair is constructed with a lightweight, yet durable aluminum frame that folds in seconds to easily fit in standard trunk sizes. The compact design is also TSA-approved and great for airline and cruise ship travel.

SAFE AND SMOOTH MANEUVERABILITY

Sturdy and stylish, the FDA-approved Vive power wheelchair provides comfortable and safe transportation at home and on the go. Compact and foldable, the electric wheelchair is equipped with powerful dual motors, flat-free front tires and large sport rear tires to safely and smoothly navigate over any terrain, including grass, pavement, gravel and more.

INTUITIVE JOYSTICK CONTROLLER

Easy to operate, the compact electric wheelchair features a 360-degree, joystick controller. The controller allows you to easily navigate around obstacles and monitor and adjust the speed of the wheelchair up to 4mph. The wheelchair controller is repositionable for use on the left or right side and includes a safety horn button, a battery charge indicator and a USB port for charging electronic devices while on the go.

COMFORTABLE SEATING

Designed for comfort, the motorized wheelchair features padded armrests and a wide, padded seat that safely supports up to 220 pounds. The cushioned seat measures 16" by 18" and is covered with a soft, breathable material that is easy to clean and includes an adjustable lap belt for added safety.

