

vive health



VIVE CORE KNEE WALKER MOB1099BLK HCPCS: E0118

The Core Vive Knee Walker is designed to provide unmatched mobility and support during your recovery, featuring large 8" wheels that ensure superior control and smooth steering for effortless navigation indoors and outdoors. Padded with soft hand grips and a cushioned knee pad, this knee walker offers ergonomic comfort, reducing strain on your knee, hands, and wrists for a pain-free experience. Built with a durable steel frame, it delivers exceptional stability and supports up to 300 pounds, giving you reliable security throughout your recovery journey. The adjustable handles and locking hand brakes allow for personalized comfort and enhanced safety, while the included storage bag keeps your essentials within easy reach.

WHAT'S INCLUDED

- 1x Vive Core Knee Walker
- 1x Storage Bag

COLOR

DURABLE STEEL FRAME FOR MAXIMUM SAFETY



www.aracent.com



(888) 550-7750



sales@aracent.com



Specifications

Dimensions	28" length, 15" width
Handle Bar Height	Adjustable in 1" increments ranging from 33" to 42"

Knee Platform Height	16.5" - 20.5"
----------------------	---------------

Wheels	8"
--------	----

Four 8-inch polyurethane wheels provide an exceptionally stable base for easy maneuverability both indoors and outdoors.

The handlebar folds down for easy storage.

Weight limit	Designed to support up to 300 pounds
--------------	--------------------------------------



BIG 8" WHEELS FOR BETTER CONTROL

Confidently navigate any environment using the Vive Knee Walker. The large 8" wheels ensure superior control and smooth steering, making it easier to maneuver indoors and outdoors.

PADDED FOR SUPERIOR COMFORT

Experience comfortable and padded support with soft hand grips and a cushioned knee pad. This ergonomic design reduces the strain and pressure placed on your knee, hands, and wrists, ensuring a more comfortable and pain-free recovery experience.

DURABLE STEEL FRAME FOR MAXIMUM SAFETY

Built from solid and durable steel, this knee walker provides a sturdy frame with exceptional stability. Designed to support up to 300 pounds, it provides dependable security throughout your recovery, giving you the confidence to move freely and safely.

ADJUSTABLE WITH CONVENIENT STORAGE

With adjustable handles, locking hand brakes, and an included storage bag, this knee walker ensures personalized comfort and enhanced safety. The knee walker is easy to assemble with detailed assembly instructions and comes with a storage bag to ensure you can keep essentials close at hand during recovery.



SR24-10R1



www.aracent.com



(888) 550-7750



sales@aracent.com