

vive health



VIVE CORE ROLLATOR

MOB1106BLK

HCPHS: E0143/E0156

This foldable rollator offers a compact and travel-ready design, easily collapsing to fit into most car trunks and tight storage spaces for hassle-free portability. Built with smooth-gliding wheels and a durable steel frame, it ensures confident navigation on both indoor and outdoor surfaces. Supporting up to 300 pounds, the sturdy steel construction provides reliable stability, while remaining lightweight at just 19 pounds for easy handling. The tool-free assembly makes setup quick and easy, and the adjustable height feature allows you to customize the rollator from 31 to 36 inches for optimal comfort. Ideal for those recovering from surgery, or injury, or anyone needing daily mobility support, this rollator also includes a convenient storage bag for carrying essentials.

SPECIFICATIONS

What's included:

➤ 1 Rollator, 1 Bag under the Seat

Product Material:

➤ Steel. Latex-free.

Dimensions:

➤ Height: 31-36"

➤ Width: 21"

➤ Length: 24"

Folded:

➤ Height: 37-41"

➤ Width: 21"

➤ Length: 13"

SEAT:

➤ 12X12"

Color:

➤ Black

Call now!

or send us an e-mail
now to find out more!



www.aracent.com



(888) 550-7750



sales@aracent.com



EASY ASSEMBLY

Put it together in minutes.



STABLE AND STURDY FRAME

Supports up to 300 pounds.



FOLDS QUICKLY AND EASILY

For easy travel and compact storage.

COMPACT & TRAVEL-READY DESIGN

This foldable rollator easily collapses to fit into most car trunks and tight storage spaces, making it perfect for travel and on-the-go use. Enjoy hassle-free portability without sacrificing stability.

SMOOTH AND SUPPORTED NAVIGATION

The core rollator is built from durable steel, with smooth gliding wheels. Confidently navigate any surface, indoors or outdoors with adjustable features and a convenient included storage bag. Perfect for those recovering from surgery or injury, or those who need added support for daily mobility.

DURABLE STEEL CONSTRUCTION

Designed to support up to 300 pounds, the durable steel frame provides enhanced support for daily use. Weighing only 19 pounds, this rollator remains portable while providing reliable stability.

TOOL-FREE ASSEMBLY & ADJUSTABLE FEATURES

Enjoy a quick and simple setup with no tools required. The adjustable design allows you to customize the rollator to your preferred height from 31 inches to 36 inches, promoting a more natural posture and stride.



(888) 550-7750



sales@aracent.com

vive
health

SR24 - 12R1