

vive health



VIVE EVERGLIDE ROLATOR MOB1057

The Vive Everglide Rollator Series-S provides stabilizing support while walking and a comfortable seat when resting. This rollator is made with a durable steel frame with height-adjustable handles, a padded backrest. The Vive Everglide Rollator Series-S folds flat for easy portability and includes an under-the-seat bin for storing everyday items.

WHAT'S INCLUDED

- The Vive Everglide Rollator
- Bag

COLORS   



A. Frame
B. Front Wheels
C. Back Wheel and Arms
D. Bolt
E. Washer

F. Lock Knob
G. Cane Holder
H. Backrest
I. Bag



www.aracent.com



(888) 550-7750



sales@aracent.com



SPECIFICATIONS

Overall Width: 24.25"

Overall Length: 27"

Armrest Lowest Height: 33"

Armrest Highest Height: 38"

Seat Width: 17.5"

Seat Depth: 10"

High Seat: 21"

Backrest adjusts from 11-16" high (from the seat)

Total weight: 18 pounds

Maximum user weight: 300 pounds



STABLE MOBILITY SUPPORT

Providing stabilizing support, the Vive rollator allows you to easily navigate indoors and out and comfortably rest when needed.

STRONG STEEL FRAME

The Vive rollator is made with a durable steel frame that requires minimal, tool-free assembly.

FOLDS FOR EASY STORAGE

Designed for easy transportation, the rollator walker folds flat for convenient storage.

LOCKING BRAKING SYSTEM

Each of the ergonomic handles is angled to position hands in a neutral position, reducing fatigue and allowing easy access to the bicycle-style hand brakes. The loop brakes also allow you to lock the wheels in place, ensuring the rollator will not move when activated.



SR22-10R2



www.aracent.com



(888) 550-7750



sales@aracent.com