

vive health

VIVE

FOLDABLE ROLLATOR SERIES T MOB1044BLU

Safely and securely navigate indoors and outdoors with the Vive Series T rollator walker. Sized to fit taller individuals up to 6' 5", the sturdy rollator features oversized wheels with a dual braking system to ensure stability on any surface. The rollator is foldable and lightweight, with removable wheels and a weight of only 18 pounds. Easily store it at home, in the trunk of your car, or anywhere while on the go. A removable storage bag is also included to keep essentials such as glasses, medications and a wallet close at hand.

WHAT'S INCLUDED

- Vive Rollator Walker Series T Frame
- 2 Removable Front Wheels
- 2 Fixed Back Wheels
- Removable Storage Bag

COLORS 

SAFELY AND SECURELY NAVIGATE INDOORS AND OUT



www.aracent.com



(888) 550-7750



sales@aracent.com



Specifications

Overall Dimensions	25" x 28"
Folded Dimensions	10" x 20" x 32"
Seat	17.5" x 9.5"
Adjustable Armrest Height	35" to 40"
Handle Height can be adjusted in 1" increments with the push pins located on the sides of the rollator	
Wheels Diameter	8"
Materials	Lightweight aluminum
Rollator Weight	19 lbs
Weight Capacity	300 lbs.
Storage Bag Dimensions	15" x 13" x 3.5"
Storage Bag Weight Capacity	11 lbs

STABLE SUPPORT FOR TALL INDIVIDUALS

Designed for taller men and women, the Vive Series T rollator walker provides stable mobility support while walking indoors or out. The tall rollator style allows you to maintain good posture while navigating any type of surface with large sport wheels. The rollator's ergonomic grip handles are height adjustable, and can be set from 35" to 40" for comfortable use by individuals 5' 11" to 6' 5" tall.

EASY TO FOLD AND STORE

Compact when folded, the collapsible rollator is designed to easily fit in the trunk of a car or even in the backseat as needed. For quick, hassle-free folding, simply pull up on the center strap located on the seat. For more compact storage, the rollator's front wheels are removable with a simple push-button release and the back wheels are hinged to fold up against the frame.

DESIGNED FOR DURABILITY AND COMFORT

Constructed with a strong, yet lightweight aluminum frame, the Vive Series T rollator safely supports up to 300 pounds. The mobility aid features dual locking handbrakes for greater security, as well as a wide fabric seat and backrest for increased comfort and better rest while on the go.

REMOVABLE STORAGE BAG

Conveniently and safely store daily necessities in the included rollator bag. The sturdy accessory bag holds a purse or wallet, medications, electronic devices and more to keep them close at hand while traveling.



SR22- 9R1



www.aracent.com



(888) 550-7750



sales@aracent.com