

# vive health

## VIVE

**UPRIGHT WALKER, SERIES T**  
**MOB1055BLU**

The Vive Mobility Upright Walker features tall adjustable handles allowing you to stand up straight, promoting proper posture. Each arm rail is height and length adjustable and features soft, contoured armrests and textured hand grips for greater control and comfort. Supporting up to 400 pounds, the foldable aluminum walker weighs approximately 25 pounds and features a wide, stable base and seat for greater stability and comfort.

### WHAT'S INCLUDED

- Upright walker just like our current model except with a wider base, larger padded seat, for users height of 5' - 6'4" and up to 400lbs

**COLORS** ■

## WALK UPRIGHT

**GREATER STABILITY  
AND COMFORT**



- |   |                         |
|---|-------------------------|
| A. 1x Base                                | E. 2x Bolts             |
| B. 2x Arm Extensions with Elbow Pads      | F. 2x Female Star Locks |
| C. 2x Front Legs                          | G. 2x Curved Washers    |
| D. 2x Back Legs Combined with Hand Brakes | H. 1x Bag               |
|   | I. 1x Cane Holder       |
|   | J. 2x Handles           |



aracent.com



(888) 550-7750



sales@aracent.com



#### Specifications

Overall Width	28"
Overall Length	30"
Armrest Lowest Height	41"
Armrest Highest Height	50"
Seat Width	17"
Seat Depth	14"
Product Weight	27 LBS
Max User Weight	400 LBS



#### WALK UPRIGHT WITH GREATER EASE

Stand tall and walk with stabilizing assistance with the Vive Mobility Series T Upright Walker. Our Series T upright walker rollator promotes a natural gait. It features a wider base for greater stability and flat-free 9.5" tires to easily navigate over most surfaces.

#### EASY HEIGHT ADJUSTMENTS

Set the armrests at the perfect height for you with simple locking push pins. The handle height is adjustable from 40" to 50" to promote proper posture and allow you to look forward and walk upright. The folding walker handles are equipped with contoured armrests to relieve fatigue and length adjustable hand grips with easy-to-use loop brakes to lock the wheels in place for stable seating.

#### WIDER, ON-THE-GO SEATING

Our Vive Mobility upright walker rollator provides comfortably padded seating everywhere you go. The stable seat measures 15" by 17" with 20" of space within the frame. Under the seat, a storage tote allows you to keep essential items close at hand.

#### LIGHTWEIGHT FOLDING FRAME:

Great for travel, our upright walker is foldable for easily transporting in car, truck or SUV. The walker is durable, safely supporting up to 400 pounds and is constructed with a lightweight aluminum frame that weighs 25 pounds.

#### VIVE GUARANTEE

One year guarantee so you can purchase now with confidence.



SR22- 8R1



[aracent.com](https://aracent.com)



(888) 550-7750



[sales@aracent.com](mailto:sales@aracent.com)