

vive health



VIVE

KNEE WALKER

MOB1007
HCPCS: E0118

The Vive knee scooter is a great alternative to crutches, providing greater independence for those recovering from an injury or surgery of the lower leg, ankle or foot, including breaks, sprains and amputations. Fitted with a contoured foam platform, the steerable scooter effectively supports the knee and lower leg to reduce painful pressure points and minimize fatigue.

Additionally, both the platform and steering handle are height adjustable for a comfortable, customized fit. The durable frame supports up to 300 pounds and is lightweight and compact, folding for easy storage and travel.

WHAT'S INCLUDED

- Knee scooter
- Storage bag

COLORS   

EASY TO MANEUVER BOTH INDOORS AND OUT



aracent.com



(888) 550-7750



sales@aracent.com



Specifications

Knee Pad Height Adjustable between 18"-23"

Handle Bar Height Adjustable in 1" increments ranging from 33" to 41"

Knee Pad Dimensions 14" x 7"
(at widest part)

Front Wheels + Frame 16" across
(widest point of scooter)

Four 8-inch polyurethane wheels provide an exceptionally stable base for easy maneuverability both indoors and outdoors.

Two fold points and the removable platform allows the scooter to be easily transported.

Weight limit 350 LBS



GREAT ALTERNATIVE TO CRUTCHES

Regain independent mobility while recovering from injury or surgery of the lower leg, ankle or foot. The knee walker is designed for safe indoor and outdoor use with 8" all-terrain tires, a dual braking system and comfort grip steering handles.

COMPACT FOLDING DESIGN

Perfect for active individuals, the Vive knee scooter features a two-point folding system and removable platform for compact storage and travel. Weighing in at 22 pounds, the knee walker is easy to lift into vehicles and move up and down stairs as necessary.

HEIGHT ADJUSTABLE

Customizable for a comfortable fit, the knee platform can be adjusted between 19"-23.5". The steering handle is also adjustable in 1" increments ranging from 33.5" to 41.5" to accommodate most individuals.

CONTOURED KNEE PLATFORM

Cradling the knee and lower leg, the high-density foam platform evenly distributes downward pressure to relieve painful pressure points and minimize fatigue in the leg, knee and hip.

SUPERIOR STABILITY AND CONTROL

Safely supporting up to 350 pounds, the steel knee walker is exceptionally durable. Creating a stable base, the steerable scooter uses large, flat-free tires and a responsive, dual-braking system for greater stability and maneuverability.

REMOVABLE STORAGE BAG

Keep essentials close at hand, on and off the scooter. Two exterior pockets allow you to quickly access electronic devices, water bottles, snacks and glasses while ensuring valuable items are securely held in one of the two interior compartments.



8955 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808

SR22- 6R1



aracent.com



(888) 550-7750



sales@aracent.com