


vive health



VIVE POWER WHEELCHAIR MOB1029L

Stylish and exceptionally durable, the Vive power wheelchair smoothly navigates most surfaces with powerful dual motors, a universal joystick controller and large, flat-free tires. The wheelchair's universal controller can be mounted on either the left or right armrest, both of which are height adjustable for greater comfort. The power wheelchair also has a travel distance of up to six miles and smoothly handles moderate inclines. Supporting up to 260 pounds, the motorized wheelchair is constructed with a durable, yet lightweight aluminum and is designed to easily fold for storage and travel.

WHAT'S INCLUDED

- Vive Mobility Power Wheelchair
- Rechargeable lithium battery
- 3-Pin charging cable

STYLISH AND EXCEPTIONALLY DURABLE



www.aracent.com



(888) 550-7750



sales@aracent.com



Specification	
Top Speed:	4 mph
Maximum Drive Range Per Charge:	9.3 mph
Weight Capacity:	265 Pounds
Total Weight with Battery included	62 lbs.
Weight without Battery:	57 lbs.
Front Tires:	7.9"
Rear Tires:	11.8"
Maximum Speed:	Up to 4.0 mph
Turning Radius:	47.2'
Length:	34"
Width:	24.5"
Height:	40"
Dimensions when folded:	35" by 25" by 14"
Seat-to-Ground Height Range:	19.7"
Seat Dimensions:	17.25" by 17.5"
Seat back height:	21.2"
Ground Clearance:	6"
Lithium Battery	

SAFE AND SMOOTH TRANSPORTATION

Easy to operate, the FDA-approved Vive power wheelchair provides comfortable, yet stylish transportation at home and on the go. The foldable power wheelchair is equipped with dual motors, flat-free tires and shock absorbers to safely and smoothly maneuver over any terrain, including grass, pavement, gravel and more.

PERFECT FOR TRAVELING

Ideal for traveling, the motorized wheelchair is constructed with a lightweight, yet durable aluminum frame that safely supports up to 265 pounds. Rechargeable, the powerful lithium battery allows the wheelchair to navigate mild inclines and travel up to six miles per charge. TSA-approved, the power wheelchair is great for airline and cruise ship travel, folding in seconds and easily fitting in standard trunk sizes.

EASY TO OPERATE

Repositionable for the left or right side, the electric power wheelchair operates with a 360-degree, joystick and speed controller. With a maximum speed of 4 mph, the power wheelchair controller allows you to easily monitor and adjust the speed of the wheelchair. The controller also includes a horn button, and a battery charge indicator.

COMFORT SEATING

Padded with dense foam, the motorized wheelchair seat is comfortably wide, measuring 17.25" wide and 17.5" deep. The cushioned seat is covered with a soft, breathable material that is easy to clean and includes an adjustable safety buckle. The armrests are also height adjustable for personalized comfort.



SR24- 3R4



www.aracent.com



(888) 550-7750



sales@aracent.com