

vive health



VIVE ELECTRIC WHEELCHAIR MODEL V MOB1054

The Vive Mobility electric wheelchair provides a smooth ride over indoor surfaces and solid outdoor surfaces such as concrete and pavement. The power chair is equipped with a powerful motor, reaching a top speed of 4mph, a dual braking system for added safety and long-lasting batteries for a range up to 15 miles per charge. Fully adjustable, the electric wheelchair features a generously sized captain's seat that also swivels for easily entering and exiting the wheelchair and repositionable armrests and footplate.

WHAT'S INCLUDED

- Vive Mobility Electric Wheelchair
- Charging cable

Body Color
Seat Color



DESIGNED FOR OPTIMAL COMFORT



aracent.com



(888) 550-7750



sales@aracent.com



COMFORTABLE MOBILITY

Our Vive Mobility electric power chair provides optimal comfort and safe, stable mobility for those with limited mobility due to injury or illness. The large electric wheelchair smoothly navigates indoors, enabling independent mobility in comfort and style.

SAFE MANEUVERABILITY

Ideal for navigating around your home, the large power wheelchair features a tight, 25" turning radius and 2.5" ground clearance to easily navigate indoors and on level outdoor surfaces such as pavement and concrete and those with a gentle incline.

EASY TO OPERATE

Featuring a 360-degree, joystick controller, our electric wheelchair is easy to operate. The intuitive joystick controller allows you to easily navigate around obstacles while monitoring and adjusting the speed of the electric wheelchair up to 4mph. The wheelchair controller is repositionable for use on the left or right side and includes a safety horn button and a battery charge indicator.

CUSTOMIZABLE SEATING

Safely supporting up to 300 pounds, the large power chair is equipped with a cushioned 18" by 20" captain's seat with a lap belt for added safety. The wide seat is height adjustable, reclines for greater comfort and swivels for easy transfers. The footplate and armrest height, width and angle are also adjustable to provide personalized seating for most adults.

Specifications

Weight capacity	300 lbs.
Top speed*	4 mph
Estimated range*	15 miles
Max. Climbing angle	6°
Turning radius	25"
Ground clearance	2.5"

Overall dimensions

Length	40"
Width	24.75"
Height	44"

Seat dimensions

Width	18 – 20"
Depth	18"
Back height	18.5"

Power

Motor	200w x 3900rpm
Controller	Dynamic shark 40a
Batteries	12v x 35ah
Battery charger	3A
Brakes	Electromagnetic
Freewheel mode	YES
Tires flat free	Flat free
Caster wheels	6" x 2"
Drive wheels	10" x 3"

Weights

Base weight 68 lbs.	68 lbs.
Battery weight (pair)	44 lbs.
Seat weight	39 lbs.
Total weight	151 lbs.