

VIVE POWER WHEELCHAIR MODEL C

MOB1094

Navigate safely and comfortably with this durable power wheelchair, designed with electromagnetic brakes and anti-tip wheels for stability on inclines of up to 6 degrees. Featuring flat-free tires and a 2.5" ground clearance, it provides smooth movement across various terrains while supporting up to 275 pounds. The ergonomic joystick allows for easy maneuvering, with a 26" turning radius for tight spaces and an automatic braking system for enhanced safety. Customize your comfort with adjustable armrests, footplate height, and a 19.5" high backrest, along with a thickly padded 17" wide seat that swivels for easy access. Weighing just 117 pounds, this portable wheelchair can be disassembled for transport, and the dual 24V batteries provide up to 8 miles of travel on a single charge, perfect for both short and long trips.

SPECIFICATIONS

What's Included:

- 1x Vive Mobility Electric Wheelchair Model C
- 1x Long-lasting battery
- 1x Quick-connect charging cable

Uses

- A power wheelchair for easy navigation on any surface.

Dimensions

- 26" x 22" by 35"

Weight

- 117 lbs

Color

- Black

Care Instructions

- Wipe down with a damp cloth and wipe dry.



Call now!

or send us an e-mail
now to find out more!



aracent.com



(888) 550-7750



sales@aracent.com



SAFE & STURDY NAVIGATION FROM A TRUSTED BRAND

Navigate safely with a durable frame, electromagnetic brakes, and anti-tip wheels for stability on inclines of up to 6 degrees. The flat-free tires and 2.5" ground clearance ensure smooth movement across various terrains. It supports up to 275 pounds.



EASY TO USE & MANEUVER

Equipped with an ergonomic joystick for effortless steering and a 26" turning radius, this wheelchair provides precise control in tight spaces, while the automatic braking system ensures safety during stops.

ADJUSTABLE COMFORT FEATURES

Enjoy customized comfort with adjustable armrest width, footplate height, and a 19.5" high backrest. The 17" wide seat—which has extra-thick padding—offers maximum support during extended use. Its seat also swivels, allowing you to more easily get on and off.



PORTABLE DESIGN WITH MODERATE WEIGHT

Weighing 117 lbs, with the heaviest part at 77 lbs, this power wheelchair can be disassembled for transport, making it suitable for occasional travel or storing when not in use.

LONG BATTERY LIFE FOR EXTENDED USE

Travel up to 8 miles on a single charge with dual 24V batteries. Ideal for both short trips and longer outings, ensuring you stay powered throughout the day.



(888) 550-7750



sales@aracent.com

vive
health

SR25- 4R1