

**NEW
ARRIVAL**

**vive
health**



VIVE MOBILITY SCOOTER SERIES C

MOB1066

The Vive Mobility Scooter Series C features an upgraded experience: Rear shock absorbers, thick padded seat, and 8" flat-free tires provide a smooth and comfortable ride. Increased ground clearance of 3.5" allows the scooter to easily pass over obstacles such as door jambs and flooring changes.

The upgraded battery gives a long range of up to 15.5 miles per charge. Adjustable speed control adjusts up to 5 miles per hour to get you where you're going fast. Larger overall dimensions provide ample legroom and a weight limit of 300lbs to accommodate most riders. The swivel-locking seat allows the rider to comfortably enter and exit from the side of the scooter.

Body Color



Seat Color



CRUISE LONG DISTANCES IN COMFORT



aracent.com



(888) 550-7750



sales@aracent.com



Specifications

Overall Dimensions (L x W x H)	43.3" x 19.6" x 36.2"
Seat Height (from ground)	23" to 26"
Seat Width	17"
Seat Depth	16"
Armrest Height (above seat)	8.2"
Seat Back Height	23.25"
Net Weight	114 lbs.
Maximum Speed	≤4.9 mph
Maximum Braking Distance	≤59"
Minimum Turning Radius	≤4.9ft
Weight Limit	300 lbs.
Maximum Travel Distance (in Theory)	15.5 mi
Static Stability	≥ 9°
Dynamic Stability	≥ 6°
Climbing Ability	≥ 9°
Transaxle	24V/250W
Battery	Lead-acid battery (24V/20AH)
Controller Maximum Output Current	45A
Controller Minimum Output Current	2A
Front Wheel Diameter	8.8" (Solid)
Rear Wheel Diameter	8.8" (Solid)

EASILY MANEUVER INDOORS AND OUT

Providing a smooth, comfortable ride, the Vive Mobility Scooter Series C is perfect for navigating in your home, at the market and on the go. With a ground clearance of 3.5", the electric scooter easily clears most door jambs and smoothly rolls over flooring changes on four 8.86" flat-free tires.

GO FARTHER ON A SINGLE CHARGE

Our Vive Mobility Scooter Series C is powered by a long-lasting, rechargeable battery pack with a range of up to 15.5 miles per charge and a top speed of 5 mph. Your mobility scooter includes a quick connect battery and charging cable.

RIDE IN COMFORT AND STYLE

Our Vive Mobility Scooter Series C is designed for comfort, safely supporting up to 300 pounds with a large platform and height-adjustable captain's seat. Comfortably padded, the wide, 16.5" seat and backrest provide gentle cushioning and support. Additionally, the scooter seat is surrounded by two padded armrests that adjust and flip upward for greater convenience when exiting or boarding.

WRAPAROUND DELTA TILLER

Equipped with a wrap around delta tiller and drive levers, the Vive Mobility Scooter Series C intuitive control panel allows you to set the maximum speed with the variable speed dial and turn on the bright LED headlight as needed when navigating in low-light areas. The tiller display also includes a color-coded battery life indicator for checking the status of the battery at a glance.



SR23-2R1



aracent.com



(888) 550-7750



sales@aracent.com