

vive health



VIVE ROLLATOR WALKER MOB1010BLK HCPCS: E0143

Oversized wheels with a dual braking system ensure stability on any surface, while easily accessible storage conveniently holds wallets, phones, drinks and more. The compact rollator folds down to easily fit in tight spaces with its hinged and removable wheels.

WHAT'S INCLUDED

- Rollator Frame
- 2 Removable Front Wheels
- 2 Fixed Back Wheels
- Removable Storage Bag

COLORS 

SAFELY AND SECURELY MANEUVER INDOORS AND OUT



aracent.com



(888) 550-7750



sales@aracent.com



Specifications

Overall Dimensions 25" x 28"

Folded Dimensions 10" x 20" x 32"

Seat 17.5" x 9.5"

Wheels Diameter 8"

Storage Bag Dimensions 15" x 13" x 3.5"

Storage Bag Weight Capacity 11 lbs

Adjustable Handle Height 31" to 36"

Handle Height can be adjusted in 1" increments with the push pins located on the sides of the rollator

Materials Lightweight Aluminum

Product Weight 19 lbs.

Weight Capacity 300 lbs.

SUPERIOR SUPPORT AND STABILITY

Lightweight and durable, the Vive rollator provides stable support, and easily maneuvers over any type of surface with large 8" sport wheels. The handle height is adjustable for a customized fit for every individual. Supporting up to 300 pounds, the lightweight aluminum frame is corrosion-resistant and folds down for easy storage and transport. Includes a wide seat and backrest for comfortably resting while on the go.

COMPACT FOLDABLE DESIGN

The collapsible rollator is compact when folded, easily fitting in the trunk of a car or even in the backseat, allowing for hassle-free folding in one easy step. Simply pull up on the center strap under the seat. For more compact storage, the front wheels are removable with a simple push button release and the back wheels are hinged to fold up against the frame.

DUAL BRAKING SYSTEM

Each of the ergonomic handles is angled to position hands in a neutral position, reducing fatigue and allowing easy access to the bicycle-style hand brakes. Wheel locking brakes ensure the rollator will not move when activated.

CONVENIENT SAFE STORAGE OPTIONS

Conveniently and safely store daily necessities in the detachable rollator bag. Safely holds a purse or wallet, medications, electronic devices and more. Also includes an integrated cane holder.



SR22- 9R1



aracent.com



(888) 550-7750



sales@aracent.com