

vive health



VIVE WHEELCHAIR ROLLATOR MOB1018 HCPCS: E0143

The versatile hybrid rollator wheelchair allows individuals to walk, ride and rest safely on any terrain, including grass, gravel, carpeting and more. Safely maneuver inclines with the ergonomic bicycle-style hand brakes and securely lock the wheels in place when seated for added stability.

Constructed with a durable, corrosion - resistant aluminum, the combination or rollator transport chair safely supports up to 300lbs.


WHAT'S INCLUDED

- Hybrid Wheelchair Rollator
- Detachable Storage Bag

COLORS



VERSATILE HYBRID ROLLATOR WHEELCHAIR



- A. Rollator Frame
- B. Brake System
- C. 2x Front Wheel Assembly
- D. 2x Foot Rest

- E. 4x Screw Knob
- F. Rollator Bag
- G. Elastic Strap loop
- H. Cane Holder Cup



aracent.com



(888) 550-7750



sales@aracent.com

Specifications

Base when open 24" by 40"

Adjustable Height 33.25" to 38.25"

Wheels 8"

Folded 8.25" by 40"

Seat dimensions 17" x 12" x 21.75"

Weight 21 pounds

Materials Lightweight aluminum

Weight limit 300 pounds



WALK, RIDE AND REST

The versatile hybrid rollator wheelchair allows individuals to walk, ride and rest safely on any terrain, including grass, gravel, carpeting and more.

COMPACT FOLDABLE DESIGN

Folding frame easily fits in the trunk of a car and in compact storage areas.

ADJUSTABLE HYBRID

Features removable footrests and adjustable handles for a comfortable fit.

DUAL BRAKING SYSTEM

Safely maneuver inclines with the ergonomic bicycle-style hand brakes and securely lock the wheels in place when seated for added stability

LIGHTWEIGHT DURABLE FRAME

Constructed with a durable, corrosion-resistant aluminum, the combination rollator transport chair safely supports up to 300 pounds.

SAFE STORAGE

Detachable shoulder bag safely stores a purse or wallet, electronic devices, medications and more.