


# MATRIX® MAC BACK SPECIFICATIONS



**H\***  
LOADED CONTOUR DEPTH

3.5"/89mm  
6"/152mm  
9"/229mm

MAC BACK SIZE	A* OVERALL BACK LENGTH	B* LOADED INNER CONTOUR WIDTH RANGE	C ADJUSTMENT RANGE TO APEX OF LATERAL	D MOUNTING SLOT LENGTH	E MOUNTING SLOT POSITION	F LATERAL SLOT LENGTH	G OVERALL LATERAL LENGTH
1618	19"/483mm	12.5-13.5"/317-343mm	5.5-11.0"/140-279mm	6.5"/165mm	3.5"/89mm	8"/203mm	11"/279mm
1620	21"/533mm	12.5-13.5"/317-343mm	5.5-11.0"/140-279mm	6.5"/165mm	3.5"/89mm	8"/203mm	11"/279mm
1622	23"/584mm	12.5-13.5"/317-343mm	5.5-12.0"/140-305mm	6.5"/165mm	4.5"/114mm	8"/203mm	11"/279mm
1716	17"/432mm	13.5-14.5"/343-368mm	5.5-8.75"/140-222mm	5.5"/140mm	3.5"/89mm	8"/203mm	11"/279mm
1718	19"/483mm	13.5-14.5"/343-368mm	5.5-11.0"/140-279mm	6.5"/165mm	3.5"/89mm	8"/203mm	11"/279mm
1720	21"/533mm	13.5-14.5"/343-368mm	5.5-11.0"/140-279mm	6.5"/165mm	3.5"/89mm	8"/203mm	11"/279mm
1722	23"/584mm	13.5-14.5"/343-368mm	5.5-12.0"/140-305mm	6.5"/165mm	4.5"/114mm	8"/203mm	11"/279mm
1816	17"/432mm	14.5-15.5"/368-394mm	5.5-8.75"/140-222mm	5.5"/140mm	3.5"/89mm	8"/203mm	11"/279mm
1818	19"/483mm	14.5-15.5"/368-394mm	5.5-11.0"/140-279mm	6.5"/165mm	3.5"/89mm	8"/203mm	11"/279mm
1820	21"/533mm	14.5-15.5"/368-394mm	5.5-11.0"/140-279mm	6.5"/165mm	3.5"/89mm	8"/203mm	11"/279mm
1822	23"/584mm	14.5-15.5"/368-394mm	5.5-12.0"/140-305mm	6.5"/165mm	4.5"/114mm	8"/203mm	11"/279mm
1916	17"/432mm	15.5-16.5"/394-419mm	5.5-8.75"/140-222mm	5.5"/140mm	3.5"/89mm	8"/203mm	11"/279mm
1918	19"/483mm	15.5-16.5"/394-419mm	5.5-11.0"/140-279mm	6.5"/165mm	3.5"/89mm	8"/203mm	11"/279mm
1920	21"/533mm	15.5-16.5"/394-419mm	5.5-11.0"/140-279mm	6.5"/165mm	3.5"/89mm	8"/203mm	11"/279mm
1922	23"/584mm	15.5-16.5"/394-419mm	5.5-12.0"/140-305mm	6.5"/165mm	4.5"/114mm	8"/203mm	11"/279mm
2016	17"/432mm	16.5-17.5"/419-445mm	5.5-8.75"/140-222mm	5.5"/140mm	3.5"/89mm	8"/203mm	11"/279mm
2018	19"/483mm	16.5-17.5"/419-445mm	5.5-11.0"/140-279mm	6.5"/165mm	3.5"/89mm	8"/203mm	11"/279mm
2020	21"/533mm	16.5-17.5"/419-445mm	5.5-11.0"/140-279mm	6.5"/165mm	3.5"/89mm	8"/203mm	11"/279mm
2022	23"/584mm	16.5-17.5"/419-445mm	5.5-12.0"/140-305mm	6.5"/165mm	4.5"/114mm	8"/203mm	11"/279mm

\*Please note measurements for soft components are approximate.



**Motion Concepts**  
www.motionconcepts.com

**USA**  
700 Ensminger Rd. Suite 112  
Tonawanda, NY 14150  
(888) 433-6818  
Fax: (888) 433-6834

**Canada**  
84 Citation Dr.,  
Unit 1-12  
Concord, ON L4K 3C1  
(866) 748-7943  
(905) 695-0134  
Fax: (905) 695-0138

©2024 Motion Concepts.  
All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Motion Concepts unless otherwise noted. Specifications are subject to change without notification.  
Form No. TRD0741 Rev A