





INVACARE® MATRX® PB DEEP BACK SPECIFICATIONS



	BACK SHELL DIMENSIONS								OVERALL DIMENSIONS (CONTACT AREA)									
	A		B		C		D		E		F		G		WEIGHT WITHOUT HARDWARE		WEIGHT CAPACITY	
	WIDTH		HEIGHT			LATERAL CONTOUR HEIGHT		HEIGHT	INNER MAXIMUM WIDTH		OUTER MAXIMUM WIDTH		HEIGHT		lb	kg	lb	kg
MODEL	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm				
MPD1412	12	31	11	28	7½	19	1½	4	12	31	14	36	12	31	3.2	1.5	300	136
MPD1416	12	31	15	38	9	23	2½	6.5	12	31	14	36	16	41	4.0	1.8	300	136
MPD1612	14	36	11	28	7½	19	1½	4	14	36	16	41	13	33	3.6	1.6	300	136
MPD1616	14	36	16	41	9	23	2½	6.5	14	36	16	41	17½	44.5	4.6	2.1	300	136
MPD1620	14	36	20	51	12	31	3	7.5	14	36	16	41	21	53	5.2	2.4	300	136
MPD1712	15	38	11	28	7½	19	1½	4	15	38	17	43	13	33	3.9	1.8	300	136
MPD1716	15	38	16	41	9	23	2½	6.5	15	38	17	43	17½	44.5	4.9	2.2	300	136
MPD1720	15	38	20	51	12	31	3	7.5	15	38	17	43	21	53	5.6	2.5	300	136
MPD1812	17	43	11	28	7½	19	1½	4	16	41	18	46	13	33	4.2	1.9	300	136
MPD1816	17	43	16	41	9	23	2½	6.5	16	41	18	46	17½	44.5	5.2	2.4	300	136
MPD1820	17	43	20	51	12	31	3	7.5	16	41	18	46	21	53	6.0	2.7	300	136
MPD2012	19	48	11	28	7½	19	1½	4	18	46	20	51	13	33	4.8	2.2	300	136
MPD2016	19	48	16	41	9	23	2½	6.5	18	46	20	51	17½	44.5	5.8	2.6	300	136
MPD2020	19	48	20	51	12	31	3	7.5	18	46	20	51	21	53	6.6	3.0	300	136

Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd.
E. Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Motion Concepts
www.motionconcepts.com

USA
700 Ensminger Rd.,
Suite 112
Tonawanda, NY 14150
(888) 433-6818
Fax: (888) 433-6834

Canada
84 Citation Dr., Unit 1-7
Concord, ON L4K 3C1
(866) 748-7943
(905)-695-0134
Fax: 905-695-0138

© 2014 Invacare Corporation. All rights reserved.
Trademarks are identified by the symbols™ and ©.
All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.
Form 12-568 Rev. 06/14 140438