



Yes, you can.®

# Account Information

Request Type: Quote                      Order

Date: \_\_\_\_\_

Account #: \_\_\_\_\_

Company: \_\_\_\_\_

SHIP TO  
Name: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_

Country: \_\_\_\_\_

Purchase Order #: \_\_\_\_\_

CONTACT  
Name: \_\_\_\_\_

Back Up Contact: \_\_\_\_\_

Phone: \_\_\_\_\_

Email: \_\_\_\_\_

Comments: \_\_\_\_\_

## Special Note

The IVC 9000 Topaz with a 700 pound weight limit.

### FRAME TYPE

9TPZ 9000 Topaz ..... \$2,979.00 K0007

### CUSTOM PACKAGES

HDPG Heavy Duty Package ..... \$321.30  
Reinforced Frame Width weight capacity of 1,000 lbs. Only available in 28" and 30" Frame Widths.

### FRAME WIDTH

W20 20" Seat Width ..... No Charge E2201  
This selection is not available with HDPG option.

W22 22" Seat Width ..... No Charge E2201  
This selection is not available with HDPG option.

W24 24" Seat Width ..... No Charge E2202  
This selection is not available with HDPG option.

W26 26" Seat Width ..... No Charge E2202  
This selection is not available with HDPG option.

W28 28" Seat Width ..... \$267.75 K0108

W30 30" Seat Width ..... \$267.75 K0108

### FRAME DEPTH

D18 18" Seat Depth ..... No Charge

D20 20" Seat Depth ..... No Charge E2203

### SEAT-TO-FLOOR HEIGHT

ADULT Adult ..... No Charge  
Based on 24" wheels and 8" casters.

H Hemi ..... No Charge  
Based on 24" wheels and 8" casters.

9TSH Super Hemi 15 1/2" ..... No Charge  
Not available with HDPG option.  
Only available with 6" casters and 20" rear wheels.



Yes, you can.®

**ARM TYPE**

Fixed height arms 9 1/2" adjustable height arms.

T18 Full Length Fxd Ht Conv .....No Charge

T19 Full Length Adj Ht Conv .....\$134.95 E0973

T28 Desk Length Fxd Ht Conv .....No Charge

T29 Desk Length Adj Ht Conv .....\$134.95 E0973

**BACK TYPE**

ADJH Adj Height Back 16", 17", 18" with 10° bend .....

.....Standard

Adjustable back height 16"-18". Factory set at 16".

**FRAME FINISH**

71P Silver Vein .....No Charge

24P Wet Black .....No Charge

**CASTER OPTIONS**

1259 6" x 2" Semi Pneu Casters ..... \$103.89

Only available with super hemi seat-to-floor,NA with HDPG option.

1228HD 8" x 2" Semi Pneu Casters Black .....No Charge

Only available with HDPG option.

1228BK 8" x 1 3/4" Semi Pneu Casters Black .....

.....No Charge

Only available with adult and hemi seat-to-floor, NA with HDPG option.

**WHEEL OPTIONS**

5763BK 20" HD Pneu with Flat Free Inserts .....

.....\$215.27 E2211+E2213

Not available with HDPG option.

Only available with super hemi seat-to-floor.

U5222SP Spoked HD Urethane Treaded Tires .....

.....No Charge

**HANDRIM OPTIONS**

AHDRM Aluminum Handrims .....Standard

CHDRM Chrome Handrims .....No Charge

**UPHOLSTERY COLORS**

U240 Black Nylon Upholstery .....Standard

U67 Black Vinyl Upholstery .....No Charge

U550 Midnight Blue Vinyl Upholstery .....No Charge

**FRONT RIGGINGS**

TZ33A Sw Footrests with Extra Large Footplates .....

.....\$205.52 K0041

TZ34A Elevating Legrests with Extra Large Footplates ...

.....\$426.36 K0195 + K0041



Yes, you can.®

**FOOTREST/LEGREST ACCESSORIES**

1337 Calf Strap ..... \$38.37 K0038

1346 Heel Loops ..... \$41.57 E0951  
Standard on swingaway footrests.

**SEAT POSITIONING STRAPS**

T1515 Topaz Push Button Style Seat Positioning Strap ....  
.....\$53.55 E0978

**BACK ACCESSORIES**

T1483 Push Bar .....\$166.01  
This accessory is shipped in a separate carton.

**SEAT ACCESSORIES**

T9120 Solid Seat Insert .....\$160.65  
This is shipped in a separate carton.  
Solid Seat Insert is manufactured in a vinyl upholstery only.

**ACCESSORIES**

The below selections are shipped in a separate carton.

T1496 IV/Oxygen Holder Package Assembly .....  
..... \$169.22 K0105 + E2208

T1322 IV Rod and Bracket, Telescopic ..... \$158.51 K0105  
This selection is not available with super hemi seat-to-floor.

T1497 Oxygen Holder Package Assembly .....  
..... \$161.72 E2208

**MISCELLANEOUS ACCESSORIES**

T1358 Rear Anti-Tippers ..... \$92.11 E0971

1305L Crutch and Cane Carrier ..... \$46.05 E2207

A1567 Pull to Lock Clamp on Whl Locks (Std on Hemi) ...  
..... No Charge

A1576 Wheel Lock Extension ..... \$41.77 E0961

1721 Hi Push to Lock Clamp on Whl Locks (Std on Adult)  
..... No Charge  
Only available with HDPG option.

1335XL Non Folding Device - Topaz ..... \$53.55