



Yes, you can.®

# Account Information

Request Type: Quote                      Order

Date: \_\_\_\_\_

Account #: \_\_\_\_\_

Company: \_\_\_\_\_

SHIP TO  
Name: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_

Country: \_\_\_\_\_

Purchase Order #: \_\_\_\_\_

CONTACT  
Name: \_\_\_\_\_

Back Up Contact: \_\_\_\_\_

Phone: \_\_\_\_\_

Email: \_\_\_\_\_

Comments: \_\_\_\_\_

## Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

Existing information (serial number) will need to be supplied to Customer Service to ensure correct configuration.

## PRE-CONFIGURED STOCK WHEELCHAIRS

### STOCK CHAIRS

9XT\_PTO\_29154 9000XT 16" x 16" Adult, Desk Length  
Fixed Ht Space Saver Arm .....\$1,355.19 K0004

9XT\_PTO\_29153 9000XT 18" x 16" Adult, Desk Length  
Fixed Ht Space Saver Arm .....\$1,355.19 K0004

9XT\_PTO\_29155 9000XT 18" x 16" Adult, Full Length  
Fixed Ht Space Saver Arm .....\$1,355.19 K0004

9XT\_PTO\_33131 9000XT 20" x 16" Adult, Desk Length  
Fixed Ht Space Saver Arm .....\$1,355.19 K0004+E2201

### STOCK FRONT RIGGINGS

T93HAP Hemi Sw. Footrests Aluminum Footplates  
w/ Heel Loops ..... \$99.13 E0951

T93HCP Hemi Sw. Footrests Composite Footplates  
w/ Heel Loops ..... \$99.13 E0951

T94HAP Hemi Sw. Elev Legrests Aluminum Footplates  
& Padded Calf Pad ..... \$191.86 K0195

T94HCP Hemi Sw. Elev Legrests Composite Footplates  
& Padded Calf Pad ..... \$191.86 K0195

AHL4A/U67 Hemi Smart Leg Articulating Legrests  
Aluminum Footplates ..... \$341.09 K0053

### STOCK OPTIONS

1497 Oxygen Holder Pkg Assy ..... \$161.72 E2208

1496 IV / Oxygen Hldr Pkg Assy ..... \$169.22 K0105+E2208

1360 Adjustable Rear Anti-Tippers- 8" Caster .....  
..... \$87.82 E0971

1360A Adjustable Rear Anti-Tippers- 6" Caster.....  
..... \$87.82 E0971

A1576 Wheel Lock Extensions ..... \$41.77 E0961



Yes, you can.®

# FRAME AND SEAT SELECTIONS

## FRAME TYPE

9XT 9000XT .....\$1,702.00 K0004

## FRAME WIDTH AND DEPTH

WD46 14" x 16" ..... No Charge E1236  
14" Seat width changes base code to E1236

WD48 14" x 18" ..... \$126.38 E1236  
14" Seat width changes base code to E1236

WD56 15" x 16" ..... \$126.38

WD58 15" x 18" ..... \$126.38

WD66 16" x 16" ..... No Charge

WD67 16" x 17" ..... \$126.38

WD68 16" x 18" ..... \$126.38

WD76 17" x 16" ..... \$126.38

WD77 17" x 17" ..... \$126.38

WD78 17" x 18" ..... \$126.38

WD86 18" x 16" ..... No Charge

WD87 18" x 17" ..... \$126.38

WD88 18" x 18" ..... \$126.38

WD96 19" x 16" ..... \$152.08  
Not available with Super Low seat-to-floor.

WD97 19" x 17" ..... \$272.03  
Not available with Super Low seat-to-floor.

WD98 19" x 18" ..... \$272.03  
Not available with Super Low seat-to-floor.

WD06 20" x 16" ..... \$152.08 E2201  
Not available with Super Low seat-to-floor.

WD07 20" x 17" ..... \$272.03 E2201  
Not available with Super Low seat-to-floor.

WD08 20" x 18" ..... \$272.03 E2201  
Not available with Super Low seat-to-floor.

WD26 22" x 16" ..... \$152.08 E2201  
Not available with Super Low seat-to-floor.

WD27 22" x 17" ..... \$272.03 E2201  
Not available with Super Low seat-to-floor.

WD28 22" x 18" ..... \$272.03 E2201  
Not available with Super Low seat-to-floor.

## SEAT-TO-FLOOR HEIGHT

A Adult (19 1/2" STF) ..... No Charge  
Based on 24" wheels and 8" casters. See seat-to-floor matrix.

H Hemi (17 1/2" STF) ..... No Charge  
Based on 24" wheels and 8" casters. See seat-to-floor matrix.

SH Super Low (14 1/2" STF) ..... \$91.04 K0056

SHE Super Hemi (15 1/2" STF) ..... No Charge K0056  
Based on 20" wheels and 6" casters. See seat-to-floor matrix.  
Requires AT902 (6" casters), U1763BK (20" composite wheels, urethane tires) AHDRM (aluminum handrims) and A1576 (wheel lock extensions).  
Seat-to-floor height is approximately 15 1/2".

# CASTERS

## CASTER OPTIONS

AT902 6" Composite Urethane 3 Spoke ..... No Charge

1228 8" x 1 3/4" Semi. Pneu Casters Lt Grey ..... \$59.98  
Not available with Super Low seat-to-floor.

AT903 8" x 1 1/4" Composite Urethane ..... Standard  
Not available with Super Low seat-to-floor.

1299 8" x 1" Solid Rubber Casters ..... No Charge  
Not available with Super Low seat-to-floor.

1292 8" x 1 1/4" Pneu Casters ..... \$93.18 E2214+E2215  
Not available with Super Low seat-to-floor.

1492 8" x 1 1/4" Pneu. with Flat Free Inserts ... \$111.38 E2217  
Not available with Super Low seat-to-floor.

## CASTER JOURNAL OPTION

Not available with Adjustable Angle or Adjustable Angle-S Back Type.  
Not available with Super Low seat-to-floor.  
Fixed caster journals are standard.

ADJC Adjustable Caster Journal ..... \$31.06

ADJCPB Adjustable Precision Bearing Caster Journal ..... \$43.91



# ARM AND BACK OPTIONS

## ARM TYPE

Arm types are removable, flip-back.

- 18 Full Length Fixed Ht Conv (10") .....No Charge
- 19 Full Length Adj. Ht. Conv (10" - 14") ..... \$134.95 E0973
- 28 Desk Length Fixed Ht. Conv (10") .....No Charge
- 29 Desk Length Adj. Ht. Conv (10" - 14") .... \$134.95 E0973
- 78 Full Length Fixed Ht. SS (10") .....No Charge  
Not available with adjustable angle back type.
- 79 Full Length Adj. Ht. SS (10" - 14") ..... \$134.95 E0973  
Not available with adjustable angle back type.

- 88 Desk Length Fixed Ht. SS (10") .....No Charge  
Not available with adjustable angle back type.
- 89 Desk Length Adj. Ht. SS (10" - 14") ..... \$134.95 E0973  
Not available with adjustable angle back type.
- 79JY Full Length Adj. Ht. SS (8" - 10") .....\$139.23 E0973  
Not available with adjustable angle back type.
- 89JY Desk Length Adj. Ht. SS (8" - 10") .....\$139.23 E0973  
Not available with adjustable angle back type.

## BACK TYPE

- ADJ Adj. Back with 10° Bend .....No Charge
- S1418 Straight Back w/ Push Handles (15" - 18") .....  
.....No Charge

- ANGLE Adj. Angle Back with 10° Bend .....\$107.10 K0108  
Not available with Super Low seat-to-floor. Standard with rear anti-tippers.
- ANGLE-S Adj. Angle Back with Straight .....\$107.10 K0108  
Not available with Super Low seat-to-floor. Standard with rear anti-tippers.

# COLOR OPTIONS

## FRAME FINISH

- 08P Deep Purple .....No Charge
- 17P Cobalt Blue .....No Charge
- 24P Wet Black .....No Charge
- 34P Satin (Flat) Black .....No Charge

- 71P Silver Vein .....No Charge
- 88P Wet Black with Twilight Sparkle .....No Charge
- 105P Midnight Blue .....No Charge

## CANDY COLORS

- 61P Electric Red .....No Charge
- 62P Electric Blue .....No Charge

- 66P Electric Purple .....No Charge
- 104P Electric Teal .....No Charge

## UPHOLSTERY COLORS

- U67 Black Vinyl .....No Charge  
Vinyl upholstery must choose WRAP or SOUPH.
- U240 Black Nylon Upholstery .....Standard

- U245 Black Sail Cloth Adj. Tension Seat ..... \$10.71  
U245 sail cloth adjustable tension seat upholstery available with U240 black nylon back upholstery only.

## NAUGAHYDE

- WRAP Wrap Around Upholstery .....\$82.47

- SOUPH Screwed On Back Upholstery .....\$95.32  
Not available with Adjustable Angle or Adjustable Angle-S Back Type.



# WHEELS

## WHEEL OPTIONS

- 1763BK 20" Composite Pneu. Tires .....\$125.31 E2211 + E2212
- 4763BK 20" Composite Pneu. Flat Free Inserts .....  
.....\$248.47 E2213 + E2211
- U1763BK 20" Composite Urethane Tires .....\$121.02
- 1764BK 22" Composite Pneu. Tires .....\$125.31 E2211 + E2212
- 4764BK 22" Composite Pneu Flat Free Inserts .....  
.....\$248.47 E2211+E2213
- 5244BK 24" HD Comp. Pneu. with FF Inserts .....  
.....\$319.16 E2211+E2213  
Not available with Super Low seat-to-floor.
- U1764BK 22" Composite Urethane Tires .....\$121.02

- 2244BK 24" Composite Pneu. Tires ...\$149.94 E2211 + E2212  
Not available with Super Low seat-to-floor.
- 4244BK 24" Composite Pneu. with FF Inserts .....  
.....\$240.98 E2211+E2213  
Not available with Super Low seat-to-floor.
- U2244BK 24" Composite Urethane Treaded Tires .....  
.....\$128.52  
Not available with Super Low seat-to-floor.
- U2222C 24" Composite Urethane Corded Tires .....  
.....Standard  
Not available with Super Low seat-to-floor.

## WHEEL OPTIONS (ONE ARM DRIVE - OAD)

The Cyclical Lever drive requires treaded rear wheels and is only available on 16", 17", and 18" deep. Hemi seat-to-floor must use space saver arms. Weight limit is 250 lbs. The above options are not available with amputee attachment or offset axle.

- 9456L LH One Arm Drive with U2222C .....\$649.03 E0958
- 9456R RH One Arm Drive with U2222C .....\$649.03 E0958
- LDR96R1 6" Right Hand Cyclical Lever One Arm Drive ...  
.....\$853.59 E0958  
Not available with Super Low seat-to-floor.  
Not available with adjustable caster journal.
- LDR96L1 6" Left Hand Cyclical Lever One Arm Drive .....  
.....\$853.59 E0958  
Not available with Super Low seat-to-floor.  
Not available with adjustable caster journal.

- LDR98R1 8" Right Hand Cyclical Lever One Arm Drive ....  
.....\$853.59 E0958  
Not available with Super Low seat-to-floor.  
Not available with adjustable caster journal.
- LDR98L1 8" Left Hand Cyclical Lever One Arm Drive .....  
.....\$853.59 E0958  
Not available with Super Low seat-to-floor.  
Not available with adjustable caster journal.

## HANDRIM OPTIONS

- AHDRM Aluminum Handrims .....No Charge  
Not available with U2222C rear wheels.
- COM Composite Handrims .....Standard  
Not available on 20" and 22" composite rear wheels.
- 1512 Plstc Coated Aluminum Handrims .....\$91.04

- 1224-8 8 Alum Oblique Proj. Hdrm, Plstc Coat ...\$269.89  
Not available on 20" and 22" composite rear wheels.
- 1225-8 8 Alum Vert. Proj. Hdrm, Plstc Coat .....\$269.89  
Not available on 20" and 22" composite rear wheels.

## REAR AXLE

- OFFAXL Offset Axle .....\$42.84 K0108

Not available with Adjustable Angle or Adjustable Angle-S Back Type.  
Not available with Super Low seat-to-floor.



### WHEEL LOCKS

The above selections are not available with Adjustable Angle or Adjustable Angle-S Back Type.

- 9690B Pull To Lock Whl Locks (Std on Hemi) .....  
.....No Charge
- 9691B Push To Lock Whl Locks (Std On Adult) .....  
.....No Charge

- A1567 Pull To lock Clamp On Whl Locks ..... \$81.40
- 1721 High Push To Lock Clamp On Whl Locks..... \$76.04

## POSITIONING STRAP OPTIONS

### SEAT POSITIONING STRAPS

- 1515 Push Button Style Seat Positioning Strap .....  
.....\$53.55 E0978

- 1321BK Hook and Loop Seat Positioning Strap .....  
.....\$42.84 E0978

### CHEST POSITIONING STRAPS

- 9515 Push Button Style Chest Positioning Strap .....  
.....\$53.55 E0960

## ACCESSORIES AND STOCK OPTIONS

### BACK ACCESSORIES

- L8601 Omit Back Upholstery .....Less \$10.00
- ADJTENS Adjustable Tension Back Upholstery .....  
.....\$119.95
- 1160 Carrying Pocket ..... \$57.83

- STC Stencil Back Upholstery..... \$31.06  
Only available with vinyl upholstery. Block letters are .75" in height.
- RB Dual Back Liner .....\$52.48
- BPO2 Back Pack .....\$53.55

### SEAT ACCESSORIES

- L8651 Omit Seat Upholstery .....Less \$10.00
- 9120 Solid Seat Insert ..... \$104.96  
Not available on 22" wide chairs.

- RS Dual Seat Liner ..... \$52.48

### ARM ACCESSORIES

- NRAS Non-Removable Arm With Screws ..... \$10.71

- 1300 Utility Tray ..... \$168.15 E0950  
Not available with 24" rear wheel set in Hemi position.

### MISCELLANEOUS ACCESSORIES

- 1305L Crutch and Cane Carrier .....\$46.05 E2207  
Not available with traditional one arm drive, but is available with CLD one arm drive.
- 1577QR Quick Release Axles .....\$92.11  
Not available with offset axle.  
Not available with any one arm drive.
- 2218 IV Rod, Telescopic .....\$160.65 K0105  
Not available with traditional one arm drive, but is available with CLD one arm drive.

- 2226 Front Anti - Tippers .....\$88.89 E0971
- 9357SH Super Low Amputee Attachment ..\$101.75 E0959  
Only available on Super Low Seat-to-Floor.
- 1360SH Super Low Adj Rear Anti - Tipper .....\$87.82 E0971  
Only available on Super Low Seat-to-Floor.
- 9357L Amputee Attachment ..... \$100.78 E0959  
Not available with offset axle.  
Not available with any one arm drive.



### STOCK OPTIONS

1497 Oxygen Holder Package Assembly ..... \$161.72 E2208  
Not available with traditional one arm drive, but is available with CLD one arm drive.

Not available with offset axle.  
For 0° and 6° seat angles. See service manual for details.

1496 IV / Oxygen Holder Package Assembly .....  
..... \$169.22 K0105+E2208

Not available with offset axle.  
Not available with any one arm drive.

1360 Adjustable Rear Anti-Tippers- 8" Caster .....  
..... \$87.82 E0971

Only available with 8" casters.  
Standard with Adjustable Angle Back.

1360A Adjustable Rear Anti - Tipper (6" Casters) .....  
..... \$87.82 E0971

Only available with 6" casters.

A1576 Wheel Lock Extensions ..... \$41.77 E0961  
Elastic cord wheel lock extensions. Hemi seat to floor come standard with non-cord wheel lock extensions.

1100 Tool Kit ..... No Charge

## FRONT RIGGING OPTIONS

### FRONT RIGGING

T93HA Hemi Sw. Footrests Aluminum Footplates  
w/ Heel Loops ..... \$99.13 E0951

T93HC Hemi Sw. Footrests Composite Footplates  
w/ Heel Loops ..... \$99.13 E0951

T94HA Hemi Sw. Elev Legrests Aluminum Footplates  
& Padded Calf Pad ..... \$191.86 K0195

T94HC Hemi Sw. Elev Legrests Composite Footplates  
& Padded Calf Pad ..... \$191.86 K0195

SH904M Super Low Sw. Elevating Legrests 11"-17" .....  
..... \$283.53 K0195

SH904MS Super Low Sw. Elevating Leg rests 14"W .....  
..... \$283.53 K0195

SH93E Super Low Sw. Ext. Footrests 11" - 17" ..... \$197.19

SH93ES Super Low Sw. Ext. Footrests 14"W ..... \$197.19

AHL4A Hemi Smart Leg Articulating Legrests  
Aluminum Footplates ..... \$341.09 K0053

AHL4C Hemi Smart Leg Articulating Legrests  
Composite Footplates ..... \$341.09 K0053

AHS4A Hemi Smart Leg Articulating Legrestst-14-15"  
W Chairs ..... \$341.09 K0053

T93HS Hemi Sw. Footrests Aluminum Footplates  
w/ Heel Loops- 14"W Chairs ..... \$95.93

T94HS Hemi Sw. Elevating Legrests Aluminum  
Footplates & Padded Calf Pad- 14"W Chairs .....  
..... \$186.53 K0195

T963A 60° Sw. Footrests Aluminum Footplates ... \$95.93

T963C 60° Sw. Footrests Composite Footplates .. \$95.93

T963S 60° Sw. Footrests Composite Footplate 14"W .....  
..... \$95.93 K0041

### FOOTPLATE OPTIONS

1350 Extra Large Aluminum Footplates ..... \$84.21 K0041

AT5543 Adjustable Angle Flip - up Footplates .....  
..... \$71.42 K0040

### FOOTREST / LEGREST ACCESSORIES

1820 Non - Removable Device ..... \$12.79

1337 Calf Strap ..... \$38.37 K0038

1346 Heel Loops ..... \$41.57 E0951

IMPACTG Impact Guards ..... \$54.36




Yes, you can.®

**Seat-to-Floor 6" Casters**

SEAT-TO-FLOOR HEIGHT	FRONT CASTER MOUNTING POSITION	REAR WHEEL SIZE	REAR WHEEL MOUNTING POSITION
15 1/2"	Top	20"	Top
16 1/2"	Middle	22"	Top
17 1/2"	Bottom	24"	Top
17 1/2"	Bottom	20"	Bottom

**Seat-to-Floor 8" Casters**

SEAT-TO-FLOOR HEIGHT	FRONT CASTER MOUNTING POSITION	REAR WHEEL SIZE	REAR WHEEL MOUNTING POSITION
17 1/2"	Top	20"	Bottom
17 1/2"	Top	24"	Top
18 1/2"	Middle	22"	Bottom
19 1/2"	Bottom	24"	Bottom



**Weight Limit 250 lbs**

Invacare recommends all wheelchairs be fitted to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.

**Seat-to-Floor 6" Casters and Offset Axle Bracket**

SEAT-TO-FLOOR HEIGHT	FRONT CASTER MOUNTING POSITION	REAR WHEEL SIZE	REAR WHEEL MOUNTING POSITION	AXLE PLATE MOUNTING ON FRAME	
				Top Lug	Bottom Lug
15 1/2"	Top Hole	22"	#5 Hole	#4 Hole	#1 Hole
15 1/2"	Top Hole	20"	#3 Hole	#4 Hole	#1 Hole
16 1/2"	Middle Hole	24"	#5 Hole	#4 Hole	#1 Hole
16 1/2"	Middle Hole	22"	#3 Hole	#4 Hole	#1 Hole
16 1/2"	Middle Hole	20"	#2 Hole	#4 Hole	#1 Hole
17 1/2"	Bottom Hole	24"	#3 Hole	#4 Hole	#1 Hole
17 1/2"	Bottom Hole	22"	#2 Hole	#4 Hole	#1 Hole
17 1/2"	Bottom Hole	20"	#3 Hole	#4 Hole	#1 Hole

For a complete list of seat-to-floor heights and angles that are configurable by the dealer refer to the service manual.

**Seat-to-Floor 8" Casters and Offset Axle Bracket**

SEAT-TO-FLOOR HEIGHT	FRONT CASTER MOUNTING POSITION	REAR WHEEL SIZE	REAR WHEEL MOUNTING POSITION	AXLE PLATE MOUNTING ON FRAME	
				Top Lug	Bottom Lug
17 1/2"	Top Hole	20"	#3 Hole	#1 Hole	#4 Hole
17 1/2"	Top Hole	22"	#2 Hole	#4 Hole	#1 Hole
17 1/2"	Top Hole	24"	#3 Hole	#4 Hole	#1 Hole
18 1/2"	Middle Hole	22"	#3 Hole	#1 Hole	#4 Hole
18 1/2"	Middle Hole	24"	#2 Hole	#4 Hole	#1 Hole
19 1/2"	Bottom Hole	20"	#3 Hole	#1 Hole	#4 Hole

For a complete list of seat-to-floor heights and angles that are configurable by the dealer refer to the service manual.