



Yes, you can.®

Account Information

Request Type: Quote Order

Date: _____

Account #: _____

Company: _____

SHIP TO
Name: _____

Address: _____

City: _____ State: _____ Zip: _____

Country: _____

Purchase Order #: _____

CONTACT
Name: _____

Back Up Contact: _____

Phone: _____

Email: _____

Comments: _____

Special Note

The IVC 9000XT Recliner is standard with attendant operated wheel locks

FRAME AND SEAT SELECTIONS

FRAME TYPE

9RC 9000 XT Recliner \$2,145.00 K0004 + E1226

FRAME WIDTH AND DEPTH

WD46 14" x 16" No Charge E1236
14" Seat width changes base code to E1236

WD48 14" x 18" \$126.38 E1236
14" Seat width changes base code to E1236

WD66 16" x 16" No Charge

WD67 16" x 17" No Charge

WD68 16" x 18" \$126.38

WD86 18" x 16" No Charge

WD87 18" x 17" No Charge

WD88 18" x 18" \$126.38

WD06 20" x 16" \$152.08 E2201

WD07 20" x 17" \$272.03

WD08 20" x 18" \$272.03 E2201

WD26 22" x 16" No Charge E2201

WD27 22" x 17" \$272.03 E2201

WD28 22" x 18" \$272.03 E2201



Yes, you can.®

SEAT-TO-FLOOR HEIGHT

R ReclinerNo Charge

ARM AND BACK OPTIONS

ARM TYPE

Fixed Height Arms 9", Adjustable Height Arms 9"-13", removable style.

18 Full Length Fxd Ht Conv No Charge

19 Full Length Adj Ht Conv\$134.95 E0973

28 Desk Length Fxd Ht ConvNo Charge

29 Desk Length Adj Ht Conv\$134.95 E0973

19JY Full Length Adj Ht Conv (8" - 10") \$139.23 E0973

29JY Desk Length Adj Ht Conv (8" - 10") ... \$139.23 E0973

BACK TYPE

Add an additional 11" for the headrest extension.

2R Reclining Back 23" High No Charge

Lift package for Reclining Back 23" High ... \$447.68 K0108

COLOR OPTIONS

FRAME FINISH

17P Cobalt Blue No Charge

61P Electric Red No Charge

24P Wet BlackNo Charge

34P Satin (flat) Black Standard

71P Silver Vein No Charge

UPHOLSTERY COLORS

All back upholstery is screwed on.

U240 Black Nylon Upholstery Standard

U245 Black Sail Cloth Adj Tension Seat \$10.71
U245 sail cloth adjustable tension seat upholstery available with U240 black nylon upholstery only.

U67 Black Vinyl No Charge

U550 Midnight Blue Vinyl No Charge

U86 Ruby Red Vinyl No Charge

CASTERS

CASTER OPTIONS

1228 8" x 1.75" Semi Pneu Casters Light Grey \$59.98

AT903 8" x 1.25" Comp Urethane Standard

1290 8" x 1" Solid Rubber Casters CharcoalNo Charge

1292 8" x 1.25" Pneu Casters\$93.18 E2214+E2215

1492 8" x 1.25" Pneu with Flat Free Inserts
..... \$111.38 E2217

1299 8" x 1" Solid Rubber Casters Dark Grey Standard



Yes, you can.®

WHEELS

WHEEL OPTIONS

- 1763BK 20" Composite Pneu Tires\$125.31 E2211 + E2212
- 4763BK 20" Composite Pneu. Flat Free Inserts
.....\$248.47 E2211 + E2213
- U1763BK 20" Composite Urethane Tires \$121.02
- 1764BK 22" Composite Pneu Tires
.....\$125.31 E2211 + E2212
- 4764BK 22" Composite Pneu. Flat Free Inserts
.....\$248.47 E2211 + E2213

- U1764BK 22" Composite Urethane Tires
..... \$121.02
- 2244BK 24" Composite Pneumatic Tires
..... \$149.94 E2211 + E2212
- 4244BK 24" Composite Pneu. Flat Free Inserts
.....\$240.98 E2211 + E2213
- U2244BK 24" Comp Urethane Treaded Tires\$128.52
- U2222C 24" Comp Urethane Corded Tires Standard

HANDRIM OPTIONS

- AHDRM Aluminum HandrimsNo Charge
Not available with U2222C rear wheels.
- COM Composite Handrims Standard
Available with 24" rear wheels only.
- 1512 Plastic Coated Aluminum Handrims \$91.04

- 1224-8 8 Alum Oblique Proj Handrims, Plst Coated
.....\$269.89
Available with 24" rear wheels only.
- 1225-8 8 Alum Vert Proj Hndrms, Plst Coated\$269.89
Available with 24" rear wheels only.

WHEEL LOCKS

- Attendant Operated Standard
Attendant operated wheel locks are standard and not available with patient operated wheel locks.
- 9696R Patient Operated Push To Lock (Recliner)
.....\$96.39

- 1329BK Toggle Wheel Lock Extension (8")\$41.77 E0961
- 1245 Patient Operated Pull to Lock (Recliner)
..... No Charge

POSITIONING STRAP OPTIONS

SEAT POSITIONING STRAPS

- 1515 Push Button Style Seat Positioning Strap
..... \$53.55 E0978

- 1321BK Hook and Loop Seat Positioning Strap
.....\$42.84 E0978

CHEST POSITIONING STRAPS

- 9515 Push Button Style Chest Positioning Strap
..... \$53.55 E0960



Yes, you can.®

ACCESSORIES AND STOCK OPTIONS

BACK ACCESSORIES

1160 Carrying Pocket \$57.83
 3046 Head Support For Recliner\$125.31 E0966

STC Stencil Back Upholstery \$31.06
 3/4" block letters. Only available with vinyl upholstery.

RB Dual Back Liner\$52.48

SEAT ACCESSORIES

9120 Solid Seat Insert\$104.96

RS Dual Seat Liner\$52.48

ARM ACCESSORIES

NRAS Non - Removable Arm With Screws\$10.71

1300 Utility Tray\$168.15 E0950

MISCELLANEOUS ACCESSORIES

2219 IV Rod, Telescopic Recliner \$117.81 K0105

1100 Tool Kit No Charge

2226 Front Anti - Tippers \$88.89

FRONT RIGGING OPTIONS

FRONT RIGGING

T93HA Hemi Sw. Footrests Aluminum Footplates
 w/ Heel Loops \$99.13

AHL4C Hemi Smart Leg Articulating Legrests
 Composite Footplates \$341.09 K0053

T93HC Hemi Sw. Footrests Composite Footplates
 w/ Heel Loops \$99.13

AHS4A Hemi Smart Leg Articulating Legrestst-
 14-15" W Chairs \$341.09 K0053

T94HA Hemi Sw. Elev Legrests Aluminum Footplates
 & Padded Calf Pad \$191.86 K0195

H93E Hemi Ext Swingaway Footrests \$197.19

T94HC Hemi Sw. Elev Legrests Composite Footplates
 & Padded Calf Pad \$191.86 K0195

T93HAB Hemi Bolt-Thru Footrests, Aluminum
 Footplates.....\$95.93

T93HS Hemi Sw. Footrests Aluminum Footplates
 w/ Heel Loops- 14"W Chairs\$95.93

T93HCB Hemi Bolt-Thru Elevating Legrests,
 Aluminum Footplates\$175.00

T94HS Hemi Sw. Elevating Legrests Aluminum
 Footplates & Padded Calf Pad- 14" W Chairs
\$186.53 K0195

AHL4A Hemi Smart Leg Articulating Legrests
 Aluminum Footplates \$341.09 K0053

FOOTPLATE OPTIONS

1350 Extra Large Aluminum Footplates\$84.21 K0041

AT5543 Adjustable Angle Flip - Up Footplates
\$71.42 K0040



Yes, you can.®

FOOTREST / LEGREST ACCESSORIES

- 1820 Non - Removable Device\$12.79
- 1337 Calf Strap\$38.37 K0038
Not available with AHL4A or AHL4C elevating legrests.
- 1348 Legrest Panel\$66.00

- 1346 Heel Loops\$41.57 E0951
Standard on footrests.
- L8689 Omit Heel LoopsLess \$10.00
- 1319 Black Non - Padded Calf Pads\$10.71
Not available with AHL4A or AHL4C elevating legrests.

Seat to Floor Matrix 8" Casters

SEAT-TO-FLOOR HEIGHT	FRONT CASTER MOUNTING POSITION	REAR WHEEL SIZE	REAR WHEEL MOUNTING POSITION
17 1/2"	Top	20"	1 Position Only
18 1/2"	Middle	22"	1 Position Only
19 1/2"	Bottom	24"	1 Position Only

Weight Limitation 250 Pounds

Invacare recommends all wheelchairs to be fitted to the consumer under the guidance of a healthcare professional knowledgeable of the particular limitations of the product.