



Invacare AVIVA FX Power Wheelchair
with Captains and Solid Seat
US PRICE LIST AND ORDER FORM - Price Effective July 15, 2025



CONTACT ARACENT HEALTHCARE CUSTOMER SERVICE: 888.550.7750 / SALES@ARACENT.COM

REQUIRED INFORMATION

☐ Quote ☐ Order

Date of Order: _____ Dealer Account #: _____

Dealer Name: _____

PO #: _____ Tag _____

Purchasing Contact: _____

Phone: _____ Fax: _____

E-mail: _____

ATP/Therapist: _____

Ship to Address: _____

City: _____ State: _____

Zip: _____ Client Gender: ☐ M ☐ F

Special Client Conditions: _____

REQUIRED MEASUREMENTS

To ensure system is accurately configured please fill in all required measurements.

Client Height:		
Client Weight:		
Use the table below to calculate your Seat Pan to Footplate Measurement. Then select the appropriate range in the legrest section.	Solid Seat	Captain Seat
Knee to Heel:		
Cushion Thickness: (Subtract)		4"
Cushion Compression:(Add)		

Notes:

HCPCS CODES	K0848, K0849
BASE WIDTH	24.25"
BASE LENGTH	42.1"
DRIVE WHEEL SIZE	14"
CASTER WHEEL SIZE	8"
BASE ONLY TURNING RADIUS	6.8"
DRIVETRAIN	4-pole Sealed Housing
RANGE ⁽¹⁾	19.7 miles

BATTERIES	Group 24
BATTERY ISOLATION SWITCH	Yes
BATTERY WEIGHT	38.5 lbs (each)
BATTERY CHARGER	8A, Off Board
WEIGHT CAPACITY	300 lbs
GROUND CLEARANCE	2.5"
MAXIMUM SPEED NARROW/HD	5.8 mph
MAXIMUM SPEED STANDARD	7.5 mph

⁽¹⁾ Actual driving range and speed may vary due to factors such as user weight, type and grade of terrain, battery condition & charge level; type and condition of the drive wheels, and variations in the drive and control systems.

SOLID SEAT



CAPTAIN SEAT



Invacare AVIVA FX Captain and Solid Seat

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

INVACARE POWER BASES

- ☐ IFX-20C AVIVA FX Base for Captain Seat
HCPCS code K0849. \$8,290
- ☐ IFX-20S AVIVA FX Base for Solid Seat Pan
HCPCS code K0848. \$8,290

SELECT SEAT TO FLOOR HEIGHT (inches)		STF Height	Captain's Seat		Solid Seat	
			No Elevate	With Elevate	No Elevate	With Elevate
☐ FLOW	☐ FLOW-S	Super Low ⁽¹⁾	19 7/8	21 5/8	16 1/8	17 7/8
☐ FLOW	☐ FLOW-S	Low	21 1/8	22 7/8	17 3/8	19 1/8
☐ FMED	☐ FMED-S	Medium	22 1/8	23 7/8	18 3/8	20 1/8

1. FLOW / FLOW-S \$250 up charge.

SEAT STYLE

- ☐ SEMIR Semi-Recline (35°) Captain's Seat STD
- ☐ SEMIRB Semi-Recline (35°) with Solid Seat \$100

SEAT COVER

- ☐ BSC Black Slip Cover. \$295

1. Includes seat, back, and headrest cover

USER WEIGHT LIMITS

- ☐ U250 Weight Capacity up to 250 lbs STD
- ☐ U300 Weight Capacity up to 300 lbs STD

PERFORMANCE MOTORS

- ☐ HSMTR Performance Motors - 7.5 mph⁽¹⁾⁽²⁾ \$1,400

1. Extended high speed use on rugged terrain is not recommended.

2. Performance/High Speed Motors may have less torque and may have less range than standard speed motors.

TRANSPORT TIE DOWN

- ☐ TRBKTS Wheelchair Transport Brackets⁽¹⁾ N/C

1. For unoccupied use only

TIRE OPTIONS

- ☐ B1431-3 14" x 3" Black Tire with Gel Foam Inserts STD
- ☐ B1430-3 14" x 3" Black Tire - Pneumatic Tire N/C

FORK OPTIONS

- ☐ AB1228-3 8" x 2.5" Black Casters w/ Solid Inlay. N/C

INSTALLATION OF BATTERIES

- ☐ M24 SLD G FT MK Group 24 FT Gel Battery (ea) Quantity()
HCPCS code E2363 \$365
- ☐ M34 SLD G MK Group 34 Gel Battery (ea) Quantity()
HCPCS code E2359 \$305

BATTERY CHARGER OPTION

- ☐ 110CHARGER 110 Volt Battery Charger. N/C
- ☐ L8900 Omit Battery Charger (\$80)

LIGHTS AND INDICATORS

- ☐ LIGHTS LED Lights and Indicators⁽¹⁾ \$700

1. Must order REM216 or REM400

FRAME FINISH

- ☐ 170P  Obsidian Black. N/C
- ☐ 171P  Ruby Red N/C
- ☐ 172P  Peacock Blue N/C
- ☐ 173P  Pink Flamingo N/C
- ☐ 174P  Purple Passion N/C
- ☐ 175P  Denim Navy N/C
- ☐ 176P  Lime Twist N/C
- ☐ 177P  Arctic Blue N/C

RIM INSERT COLORS - MUST PICK ONE

- ☐ 170PR  Obsidian Black. N/C
- ☐ 171PR  Ruby Red N/C
- ☐ 172PR  Peacock Blue N/C
- ☐ 173PR  Pink Flamingo N/C
- ☐ 174PR  Purple Passion N/C
- ☐ 175PR  Denim Navy N/C
- ☐ 176PR  Lime Twist N/C
- ☐ 177PR  Arctic Blue N/C
- ☐ 178PR  Banana Split N/C
- ☐ 179PR  Orange Marmalade N/C
- ☐ 180PR  Morning Mist N/C

ELEVATING SEAT

- ☐ ESM12C Elevate lift 12" \$2,995
Weight capacity up to 300lbs
- ☐ RDE Reduced Drive for Elevating Seat⁽¹⁾ \$400

1. If reduced drive is not selected the system will be in drive lock-out when elevated.

ELECTRONICS FOR ELEVATING SEAT

- ☐ SFCB Single Function Electronics⁽¹⁾ STD
- Select Switch:
- ☐ SPBSO Single Push Button
- ☐ STO Single Toggle Switch
- ☐ SFCBA Single/Dual Function Drive Control⁽¹⁾
HCPCS code E2310 \$1,795

Invacare AVIVA FX Captain and Solid Seat

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

ADDITIONAL SWITCH FOR POWER SEATING

☐ SPBSA	Single Push Button	\$150
☐ STA	Single Toggle Switch	\$150
☐ ASL320R	ASL Pulse Mechanical Switch - Red	\$125
☐ ASL320G	ASL Pulse Mechanical Switch - Green	\$125
☐ ASL320BK	ASL Pulse Mechanical Switch - Black	\$125
☐ SESG	Stealth Egg Switch - Green	\$135
☐ SESB	Stealth Egg Switch - Black	\$135

MOUNTING FOR POWER SEATING SWITCH

☐ ESML	Left Side (Armpad)
☐ ESMR	Right Side (Armpad)
☐ ESMJL	Joystick Left (Mounts to Joystick Bracket)
☐ ESMJR	Joystick Right (Mounts to Joystick Bracket)
☐ ESMD	Do not mount

SEAT WIDTH SETTING

☐ W18	18" Wide (Actual Width 17.7" Solid / 19" Captain)	N/C
☐ W20	20" Wide (Actual Width 18.5" Solid / 20" Captain)	N/C
☐ W22	22" Wide (Actual Width 20.5" Solid / 22" Captain)	N/C

SEAT DEPTH SETTING

☐ D16	16" Deep	N/C
☐ D17	17" Deep	N/C
☐ D18	18" Deep	N/C
☐ D19	19" Deep	N/C
☐ D20	20" Deep	N/C
☐ D21	21" Deep	N/C

Seat Width and Depth Compatibility

	16D	17D	18D	19D	20D	21
18W	✓	✓	✓	✓	N/A	N/A
20W	N/A	N/A	N/A	✓	✓	✓
22W	N/A	N/A	N/A	✓	✓	✓

NON-EXPANDABLE SINGLE DRIVE JOYSTICKS

 REM110S	LED Remote ⁽¹⁾⁽²⁾N/C <i>Control of seating through separate switch.</i>
 REM210S	LED Non-Expandable Remote/Joystick ⁽¹⁾⁽³⁾N/C <i>Control of power seating through drive control</i>

1. Not available with lights and indicators.

2. Select switch in Switch For Power Seating.

3. Must select SFCBAC for through joy operation.

EXPANDABLE JOYSTICKS AND CONTROLLER

☐ EXPC	Expandable Controller HCPCS code E2377	\$700
☐ PWH	Harness Required for Expandable System HCPCS code E2313	\$450
☒ GTRAC-LX	LiNX G-TRAC Module	STD
☐ REM211	LED Expandable Remote ⁽¹⁾	N/C
☐ REM216	LED Expandable Remote with Lights ⁽¹⁾⁽²⁾	\$350
☐ REM400	Color 3.5" Touch Screen Remote ⁽¹⁾⁽²⁾	\$1,000
☐ SW400TGL	Toggle Kit for REM400 ⁽³⁾	\$395

1. Must choose EXPC Expandable Controller and PWH Harness.

2. Can control lighting, must choose lighting package separately.

3. Adds 2 toggle switches to REM400.

MOUNTING FOR JOYSTICKS

☐ FR	Fixed Joystick Mounting - Right	STD
☐ FL	Fixed Joystick Mounting - Left	STD
☐ MRR	Maxx Resolve Swing-Away Joystick Mount (Right Hand) HCPCS code E1032	\$395
☐ MRL	Maxx Resolve Swing-Away Joystick Mount (Left Hand) HCPCS code E1032	\$395
☐ SAQR	Motion Height Adjustable Swing-Away Quad Link (Right Hand) HCPCS code E1032	\$235
☐ SAQL	Motion Height Adjustable Swing-Away Quad Link (Left Hand) HCPCS code E1032	\$235

LINX ELECTRONIC ACCESSORIES

☐ LAK	LiNX Access Key ⁽¹⁾	\$100
☐ MUSB	Motion Dual USB Charger	N/C

1. Required for programming

JOYSTICK TOPS⁽¹⁾

HCPCS code E2323

☐ PC101A	Bodypoint U Shaped Handle 3"	\$120
☐ PC102A	Bodypoint U Shaped Handle 4"	\$120
☐ PC107A	Bodypoint Rubber Dome	\$120
☐ PC110A	Grooved Mushroom Joystick Handle	\$120

1. To remove/disengage the joystick knob from the REM400 Remote pull straight up on the joystick (DO NOT TWIST) otherwise damage may occur and may void the warranty.

LAP BELTS

☐ BELT48	Push Button Style Lap Belt, Length 48"	STD
☐ BELT60	Push Button Style Lap Belt, Length 60"	N/C
☐ BELT71	Push Button Style Lap Belt, Length 71"	N/C
☐ PLB2	Bodypoint Padded Lap Belt 2 Point ⁽¹⁾	\$190
☐ PLB4	Bodypoint Padded Lap Belt 4 Point ⁽¹⁾	\$235

1. 67" length only option

ARM RESTS

HCPCS code E0973

☒ LVAHF	Left Full Length, Adjustable Height	\$115
☒ RVAHF	Right Full Length, Adjustable Height	\$115

Invacare AVIVA FX
Captain and Solid Seat

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

MATRIX SEAT CUSHIONS - FOR SOLID SEAT SYSTEMS

Matrix PS Seat Cushion 14" - 20" wide/reversable cover
 HCPCS code E2605. \$345

☐ PS1816 ☐ PS1818 ☐ PS2018 ☐ PS2020

Matrix PS Seat Cushion 21"-22" wide
 HCPCS code E2606. \$545

☐ PS2218 ☐ PS2220

Matrix PSP Seat Cushion 14" - 20" wide/reversable cover
 HCPCS code E2607. \$395

☐ PSP1816 ☐ PSP1818 ☐ PSP2018 ☐ PSP2020

Matrix PSP Seat Cushion 21"-22" wide
 HCPCS code E2608. \$595

☐ PSP2218

Matrix Vi Seat Cushion 14" - 20" wide
 HCPCS code E2607. \$425

☐ MA1816-VI ☐ MA1818-VI ☐ MA2018-VI ☐ MA2020-VI

Matrix Vi Seat Cushion 21" - 22" wide
 HCPCS code E2608. \$475

☐ MA2218-VI ☐ MA2220-VI

Matrix Libra Seat Cushion 14" - 20" wide
 HCPCS code E2624. \$485

☐ LC1816 ☐ LC1818 ☐ LC2018 ☐ LC2020

Matrix Libra HD Seat Cushion 21" wide
 HCPCS code E2624. \$525

☐ LC2218 ☐ LC2220

Matrix Libra Adjust Seat Cushion 14" -20" wide
 HCPCS code E2624. \$625

☐ LC1816-ADJ ☐ LC1818-ADJ ☐ LC2018-ADJ ☐ LC2020-ADJ

Matrix Libra Adjust Seat Cushion 22" deep
 HCPCS code E2624. \$735

☐ LC2218-ADJ ☐ LC2220-ADJ

Matrix Multi Cushion 14" - 20" wide w/reversable cover
 HCPCS code E2607. \$445

☐ MM1816 ☐ MM1818 ☐ MM2018 ☐ MM2020

Matrix Multi Cushion 22" wide w/reversable cover
 HCPCS code E2608. \$595

☐ MM2218 ☐ MM2220 ☐ MM2222

☐ OCCC "ON-CHAIR" Cushion Credit⁽¹⁾ (\$90)

1. Not available on Matrix PS or PSP Cushions.

STANDARD FRONT RIGGING


☐ FCMP Fixed Center Mount Platform⁽¹⁾ N/C
☐ FIFP Rubber Coated Footplate (11.5" W x 10") N/C
☐ FBIFP Large Rubber Coated Footplate (17" w x 13"d) ... N/C

Select seat pan to footplate:

☐ FSPTF3 10.25" - 14.25" N/C
☐ FSPTF4 14.25" - 20.25" N/C

Set at:

☐ FSA70 Set at 70° N/C
☐ FSA83 Set at 83° N/C
☐ FSA90 Set at 90° N/C
☐ FSA97 Set at 97° N/C



1. 3" ground clearance required.

MANUAL ELEVATING CENTER MOUNT LEGREST

☐ MECML Manual Elevating Center Mount Legrest
 With Plantar Flex \$595

Select legrest set up position:

☐ 90SU Set at 90° (elevating range 90-40°) STD
☐ 110TU Tucked at 110° (elevating range 110-60°)⁽¹⁾ N/C

Rubber Coated Footplate:

☐ MRCP Rubber Coated Footplate 11.5" W x 10" D \$165

Select seat pan to footplate for One Piece Foot Plate

☐ CSPTF1 12.75" - 16.25" N/C
☐ CSPTF2 15.25" - 20.25" N/C

Individual Footplates:

☐ IFSL Small - (5" x 7.5" deep) - Left \$165
☐ IFSR Small - (5" x 7.5" deep) - Right \$165
☐ IFML Medium - (5.5" x 9.5" deep) - Left \$165
☐ IFMR Medium - (5.5" x 9.5" deep) - Right \$165
☐ IFLL Large - (6" x 11.5" deep) - Left \$165
☐ IFLR Large - (6" x 11.5" deep) - Right \$165


Select seat pan to footplate for Individual Footplates

☐ CSPTF3 10" - 15" N/C
☐ CSPTF4 12.5" - 19.5" N/C

Optional:

☐ CICP Individual calf pads \$200

1. Not available on seat depths 16", 17".



Invacare AVIVA FX
Captain and Solid Seat

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

MAXX STYLE FRONT RIGGING

☐ MSFCMFP Maxx Style Fixed Center Mount
Foot Platform⁽¹⁾⁽²⁾\$395

☐ MIFP Comes with a rubber coated footplate
11.5" W x 10" DN/C

☐ MICP Comes standard with individual calf padsN/C

Select seat pan to footplate:

☐ MSPTF3 10.25" - 14.25"N/C

☐ MSPTF14 14.25" -20.25"N/C


Set at:

☐ MSA70 Set at 70°N/C

☐ MSA83 Set at 83°N/C

☐ MSA90 Set at 90°N/C

☐ MSA97 Set at 97°N/C



1. 3" ground clearance required.

2. Fixed Center Mount will be set at factory seating to prevent castor interference.

FOOT PLATFORM OPTIONS ⁽¹⁾



☐ SFFP Small Foot Platform (12"w x 8.5"d)\$180

☐ MFFP Medium Platform (12"w x10.75"d)\$180

☐ LFFP Large Platform (14.5"w x13"d)\$250

☐ XLFFP Extra Large Platform (17"w x13"d)\$250

☐ XXLFFP Extra Extra Large Platform (21"w x13"d)\$250

1. May need to adjust leg position to avoid interference with front casters

ACCESSORIES

☐ PH X-Grip Phone Holder\$95

☐ UCH Universal Cup Holder ⁽¹⁾\$35

☐ CCCC Cane and Crutch Holder\$395

☐ CSO2 Non-Articulating O2 Gas Holder\$250

☐ CFH Flag Holder\$65

☐ CRVMR Rear View Mirror Right Mounted\$45

☐ CRVML Rear View Mirror Left Mounted\$45

☐ TAG571850 Captain Seat Wheelchair Cover\$105

☐ AMC Accessory Mount with Clip\$50

☐ CSWH Captain Seat Walker Holder ⁽²⁾\$200

1. Not available on the joystick side

2. Maximum load is 20 lbs

MILITARY DECAL ⁽¹⁾

☐ DECALA Army Decal 3" Diameter\$50

☐ DECALAF Air Force Decal 3" Diameter\$50

☐ DECALN Navy Decal 3" Diameter\$50

☐ DECALM Marines Decal 3" Diameter\$50



1. Comes uninstalled off chair.

NOTES

RETAIL \$ _____

SUBMIT

PRINT FORM

SAVE

CLEAR FORM