



**TURNING,
DRIVING
& RIDE
COMFORT!**

ULTRA LOW MAXX CG POWER POSITIONING SYSTEM ON STORM TORQUE RX20 BASE

US PRICE LIST AND ORDER FORM - Price Effective May 1, 2025



PRINT FORM

SAVE

CLEAR FORM

Aracent Healthcare Customer Service:

1-888-550-7750

sales@aracent.com

REQUIRED INFORMATION

☐ Quote ☐ Order

Date of Order: _____ Dealer Account #: _____

Dealer Name: _____

PO #: _____ Tag: _____

Purchasing Contact: _____

Phone: _____ Fax: _____

E-mail: _____

ATP/Therapist: _____

Ship to Address: _____

City: _____ State: _____


Zip: _____ Client Gender: ☐ M ☐ F

Special Client Conditions: _____

REQUIRED MEASUREMENTS

To ensure system is accurately configured please fill in all **REQUIRED*** measurements

Client Height*	
Client Weight (lbs)*	
A. Seat to Shoulder	
B. Trunk Depth	
C. Chest Width	
D. Knee to Back	
E. Seat to Top of Head	
F. Back to Joystick	
G. Seat to Elbow	
H. Hip Width	
I. Knee to Heel*	
Cushion Thickness*	



BASE OPTIONS

NEW!

☐ SRX-20R Storm **Torque** RX20 Base for Rehab Seat System⁽¹⁾⁽²⁾
HCPCS code K0848 \$10,995

☐ SRX-20SP Storm **Torque** RX20 Base for Single Actuator System⁽¹⁾
HCPCS code K0856 \$7,995

☐ SRX-20MP Storm **Torque** RX20 Base for Multiple Actuator System⁽¹⁾
HCPCS code K0861 \$7,995



-Enhanced Front Caster Housing and Bearings
-New Front Caster Design
-Improved double-sided Front Forks
-Upgraded Rear Suspension
-Optimized Linx Parameters

1. Comes Standard with Lighting Package
2. Must select Ultra Low Maxx Rehab Seat (ULMRSN) system

BASE WIDTH

☐ NARROWFRM Narrow Base 22.4" (MK 50 Batteries) N/C
☐ STDFRM Standard Base 24.4" (MK 24 FT Batteries) N/C

SEAT-TO-FLOOR HEIGHT⁽¹⁾

☐ SLOW 16.75" Seat-To-Floor⁽²⁾ \$250
☐ LOW 17.25" Seat-To-Floor N/C
☐ MED 18.25" Seat-To-Floor N/C
☐ TALL 19.25" Seat-To-Floor N/C

1. Seat-to-floor height is measured at the front of the seat pan based on a standard seat size. Seat-to-floor height may vary +/- .25" dependent on condition of the power base and seat depths greater than 18"
2. This seat to floor height is only available with the Group 24 Battery Tray (24TRY) and on 17" seat depths and greater.

USER WEIGHT LIMIT

☐ UW250 Weight Capacity up to 250lbs. STD
☐ U300 Weight Capacity up to 300lbs STD

1. Required for Elevate Systems and/or NARROWFRM

WHEELCHAIR TRANSPORT BRACKETS

☐ TRBKTS Wheelchair Transport Brackets⁽¹⁾ STD
☐ TRRO Transport Ready Option⁽²⁾⁽⁴⁾ \$275
☐ TRRO-E Transport Ready Option for Elevate⁽²⁾⁽³⁾⁽⁴⁾ \$275

1. For unoccupied use only
2. TRRO meets ISO-7176 Part 19
3. Required for Elevate Systems to meet ISO-7176 Part 19, includes two additional tie-down loops on rail.
4. Must select a Motion Concepts or Matrix Elan headrest with hardware and Motion Concepts Standard Rehab, High Back or Matrix Upgradeable back. Comes with pelvic safety belt standard, please select TROBELT on page 5.

TIRE OPTIONS

☐ AB1431-3 14" x 3" Black Tire with Gel Foam Inserts STD
☐ AB1430-3 14" x 3" Black Tire - Pneumatic Tire N/C

CASTER OPTION

NEW!

☐ AB1228-3 8" x 2.5" Black Ribbed Casters with Solid Inlay .. STD

BATTERY CHARGER OPTION

☐ 110CHARGER 110 Volt Battery Charger N/C
☐ L8900 Omit Battery Charger (\$80)

HIGHSPEED MOTORS AND ENHANCED DRIVING

☐ HSMTR Performance Motors 7.5 mph⁽¹⁾⁽²⁾⁽³⁾ \$1,400
☐ GTRAC-LX LINX G-TRAC Module⁽⁴⁾
HCPCS code K0108 \$995

1. Not available on the NARROWFRM.
2. Extended high speed use on rugged terrain is not recommended.
3. Performance/High Speed Motors may have less torque and may have less range than standard speed motors.
4. Recommended for use with Specialty Drive Controls

Ultra Low Maxx CG Power Positioning System

Storm Torque RX20

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

FRAME FINISH

☐ 170P		Obsidian Black.....	N/C
☐ 171P		Ruby Red.....	N/C
☐ 172P		Peacock Blue.....	N/C
☐ 173P		Pink Flamingo.....	N/C
☐ 174P		Purple Passion.....	N/C
☐ 175P		Denim Navy.....	N/C
☐ 176P		Lime Twist.....	N/C
☐ 177P		Arctic Blue.....	N/C

RIM INSERT COLORS - MUST PICK ONE

☐ 170PR		Obsidian Black.....	N/C
☐ 171PR		Ruby Red.....	N/C
☐ 172PR		Peacock Blue.....	N/C
☐ 173PR		Pink Flamingo.....	N/C
☐ 174PR		Purple Passion.....	N/C
☐ 175PR		Denim Navy.....	N/C
☐ 176PR		Lime Twist.....	N/C
☐ 177PR		Arctic Blue.....	N/C
☐ 178PR		Banana Split.....	N/C
☐ 179PR		Orange Marmalade.....	N/C
☐ 180PR		Morning Mist.....	N/C

NON-EXPANDABLE SINGLE/DUAL JOYSTICK

☐ REM215N	LED Non-Expandable Remote with Lights.....	N/C
<i>Three drive profiles, control of two seat functions and light operation</i>		

EXPANDABLE JOYSTICKS⁽¹⁾

☐ REM216N	LED Expandable Remote/Joystick with Lights ⁽³⁾	N/C
<i>Three drive profiles, control of up to 6 actuators and light operation. Not for use with specialty controls</i>		
☐ REM400	Color 3.5" Touch Screen Remote/Joystick ⁽²⁾	\$1,000
<i>Four drive profiles, control of up to 6 actuators through the drive control and light operation.</i>		
☐ SW400TGL	Toggle Kit for REM400 ⁽³⁾	\$200

1. Must choose EXPC Expandable Controller and PWH Harness.
2. REM400 or REM500 is required when using specialty controls, choose specialty controls purchased separately.
3. Adds 2 toggle switches to REM400.

EXPANDABLE CONTROLLER AND HARNESS

☐ EXPC	Expandable Controller ⁽¹⁾ HCPCS code E2377.....	\$700
☐ PWH	Harness for Expandable Controller HCPCS code E2313.....	\$450

1. This item is required when ordering power tilt, power recline, and power elevating legrest together OR when ordering an Alternative Drive Control.

EXPANDABLE SPECIALTY DISPLAY

☐ REM500	LiNX Color 3.5" Touch Screen Display ⁽¹⁾⁽²⁾⁽³⁾	\$1,500
---------------------------------	---	---------

Select mounting position for REM500

☐ R5MR	Right.....	N/C
☐ R5ML	Left.....	N/C
☐ R5RPB	Rear Push Bar.....	\$140

1. Must choose EXPC Expandable Controller and PWH Harness.
2. Driver Control not included.
3. REM400 or REM500 is required when using specialty controls, choose specialty controls purchased separately.

SPECIALTY PROPORTIONAL CONTROLS⁽¹⁾⁽²⁾

☐ CREM	LiNX Compact Remote HCPCS code E2373...	\$1,500
☐ CREM-LF	LiNX Compact Remote Low Force HCPCS code E2373.....	\$1,500

1. Requires REM400 or REM500
2. Select Mount in Mounting For Driver Control section

ATTENDANT DRIVE CONTROL CHOICES

☐ ACU	Proportional Attendant Control ⁽¹⁾ HCPCS code E2331.....	\$900
------------------------------	--	-------

Select mounting position for ACU

☐ RPBMR	Right
☐ RPBML	Left
☐ RPBMC	Center

1. Requires Expandable Controller.

THIRD PARTY SPECIALTY CONTROLS

Motion Concepts Install recommended

Installed at Motion Concepts: ☐ YES ☐ NO

☐ Type: _____

☐ Specialty Control Mount:	☐ Right	☐ Left
☐ Display Mount:	☐ Right	☐ Left
☐ Mode Switch Mount:	☐ Right	☐ Left
☐ On/Off Mount:	☐ Right	☐ Left

Notes

Ultra Low Maxx CG Power Positioning System

Storm Torque RX20

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

JOYSTICK MOUNT

Arm Mounted

- ☐ FR Fixed Joystick Mounting - Right STD
- ☐ FL Fixed Joystick Mounting - Left STD
- ☐ MRR Maxx Resolve Swing-Away Joystick Mount
(Right Hand) HCPCS code E1028\$395
- ☐ MRL Maxx Resolve Swing-Away Joystick Mount
(Left Hand) HCPCS code E1028\$395
- ☐ SAQR Motion Height Adjustable Swing-Away
Quad Link (Right Hand)
HCPCS code E1028\$235
- ☐ SAQL Motion Height Adjustable Swing-Away
Quad Link (Left Hand)
HCPCS code E1028\$235

Rear Push Bar Mounted

- ☐ RPBMR Right\$140
- ☐ RPBML Left\$140
- ☐ RPBMC Center\$140

JOYSTICK KNOB - OPTIONAL

HCPCS code E2323

- ☐ PC101A-20 BodyPoint 3" U-Shaped\$120
- ☐ PC102A-20 BodyPoint 4" U-Shaped\$120
- ☐ PC104A-20 BodyPoint Straight Handle\$120
- ☐ PC107A-20 BodyPoint Domed Rubber\$120
- ☐ PC110A-20 BodyPoint Grooved Mushroom\$120

1. To remove the joystick knob from the REM400, pull straight up on the joystick (DO NOT TWIST) otherwise damage may occur and may void the warranty.

ELECTRONICS - EXPANDABLE MODULES

- ☐ OUT500 Output Module\$500
- ☐ IN200 LiNX Switch Input Module with
Power ON/OFF for REM200 Series⁽¹⁾⁽²⁾\$500
- ☐ IN500H LiNX 9 Pin Input Module for
Alternative Driver Controls⁽²⁾⁽³⁾⁽⁴⁾⁽⁶⁾\$1,500
- ☐ IN500 LiNX 9 Pin Input Module with
Sip-N-Puff Nozzle⁽²⁾⁽³⁾⁽⁴⁾⁽⁶⁾ HCPCS code E2325.\$1,500
- ☐ PKG32666 Therafin Sip-N-Puff Breath Tube Kit\$380
HCPCS code E2326
- ☐ COVER-L Large Recline Shroud Cover⁽⁵⁾N/C

1. Not for use with Alternative Drive Controls. (Uses switch ports).

2. Only one input module per chair.

3. Input module also required for any 9-pin specialty control.

4. 9 pin specialty control and Sip-N-Puff cannot both be used together.

5. May be required when choosing OUT500, IN200 or IN500 options.

6. Requires REM400 or REM500 for compatibility.

LINX ELECTRONIC ACCESSORIES

- ☐ LAK LiNX Access Key⁽¹⁾\$100

1. Required for programming

NOTES

Ultra Low Maxx CG Power Positioning System

Storm Torque RX20

NOTE: Systems Include Seat Pan, Back Pan, and Standard Armrests. They are designed for use with a 2" - 3" cushion

REHAB SEAT SYSTEM⁽¹⁾

☐ ULMRSN Ultra Low Maxx Rehab SeatN/C

1. Must order Rehab Seat Base (SRX-20R)

POWER POSITIONING SYSTEMS⁽¹⁾

☐ ULMT Ultra Low Maxx CG Tilt (50°) Only System⁽²⁾
HCPCS code E1002\$6,595

☐ ULMTRE Ultra Low Maxx CG Tilt (50°) & Power Recline (170°) w/Extended Shear Reduction(ESR)⁽²⁾
HCPCS code E1007\$10,795

☐ ULMTRP Ultra Low Maxx CG Tilt (50°) & Power Recline (170°) w/Power Sliding Back(PSB)⁽²⁾⁽³⁾
HCPCS code E1008\$10,995

☐ ULMRE Ultra Low Maxx Power Recline (170°) w/Extended Shear Reduction(ESR)
HCPCS code E1004\$5,450

☐ ULMRP Ultra Low Maxx Power Recline (170°) w/Power Sliding Back(PSB)⁽³⁾
HCPCS code E1005\$6,250

☐ ULME Ultra Low Maxx 12 " Elevating Lift Seat Only⁽⁴⁾ HCPCS code E2298\$4,995
Weight capacity up to 250lbs

1. Not available with Rehab Seat Base (SRX-20R)

2. Tilt will be reduced to 45° if adding an Elevating Lift Seat Module

3. Requires Multi-Function M610 through Drive Control (MFCBD). Not available on Multi Function Act4-LiNX Electronics.

4. Not available on seat sizes greater than 20" wide by 20" deep.

ADD: POWER ADJUSTABLE SEAT HEIGHT MODULE

☐ ESM12 12" Power Adjustable Seat Height Module⁽¹⁾⁽²⁾⁽³⁾
HCPCS code E2298\$2,495
Weight capacity up to 250lbs

1. Uses an integrated Tilt/Lift mechanism, tilt is 45° and the seat to floor height is not affected by the Elevating Lift Seat.

2. Not available on seat sizes greater than 20" wide by 20" deep

3. For weights greater than 175lbs tilt while elevated is not available.

4. Requires Multi-Function M610 through Drive Control (MFCBD). Not available on Multi Function Act4-LiNX Electronics.

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

MODULES - MUST SELECT SYSTEM TO ADD MODULES

☐ PRETILT 5° Fixed Pre-Tilt⁽¹⁾N/C
☐ PRM8 8° Precline Module⁽²⁾N/C
☐ PRM30 30° Precline Module⁽²⁾
HCPCS code K0108\$500
☐ ERSSP Programable Electronic Recline Function
Start & Stop⁽³⁾\$400

Specify:

Start: _____degrees Stop: _____degrees

1. This Fixed Pre-tilt option adds approximately 2" to the front seat to floor height.

2. Total degrees of recline is 170° ie. if 30° of precline is used, then 140° of recline will be available. Must order Recline System

3. Requires Multi-Function M610 through Drive Control (MFCBD). Set range between 90° and 170°. Will be set to factory settings unless specified.

ADULT SEAT WIDTH SETTINGS

☐ SW116 16" Wide (Adjustable from 16" to 20")N/C
☐ SW117 17" Wide (Adjustable from 16" to 20")N/C
☐ SW118 18" Wide (Adjustable from 16" to 20")N/C
☐ SW119 19" Wide (Adjustable from 16" to 20")N/C
☐ SW120 20" Wide (Adjustable from 16" to 20")N/C
☐ SW219 19" Wide (Adjustable from 19" to 22")N/C
☐ SW220 20" Wide (Adjustable from 19" to 22")N/C
☐ SW221 21" Wide (Adjustable from 19" to 22")N/C
☐ SW222 22" Wide (Adjustable from 19" to 22")N/C

ADULT ULTRA RAIL SEAT DEPTH SETTINGS

☐ SD116 16" Deep (Adjustable from 16" to 21")N/C
☐ SD117 17" Deep (Adjustable from 16" to 21")N/C
☐ SD118 18" Deep (Adjustable from 16" to 21")N/C
☐ SD119 19" Deep (Adjustable from 16" to 21")N/C
☐ SD120 20" Deep (Adjustable from 16" to 21")N/C
☐ SD121 21" Deep (Adjustable from 16" to 21")N/C
☐ SD219 19" Deep (Adjustable from 19" to 23")N/C
☐ SD220 20" Deep (Adjustable from 19" to 23")N/C
☐ SD221 21" Deep (Adjustable from 19" to 23")N/C
☐ SD222 22" Deep (Adjustable from 19" to 23")N/C
☐ SD223 23" Deep (Adjustable from 19" to 23")N/C

SMALL ADULT SEAT WIDTH SETTINGS⁽¹⁾

☐ SW314 14" Wide (Adjustable from 14" to 15")N/C
☐ SW315 15" Wide (Adjustable from 14" to 15")N/C

1. Weight limit limited to 175 lbs. If recline is selected it may be limited to 165° of recline. You must select both small adult seat width and depth together.

SMALL ADULT SEAT DEPTH SETTINGS⁽¹⁾

☐ SD314 14" Deep (Adjustable from 14" to 18")N/C
☐ SD315 15" Deep (Adjustable from 14" to 18")N/C
☐ SD316 16" Deep (Adjustable from 14" to 18")N/C
☐ SD317 17" Deep (Adjustable from 14" to 18")N/C
☐ SD318 18" Deep (Adjustable from 14" to 18")N/C

1. Weight limit limited to 175 lbs. If recline is selected it may be limited to 165° of recline. You must select both small adult seat width and depth together.

Ultra Low Maxx CG Power Positioning System

Storm Torque RX20

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

MATRIX SEAT CUSHION OPTIONS

Matrix PS Seat Cushion 14" - 20" wide/reversable cover

HCPCS code E2605 \$345

☐ PS1414 ☐ PS1416 ☐ PS1616 ☐ PS1618 ☐ PS1620
☐ PS1716 ☐ PS1718 ☐ PS1720 ☐ PS1816 ☐ PS1818
☐ PS1820 ☐ PS2016 ☐ PS2018 ☐ PS2020

Matrix PS Seat Cushion 21"-22" wide

HCPCS code E2606 \$545

☐ PS2218 ☐ PS2220 ☐ PS2222

Matrix PSP Seat Cushion 14" - 20" wide/reversable cover

HCPCS code E2607 \$395

☐ PSP1414 ☐ PSP1416 ☐ PSP1616 ☐ PSP1618 ☐ PSP1620
☐ PSP1716 ☐ PSP1718 ☐ PSP1720 ☐ PSP1816 ☐ PSP1818
☐ PSP1820 ☐ PSP2016 ☐ PSP2018 ☐ PSP2020

Matrix PSP Seat Cushion 21"-22" wide

HCPCS code E2608 \$595

☐ PSP2218 ☐ PSP2220 ☐ PSP2222

Matrix Vi Seat Cushion 14" - 20" wide

HCPCS code E2607 \$425

☐ KA1414-VI ☐ MA1415-VI ☐ MA1417-VI ☐ MA1616-VI ☐ MA1618-VI
☐ MA1620-VI ☐ MA1716-VI ☐ MA1718-VI ☐ MA1720-VI ☐ MA1816-VI
☐ MA1818-VI ☐ MA1820-VI ☐ MA2016-VI ☐ MA2018-VI
☐ MA2020-VI

Matrix Vi Seat Cushion 21" - 22" wide

HCPCS code E2608 \$475

☐ MA2218-VI ☐ MA2220-VI

Matrix Vi HD Seat Cushion 22" wide

HCPCS code E2608 \$625

☐ MA2222-VI

Matrix Libra Seat Cushion 14" - 20" wide

HCPCS code E2624 \$485

☐ LC1414 ☐ LC1416 ☐ LC1516 ☐ LC1518 ☐ LC1520
☐ LC1616 ☐ LC1618 ☐ LC1620 ☐ LC1716 ☐ LC1718
☐ LC1720 ☐ LC1816 ☐ LC1818 ☐ LC1820 ☐ LC1916
☐ LC1918 ☐ LC1920 ☐ LC2016 ☐ LC2018 ☐ LC2020

Matrix Libra HD Seat Cushion 21" wide

HCPCS code E2624 \$525

☐ LC2118 ☐ LC2120

Matrix Libra Seat Cushion 22" deep

HCPCS code E2624 \$595

☐ LC1622 ☐ LC1722 ☐ LC1822 ☐ LC1922 ☐ LC2022
☐ LC2122

Matrix Libra HD Seat Cushion 22" wide

HCPCS code E2625 \$595

☐ LC2218 ☐ LC2220 ☐ LC2222

MATRIX SEAT CUSHION OPTIONS

Matrix Libra Adjust Seat Cushion 14" -20" wide

HCPCS code E2624 \$625

☐ LC1414-ADJ ☐ LC1416-ADJ ☐ LC1516-ADJ ☐ LC1518-ADJ
☐ LC1520-ADJ ☐ LC1616-ADJ ☐ LC1618-ADJ ☐ LC1620-ADJ
☐ LC1716-ADJ ☐ LC1718-ADJ ☐ LC1720-ADJ ☐ LC1816-ADJ
☐ LC1818-ADJ ☐ LC1820-ADJ ☐ LC1916-ADJ ☐ LC1918-ADJ
☐ LC1920-ADJ ☐ LC2016-ADJ ☐ LC2018-ADJ ☐ LC2020-ADJ

Matrix Libra Adjust Seat Cushion 22" deep

HCPCS code E2624 \$735

☐ LC1622-ADJ ☐ LC1722-ADJ ☐ LC1822-ADJ ☐ LC1922-ADJ
☐ LC2022-ADJ ☐ LC2122-ADJ

Matrix Libra Adjust HD Seat Cushion 22" wide

HCPCS code E2625 \$735

☐ LC2118-ADJ ☐ LC2120-ADJ

Matrix Multi Cushion 14" - 20" wide w/reversable cover

HCPCS code E2607 \$445

☐ MM1414 ☐ MM1416 ☐ MM1616 ☐ MM1618 ☐ MM1620
☐ MM1716 ☐ MM1718 ☐ MM1720 ☐ MM1816 ☐ MM1818
☐ MM1820 ☐ MM2016 ☐ MM2018 ☐ MM2020

Matrix Multi Cushion 22" wide w/reversable cover

HCPCS code E2608 \$595

☐ MM2216 ☐ MM2218 ☐ MM2220 ☐ MM2222

☐ OCCC "ON-CHAIR" Cushion Credit⁽¹⁾ (\$90)

1. Not available on Matrix PS or PSP Cushions.

FABRIC OPTIONS FOR CUSHIONS

Cushion	Meshtex	Startex	Infection Control	Spacetex	Onyx	Fusion
Libra/Adjust	N/A	STD	\$50	N/A	N/A	N/A
PS/PSP	N/A	STD Reversible	STD Reversible	N/A	N/A	N/A
Vi	STD	\$50	\$50	N/A	\$100	\$150

☐ LC-SP Spacetex Fabric Option for
Libra Cushion N/C

☐ LC-IC Infection Control Fabric Upgrade for
Libra Cushion \$50

☐ FSU-MA Fusion Fabric Upgrade for Matrix Vi Cushion \$150

☐ OCU-MA Onyx Fabric Upgrade for Matrix Vi Cushion \$100

☐ STARTEX-S Smooth side Startex Fabric Upgrade for
Matrix Vi Cushion \$50

☐ STARTEX-F Fabric side Startex Fabric Upgrade for
Matrix Vi Cushion \$50

Ultra Low Maxx CG Power Positioning System

Storm Torque RX20

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

LAP BELTS

- ☐ BELT48 Push Button Style Lap Belt, Length 48"⁽¹⁾ STD
- ☐ BELT60 Push Button Style Lap Belt, Length 60"⁽¹⁾ STD
- ☐ BELT71 Push Button Style Lap Belt, Length 71"⁽¹⁾ STD
- ☐ PLB2 Bodypoint Padded Lap Belt 2 Point⁽²⁾
HCPCS code E0978. \$190
- ☐ PLB4 Bodypoint Padded Lap Belt 4 Point⁽²⁾
HCPCS code E0978. \$235
- ☐ TROBELT Safety Belt for use only with TRRO & TRRO-E N/C

1. Not available with TRRO and TRRO-E.

2. 67" length only option

MATRIX HEADREST PADS⁽¹⁾

HCPCS code E0955

Standard Pads

- ☐ ESP6 Matrix Elan Standard Pad 6" \$235
- ☐ ESP10 Matrix Elan Standard Pad 10" \$235
- ☐ ESP14 Matrix Elan Standard Pad 14" \$235
- ☐ ESP6-HUG Matrix Elan Hug Standard Pad 6" \$235
- ☐ ESP10-HUG Matrix Elan Hug Standard Pad 10" \$235
- ☐ ESP14-HUG Matrix Elan Hug Standard Pad 14" \$235
- ☐ SNUG14 Matrix SNUG Standard Pad 14" \$415
- ☐ SNUG14-A Matrix SNUG Asymmetrical Pad 14" \$415
- ☐ SNUG18 Matrix SNUG Standard Pad 18" \$415
- ☐ SNUG18-A Matrix SNUG Asymmetrical Pad 18" \$415

Standard Pads with Infection Control Cover

- ☐ ESP6-IC Matrix Elan Infection Control Pad 6" \$235
- ☐ ESP10-IC Matrix Elan Infection Control Pad 10" \$235
- ☐ ESP14-IC Matrix Elan Infection Control Pad 14" \$235
- ☐ ESP6-IC-HUG Matrix Elan Hug Infection Control Pad 6" \$235
- ☐ ESP10-IC-HUG Matrix Elan Hug Infection Control Pad 10" \$235
- ☐ ESP14-IC-HUG Matrix Elan Hug Infection Control Pad 14" \$235
- ☐ SNUG14-IC Matrix SNUG Standard
Infection Control Pad 14" \$445
- ☐ SNUG14-A-IC Matrix SNUG Asymmetrical
Infection Control Pad 14" \$445
- ☐ SNUG18-IC Matrix SNUG Standard
Infection Control Pad 18" \$445
- ☐ SNUG18-A-IC Matrix SNUG Asymmetrical
Infection Control Pad 18" \$445

Occipital /4-Point Pads

- ☐ EOP9 Matrix Elan Occipital Pad 9" \$235
- ☐ EOP12 Matrix Elan Occipital Pad 12" \$235
- ☐ E4PTSM Matrix Elan 4-Point Pad 9"W x 7"H \$325
- ☐ E4POINT Matrix Elan 4-Point Pad 11"W x 10"H \$325

Occipital / 4-Point Pads with Infection Control Cover

- ☐ EOP9-IC Matrix Elan Occipital Pad 9"
with Infection Control. \$235
- ☐ EOP12-IC Matrix Elan Occipital Pad 12"
with Infection Control. \$235
- ☐ E4PTSM-IC Matrix Elan 4-Point Pad 9"W x 7"H
with Infection Control. \$325
- ☐ E4POINT-IC Matrix Elan 4-Point Pad 11"W x 10"H
with Infection Control \$325

1. If mounting on PB Elite must order removable hardware

MOTION CONCEPTS HEADREST PAD

HCPCS code E0955

- ☐ MCSH Motion Concepts Standard \$295
- ☐ MCOH Motion Concepts Auto Style⁽¹⁾ \$350

1. Comes Standard with Fixed Mounting hardware; for use on Matrix PB or Elite Backs you must select the removable hardware upgrade below.

HEADREST MOUNTING HARDWARE

HCPCS code E1028


- ☐ FM Fixed Mounting⁽¹⁾ N/C
- ☐ LXSHW Loxx Standard Hardware 15" \$215
- ☐ LXMMH Loxx Mini Hardware 11" \$215
- ☐ MEHW Elan Hardware 15" \$220
- ☐ MEMH Elan Mini Hardware 11" \$225
- ☐ MEHW-EXT Elan Extended Hardware \$295
- ☐ FPLXSHW FLIPP Loxx Standard Hardware 15" \$515
- ☐ FPLXSHW-EX FLIPP Loxx Extended Hardware \$525
- ☐ FPLXMMH FLIPP Loxx Mini Hardware 11" \$515
- ☐ FPLXMMH-EX FLIPP Loxx Mini Extended Hardware \$525
- ☐ FPMEHW FLIPP Elan Hardware 15" \$525
- ☐ FPMMH FLIPP Elan Mini Hardware 11" \$525
- ☐ FPMEHW-EXT FLIPP Elan Extended Hardware \$605

- ☐ OHR OMIT: Headrest N/C

1. HCPCS code does not apply.

TILT ONLY BACK CANE/SPREADER BAR OPTION⁽¹⁾

- ☐ MSBSB Maxx Style Back Adjustable Spreader Bar STD
Adjust height easily on back canes
- ☐ CSBSB Classic Style Back Adjustable Spreader Bar(2) .. \$100
- ☐ OMR Optional Multi-Rail \$150
Standard Rehab Back with Classic Style Spreader Bar Only.



Ultra Low Maxx CG Power Positioning System

Storm Torque RX20

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

STANDARD BACK CANES - FOR TILT SYSTEMS

☐ SBC18	18" Straight Back Canes	STD
☐ SBC20	20" Straight Back Canes	STD
☐ SBC22	22" Straight Back Canes	STD
☐ SBC24	24" Straight Back Canes	STD

Select Cane Angle

☐ BA75 75°	☐ BA79 79°	☐ BA85 85°	☐ BA90 90°	☐ BA95 95°
☐ BA101 101°	☐ BA106 106°	☐ BA112 112°	☐ BA116 116°	

☐ ABC20	20" Angled (5°) Back Canes	N/C
☐ ABC22	22" Angled (5°) Back Canes	N/C
☐ ABC24	24" Angled (5°) Back Canes	N/C

Select Cane Angle

☐ BA81 81°	☐ BA84 84°	☐ BA91 91°	☐ BA96 96°	☐ BA101 101°
☐ BA106 106°	☐ BA111 111°	☐ BA117 117°	☐ BA121 121°	

☐ EBC22	22" Endomorph Back Canes ⁽¹⁾	\$350
☐ EBC24	24" Endomorph Back Canes ⁽¹⁾	\$350

Please Select

☐ EBCN1	1" Narrower	N/C
☐ EBCN2	2" Narrower	N/C
☐ EBCN4	4" Narrower	N/C

☐ TBC18	18" 10° Mid-Angle Tilt Back Canes ⁽²⁾	\$250
☐ TBC20	20" 10° Mid-Angle Tilt Back Canes ⁽²⁾	\$250
☐ TBC22	22" 10° Mid-Angle Tilt Back Canes ⁽²⁾	\$250
☐ TBC24	24" 10° Mid-Angle Tilt Back Canes ⁽²⁾	\$250

☐ FDBC	Fold Down Back Canes w/ Pull Pin ⁽⁶⁾	\$250
-------------------------------	---	-------

1. Comes standard with Cantilever Armrest Cane Mounted.
2. Quick Height Adjustment on Backrest and Maxx Spreader Bars not available.
3. Minimum finished back widths of 16".
4. Handle measures 6.5" high and has .75" height adjustment with angle adjustment.
5. 4" of total vertical mounting adjustment.
6. 250 lbs weight capacity. Must select canes above.

CANE OPTIONS - FOR RECLINE SYSTEMS

☐ BFS2	Recline Body Fit 2" Narrower ⁽¹⁾	\$995
☐ BFS4	Recline Body Fit 4" Narrower ⁽¹⁾	\$995
☐ AHH	Adjustable Hooking Handles ⁽²⁾ (pair)	\$450
☐ RPHS	Push Handles for Recline - Straight ⁽³⁾ (pair)	\$350
☐ RPHA	Push Handles for Recline - Angled ⁽³⁾ (pair)	\$350

1. Minimum finished back widths of 16".
2. Handle measures 6.5" high and has .75" height adjustment with angle adjustment.
3. 4" of total vertical mounting adjustment.

FINISHED BACK HEIGHT⁽¹⁾⁽²⁾ - FOR TILT ONLY/ELEVATE ONLY

☐ FBT18 18"	☐ FBT19 19"	☐ FBT20 20"	☐ FBT21 21"
☐ FBT22 22"	☐ FBT23 23"	☐ FBT24 24"	☐ FBT25 25"

1. 4" Gap from the seat pan to the bottom of back pan is recommended.
2. Finished back height is measured from seat pan to top of backrest.

FINISHED BACK HEIGHT⁽¹⁾⁽²⁾ - FOR RECLINE OR TILT/RECLINE

☐ FBR20 20"	☐ FBR21 21"	☐ FBR22 22"	☐ FBR23 23"
☐ FBR24 24"	☐ FBR25 25"	☐ FBR26 26"	☐ FBR27 27"

1. 6" Gap from the seat pan to the bottom of back pan is recommended.
2. Finished back height is measured from seat pan to top of backrest.

STANDARD BACK OPTIONS

☐ SRB	Standard Rehab Back & Cushion	N/C
☐ SS	Additional Super Soft/HR Foam ⁽¹⁾	\$450
☐ MFU	Matrx Fabric Upcharge	\$450
☐ OBC	OMIT: Back Cushion	(\$25)
☐ PIPF	Flat Planar Interface/I-Back Plate ⁽²⁾⁽³⁾⁽⁵⁾	\$125
☐ PIPR	Recessed Planar/I-Back Interface Plate ⁽²⁾⁽³⁾⁽⁵⁾⁽⁶⁾	\$125
☐ ORBC	OMIT: Rigid Back & Cushion ⁽⁴⁾	(\$75)

1. Additional foam added may impact seat depth
2. Not available on tilt only or tilt/elevate systems, for use with Recline Module.
3. Required with all Recline Systems when no back is selected. May also be used as an interface for after market backs.
4. Not available with Recline/ESR or Recline/PSB.
5. Not available with TRRO and TRRO-E.
6. Fits Comfort Company, Freedom, Pindot and Moulded Backs.

MOTION CONCEPTS SINGLE POST BACKREST⁽¹⁾⁽²⁾⁽³⁾

☐ SPBMI	Single Post Angle Adjustable Fold Down Back interface-Matrx Back.	\$750 (Folds to 25° forward)
☐ SPBCI	Single Post Angle Adjustable Fold Down Back for Contour/High Back ⁽⁴⁾	\$750 (Folds to 25° forward)

1. Single Post Back is an angle adjustable back between 73° and 109° in 6° increments from 17° forward to 19° back. The back folds down and is only available for 16"-20" seat widths & 21"-24" back height. A 4" gap from bottom of back to seat pan is required.
2. This option is not available with Recline and/or a vent tray and has a 250lb weight capacity.
3. Cane mounted armrest are not available with this option.
4. Not available with TRRO/TRRO-E

MOTION CONCEPTS REHAB SEAT SLING BACK⁽¹⁾

☐ SLINGB-TT	Sling Back with Tension Fabric	\$950
☐ SLINGB-TA	Sling Back with Adjustable Fabric	\$1,250

1. Only available on Rehab Seat System. Not available with Elevate. Cane height must be 2" taller than Finished Back Height. Not available with TRRO/TRRO-E.

FREEDOM BACK AND SEAT INTERFACE HARDWARE⁽¹⁾

Freedom Back Interface Hardware for Tilt Only

☐ FBIFEM	Flush Mount	\$125 Back can be seat width or up to 4" narrower
☐ FBITRM	Recessed Mount	\$125 Back must be 4" less than seat width

Freedom Back Interface Hardware with Recline

☐ FBIRFM	Flush Mount	\$125 Back can be seat width or up to 5" narrower
☐ FBIRRM	Recessed Mount	\$125 Back must be 5" less than seat width

☐ FSIP	Freedom Seat Interface Hardware	\$125 Seat can be up to 3" less than system width
-------------------------------	---------------------------------	--

1. Not available with TRRO and TRRO-E.

RIDE STYLE MULTI-AXIS INTERFACE⁽¹⁾

☐ RSMI	Ride Style Multi-Axis Interface for Recline Systems 14" to 20" wide	\$250
-------------------------------	--	-------

1. Only available with TRRO and TRRO-E when using with ELITE back.

Ultra Low Maxx CG Power Positioning System
Storm Torque RX20

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

UPGRADABLE BACK OPTIONS

Contour Back 15" wide - 20" wide

HCPCS code E2613 \$495

- | | | | |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| ☐ CB1512 | ☐ CB1515 | ☐ CB1517 | ☐ CB1612 |
| ☐ CB1615 | ☐ CB1617 | ☐ CB1621 | ☐ CB1712 |
| ☐ CB1715 | ☐ CB1717 | ☐ CB1721 | ☐ CB1812 |
| ☐ CB1815 | ☐ CB1817 | ☐ CB1821 | ☐ CB1912 |
| ☐ CB1915 | ☐ CB1917 | ☐ CB1921 | ☐ CB2012 |
| ☐ CB2015 | ☐ CB2017 | ☐ CB2021 | |

Contour Back 21"-22" wide

HCPCS code E2613/E2614 \$595

- | | | | |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| ☐ CB2115 | ☐ CB2117 | ☐ CB2121 | ☐ CB2215 |
| ☐ CB2217 | ☐ CB2221 | | |

High Back 16" - 22" W

HCPCS code E2620 \$595

- | | | | |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| ☐ HB1616 | ☐ HB1618 | ☐ HB1620 | ☐ HB1716 |
| ☐ HB1718 | ☐ HB1720 | ☐ HB1816 | ☐ HB1818 |
| ☐ HB1820 | ☐ HB1916 | ☐ HB1918 | ☐ HB1920 |
| ☐ HB2016 | ☐ HB2018 | ☐ HB2020 | ☐ HB2116 |
| ☐ HB2118 | ☐ HB2120 | ☐ HB2216 | ☐ HB2218 |
| ☐ HB2220 | | | |

Matrix Elite Back 14" - 20" W (3" Contour)

HCPCS code E2620 \$625

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ PBE1414 | ☐ PBE1416 | ☐ PBE1418 | ☐ PBE1420 |
| ☐ PBE1514 | ☐ PBE1516 | ☐ PBE1518 | ☐ PBE1520 |
| ☐ PBE1614 | ☐ PBE1616 | ☐ PBE1618 | ☐ PBE1620 |
| ☐ PBE1714 | ☐ PBE1716 | ☐ PBE1718 | ☐ PBE1720 |
| ☐ PBE1814 | ☐ PBE1816 | ☐ PBE1818 | ☐ PBE1820 |
| ☐ PBE1916 | ☐ PBE1918 | ☐ PBE1920 | ☐ PBE2014 |
| ☐ PBE2016 | ☐ PBE2018 | ☐ PBE2020 | |

Matrix Elite HD Back 21" - 22" W (3" Contour)

HCPCS code E2620 for 21" E2621 for 22" \$795

- | | | | |
|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|
| ☐ PBE2116-HD | ☐ PBE2120-HD | ☐ PBE2216-HD | ☐ PBE2220-HD |
|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|

Matrix Elite TR Back 15" - 20" W (3" Contour)

HCPCS code E2620 \$650

- | | | | |
|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|
| ☐ PBE1516-TR | ☐ PBE1518-TR | ☐ PBE1520-TR | ☐ PBE1616-TR |
| ☐ PBE1618-TR | ☐ PBE1620-TR | ☐ PBE1716-TR | ☐ PBE1718-TR |
| ☐ PBE1720-TR | ☐ PBE1816-TR | ☐ PBE1818-TR | ☐ PBE1820-TR |
| ☐ PBE1916-TR | ☐ PBE1918-TR | ☐ PBE1920-TR | ☐ PBE2016-TR |
| ☐ PBE2018-TR | ☐ PBE2020-TR | | |

Matrix Elite TR HD Back 21" - 22" W (3" Contour)

HCPCS code E2620 for 21" E2621 for 22" \$845

- | | |
|---------------------------------------|---------------------------------------|
| ☐ PBE2120-TRHD | ☐ PBE2220-TRHD |
|---------------------------------------|---------------------------------------|

UPGRADABLE BACK OPTIONS

Matrix Elite Deep Back 14" - 20" W (6" Contour)

HCPCS code E2620 \$760

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ EDB1416 | ☐ EDB1418 | ☐ EDB1420 | ☐ EDB1514 |
| ☐ EDB1516 | ☐ EDB1518 | ☐ EDB1520 | ☐ EDB1614 |
| ☐ EDB1616 | ☐ EDB1618 | ☐ EDB1620 | ☐ EDB1714 |
| ☐ EDB1716 | ☐ EDB1718 | ☐ EDB1720 | ☐ EDB1814 |
| ☐ EDB1816 | ☐ EDB1818 | ☐ EDB1820 | ☐ EDB1914 |
| ☐ EDB1916 | ☐ EDB1918 | ☐ EDB1920 | ☐ EDB2014 |
| ☐ EDB2016 | ☐ EDB2018 | ☐ EDB2020 | |

Matrix Elite Deep HD Back 21" - 22" W (6" Contour)

HCPCS code E2620 for 21" E2621 for 22" \$930

- | | | |
|-------------------------------------|-------------------------------------|-------------------------------------|
| ☐ EDB2116-HD | ☐ EDB2118-HD | ☐ EDB2120-HD |
| ☐ EDB2216-HD | ☐ EDB2218-HD | ☐ EDB2220-HD |

Matrix Elite E2S Back 14" - 20" W x 14"/16"H (3" Contour)

HCPCS code E2615 \$625

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ E2S1414 | ☐ E2S1416 | ☐ E2S1514 | ☐ E2S1516 |
| ☐ E2S1614 | ☐ E2S1616 | ☐ E2S1714 | ☐ E2S1716 |
| ☐ E2S1814 | ☐ E2S1816 | ☐ E2S1914 | ☐ E2S1916 |
| ☐ E2S2014 | ☐ E2S2016 | | |

Matrix Elite E2S Back 14" - 20" W x 18"/20"H (3" Contour)

HCPCS code E2615 \$645

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ E2S1418 | ☐ E2S1420 | ☐ E2S1518 | ☐ E2S1520 |
| ☐ E2S1618 | ☐ E2S1620 | ☐ E2S1718 | ☐ E2S1720 |
| ☐ E2S1818 | ☐ E2S1820 | ☐ E2S1918 | ☐ E2S1920 |
| ☐ E2S2018 | ☐ E2S2020 | | |

Matrix Elite E2 Lateral Back

14" - 20" W x 14"/16"H (3" Contour) HCPCS code E2620 \$725

- | | | | |
|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|
| ☐ E2LB1414 | ☐ E2LB1416 | ☐ E2LB1514 | ☐ E2LB1516 |
| ☐ E2LB1614 | ☐ E2LB1616 | ☐ E2LB1714 | ☐ E2LB1716 |
| ☐ E2LB1814 | ☐ E2LB1816 | ☐ E2LB1914 | ☐ E2LB1916 |
| ☐ E2LB2014 | ☐ E2LB2016 | | |

Matrix Elite E2 Lateral Back

14" - 20" W x 18"/20"H (3" Contour) HCPCS code E2620 \$745

- | | | | |
|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|
| ☐ E2LB1418 | ☐ E2LB1420 | ☐ E2LB1518 | ☐ E2LB1520 |
| ☐ E2LB1618 | ☐ E2LB1620 | ☐ E2LB1718 | ☐ E2LB1720 |
| ☐ E2LB1818 | ☐ E2LB1820 | ☐ E2LB1918 | ☐ E2LB1920 |
| ☐ E2LB2018 | ☐ E2LB2020 | | |

Matrix Elite E2 Deep Back

14" - 20" W x 14"/16"H (6" Contour) HCPCS code E2620 \$755

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ E2D1414 | ☐ E2D1416 | ☐ E2D1514 | ☐ E2D1516 |
| ☐ E2D1614 | ☐ E2D1616 | ☐ E2D1714 | ☐ E2D1716 |
| ☐ E2D1814 | ☐ E2D1816 | ☐ E2D1914 | ☐ E2D1916 |
| ☐ E2D2014 | ☐ E2D2016 | | |

Ultra Low Maxx CG Power Positioning System

Storm Torque RX20

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

UPGRADABLE BACK OPTIONS

Matrix Elite E2 Deep Back
14" - 20" W x 18"/20"H (6" Contour) HCPCS code E2620.\$755

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ E2D1418 | ☐ E2D1420 | ☐ E2D1518 | ☐ E2D1520 |
| ☐ E2D1618 | ☐ E2D1620 | ☐ E2D1718 | ☐ E2D1720 |
| ☐ E2D1818 | ☐ E2D1820 | ☐ E2D1918 | ☐ E2D1920 |
| ☐ E2D2018 | ☐ E2D2020 | | |

Matrix PB Posture Back
14" - 20" W (5" Contour) HCPCS code E2615.\$625

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ MPB1412 | ☐ MPB1612 | ☐ MPB1616 | ☐ MPB1620 |
| ☐ MPB1712 | ☐ MPB1716 | ☐ MPB1720 | ☐ MPB1812 |
| ☐ MPB1816 | ☐ MPB1820 | ☐ MPB2012 | ☐ MPB2016 |
| ☐ MPB2020 | | | |

Matrix PB Heavy Duty Back⁽¹⁾
20" - 22" W (5" Contour) HCPCS code E2615 for 21" E2616 for 22" \$835

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ HDB2016 | ☐ HDB2020 | ☐ HDB2216 | ☐ HDB2220 |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|

Matrix PB Posture Back Deep⁽¹⁾
14" - 20" W (7" Contour) HCPCS code E2615 \$645

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ PBD1412 | ☐ PBD1612 | ☐ PBD1616 | ☐ PBD1620 |
| ☐ PBD1712 | ☐ PBD1716 | ☐ PBD1720 | ☐ PBD1812 |
| ☐ PBD1816 | ☐ PBD1820 | ☐ PBD2012 | ☐ PBD2016 |
| ☐ PBD2020 | | | |

Matrix PB Heavy Duty Deep Back⁽¹⁾
20"-22" W (5" Contour)\$935

- | | | | |
|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|
| ☐ DHDB2016 | ☐ DHDB2020 | ☐ DHDB2216 | ☐ DHDB2220 |
|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|

Matrix MAC Back⁽²⁾ HCPCS code E2620.\$955

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ MAC1616 | ☐ MAC1618 | ☐ MAC1620 | ☐ MAC1622 |
| ☐ MAC1716 | ☐ MAC1718 | ☐ MAC1720 | ☐ MAC1722 |
| ☐ MAC1816 | ☐ MAC1818 | ☐ MAC1820 | ☐ MAC1822 |
| ☐ MAC1916 | ☐ MAC1918 | ☐ MAC1920 | ☐ MAC1922 |
| ☐ MAC2016 | ☐ MAC2018 | ☐ MAC2020 | ☐ MAC2022 |

Matrix MAC Back w/ HUG⁽²⁾ HCPCS code E2620. \$1,055

- | | | |
|--------------------------------------|--------------------------------------|--------------------------------------|
| ☐ MAC1616-HUG | ☐ MAC1618-HUG | ☐ MAC1620-HUG |
| ☐ MAC1622-HUG | ☐ MAC1716-HUG | ☐ MAC1718-HUG |
| ☐ MAC1720-HUG | ☐ MAC1722-HUG | ☐ MAC1816-HUG |
| ☐ MAC1818-HUG | ☐ MAC1820-HUG | ☐ MAC1822-HUG |
| ☐ MAC1916-HUG | ☐ MAC1918-HUG | ☐ MAC1920-HUG |
| ☐ MAC1922-HUG | ☐ MAC2016-HUG | ☐ MAC2018-HUG |
| ☐ MAC2020-HUG | ☐ MAC2022-HUG | |

Matrix MAC Lateral Contour

- | | | | |
|---------------------------------|-----------|---------------------------------|------------|
| ☐ MLCSL | Left 3.5" | ☐ MLCSR | Right 3.5" |
| ☐ MLCDL | Left 6" | ☐ MLCDR | Right 6" |
| ☐ MLCXDL | Left 9" | ☐ MLCXDR | Right 9" |

- ☐ OCBC "ON-CHAIR" Matrix Back Credit(\$125)

1. The Heavy Duty PB Back and PB Deep Back are not available on Recline systems with seat widths of 21" and 22"

2. Must select Pad Contour in Matrix MAC Lateral Contour section.

OTHER UPGRADABLE BACK OPTIONS

Cushion	Meshtex	4-Way Stretch Polyester	Infection Control	Onyx	Fusion
Contour	N/A	STD	N/A	\$100	N/A
High	STD	N/A	N/A	\$100	\$100
E2, Elite, Elite HD, Elite Deep	STD	N/A	\$100	N/A	N/A
E2, Elite TR, Elite TR HD	STD	N/A	\$100	N/A	N/A

- | | | |
|----------------------------------|--|-------|
| ☐ EBSS | Super Soft Foam for Elite Back Upgrade. | \$100 |
| ☐ FFUHB | Fusion Fabric Upgrade for High Backs. | \$100 |
| ☐ OCUCB | Onyx Fabric Upgrade for Contour Backs. | \$100 |
| ☐ OCUHB | Onyx Fabric Upgrade for High Backs. | \$100 |
| ☐ CHIEB | Chest Harness Interface Kit for Elite Back ⁽¹⁾⁽²⁾ | \$95 |
| ☐ CHIPB | Chest Harness Interface Kit for PB Back ⁽¹⁾⁽²⁾ | \$45 |
| ☐ E2CHI-S | Chest Harness Interface Kit for E2 Back for 14"-15" Backs ⁽¹⁾⁽²⁾ | \$85 |
| ☐ E2CHI-M | Chest Harness Interface Kit for E2 Back for 16"-20" Backs ⁽¹⁾⁽²⁾ | \$85 |
| ☐ E2BKB | Contoured anatomical PSIS/Sacrolumbar Support. | \$55 |

1. Does not include the chest strap.

2. This Interface Kit is only Required for 4 point harness

CHEST STRAPS - ABDOMINAL SUPPORT HCPCS code E0960

- | | | |
|----------------------------------|---|-------|
| ☐ MFUS | BP Monoflex, Underarm RLS S(3" x 15.5") | \$195 |
| ☐ MFUM | BP Monoflex, Underarm RLS M(4" x 18.25") | \$195 |
| ☐ MFUL | BP Monoflex, Underarm RLS L(4.5" x 21.5") | \$195 |
| ☐ MFCS | BP Monoflex, Center RLS S(2.5" x 15.5") | \$195 |
| ☐ MFCM | BP Monoflex, Center RLS M(3.25" x 18.25") | \$195 |
| ☐ MFCL | BP Monoflex, Center RLS L(3.5" x 21.5") | \$195 |
| ☐ PFLEX-M | PostureFlex Abdominal Support Height 6.5", length 16" to 23" ⁽¹⁾ | \$245 |
| ☐ PFLEX-L | PostureFlex Abdominal Support Height 7", length 20" to 32" ⁽¹⁾ | \$245 |
| ☐ QRTAB | Quad release tab for easy independent application and removal of PostureFlex | \$35 |

1. Only available on the E2 and E2 Deep Backs.

Ultra Low Maxx CG Power Positioning System

Storm Torque RX20

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

TILT ARMRESTS

Adult Dual Post

☐ DP1	9.75" - 12.75"	STD
☐ DP2	12.75" - 15.75"	STD
☐ DP3	15.75" - 18.75"	STD
☐ FP	Flat Pouch	STD
☐ GB	Glove Boxes	N/C
☐ OAL	Outback Arm Option Left	\$195
☐ OAR	Outback Arm Option Right	\$195

Ultra Rail Seat Mount Flip Back Cantilever MAXX Tilt Arm⁽¹⁾

☐ RMCAN	9" - 12.5"	\$488
☐ RMCAN2	11" - 14"	\$488
☐ DPM	Dual Post Module to RMCAN (pair)	\$250

Outback Cantilever Tilt Armrest

☐ OCTA1	9"-12"	\$488
☐ OCTA2	12"-15"	\$488
☐ OCTA2	6"-12"	\$488

Cantilver Arms for Endomorph Back Canes⁽¹⁾

☐ CAFEB	See page 4	STD
--------------------------------	------------	-----

1. Not available with Small Adult Sizes.

RECLINE ARMRESTS

Reclining Adjustable Height

☐ RC1	9.5"-13"	STD
☐ RC2	13"-16"	STD
☐ OAL	Outback Arm Option Left	\$195
☐ OAR	Outback Arm Option Right	\$195

Maxx Style Cane Mounted Height, Width, Angle Adj.

Cantilever Flip Back Armrests⁽¹⁾

☐ MCAN1	9.5"-13"	\$488
☐ MCAN2	12.5"-16"	\$488

Outback Cantilever Recline Armrest

☐ OCRA1	9"-12"	\$488
☐ OCRA2	12"-15"	\$488
☐ OCRA3	7" - 9"	\$488

Cane Mounted Height and Angle Adj.

Cantilver Flip Back Armrest⁽¹⁾

☐ CANQ1	9.5"-12"	\$488
☐ CANQ2	10"-13"	\$488
☐ CANQ3	13"-16"	\$488

Adult Dual Post On Recline

☐ DPOR1	9.75" - 12.75"	\$250
☐ DPOR2	12.75" - 15.75"	\$250
☐ DPOR3	15.75" - 18.75"	\$250
☐ FP	Flat Pouch	STD
☐ GB	Glove Boxes	N/C

1. Not available with Small Adult Sizes.

ARMREST TUBE LENGTH - PLEASE SPECIFY

☐ ATFL	Left Full	STD
☐ ATDL	Left Desk	STD
☐ ATFR	Right Full	STD
☐ ATDR	Right Desk	STD

ARMPADS⁽⁵⁾

Standard

☐ SAPFL	Left Full (2.25" x 14")	STD
☐ SAPDL	Left Desk (2.25" x 10")	STD
☐ SAPML	Left Mini (2.25" x 7")	STD
☐ SAPFR	Right Full (2.25" x 14")	STD
☐ SAPDR	Right Desk (2.25" x 10")	STD
☐ SAPMR	Right Mini (2.25" x 7")	STD

Waterfall Arm Pads⁽¹⁾

☐ WAPFL	Left Full (13")	\$50
☐ WAPDL	Left Desk (9")	\$50
☐ WAPFR	Right Full (13")	\$50
☐ WAPDR	Right Desk (9")	\$50

Modular Arm Pad - Left⁽²⁾

☐ MAPTFL	Left Tray Full (14")	\$55
☐ MAPTDL	Left Tray Desk (10")	\$55

Select Pad Insert below:

☐ ISSFL	Iskin - Integrated soft skin foam
☐ SCVFL	Startex Covered Visco foam

Modular Arm Pad - Right⁽²⁾

☐ MAPTFR	Right Tray Full (14")	\$55
☐ MAPTDR	Right Tray Desk (10")	\$55

Select Pad Insert below:

☐ ISSFR	Iskin - Integrated soft skin foam
☐ SCVFR	Startex Covered Visco foam

Flat Multi Position Arm Pads

☐ MPAFL	Left Full (14")	\$55
☐ MPADL	Left Desk (10")	\$55
☐ MPAML	Left Mini (7")	\$55
☐ MPAFR	Right Full (14")	\$55
☐ MPADR	Right Desk (10")	\$55
☐ MPAMR	Right Mini (7")	\$55

Gel Pads

☐ GAPNL	Narrow Left Desk (2" x 12")	\$75
☐ GAPNR	Narrow Right Desk (2" x 12")	\$75
☐ GAPWL	Wide Left Desk (3.5" x 12")	\$75
☐ GAPWR	Wide Right Desk (3.5" x 12")	\$75
☐ GAPLL	Long Left Full (3.5" x 14")	\$90
☐ GAPLR	Long Right Full (3.5" x 14")	\$90

Arm Troughs

☐ EATL	Ergonomic Arm Trough - Left HCPCS code E2209	\$85
☐ EATR	Ergonomic Arm Trough - Right HCPCS code E2209	\$85
☐ FHPL	Flat Hand Pad - Left ⁽³⁾ HCPCS code K0108	\$55
☐ FHPR	Flat Hand Pad - Right ⁽³⁾ HCPCS code K0108	\$55

Arm Pad Storage/Cup Holder⁽⁴⁾

☐ APSL	Left	\$295
☐ APSR	Right	\$295
☐ OAP	Omit Armpads(7/8" diameter tube)	(\$25)

1. Full length pad not compatible with the desk length armrest tube if using Ultra Rail mounted Tilt Arms.
2. Requires interchangeable pad insert.
3. Must order with ergonomic arm trough.
4. Not available on the joystick side, or with Outback Cantilever arms.
5. Full length arm pads used on the same side as the Joystick result in a "Back to Joystick" measurement of 20" minimum (24" including switch mounting)

Ultra Low Maxx CG Power Positioning System

Storm Torque RX20

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

ARM ADDUCTOR PAD⁽¹⁾⁽²⁾

☐ AADSL	Left Small (6" x 2.25")	\$100
☐ AADSR	Right Small (6" x 2.25")	\$100
☐ AADLL	Left Large (8.5" x 2.25")	\$100
☐ AADLR	Right Large (8.5" x 2.25")	\$100

1. Only available with Modular or Flat Multi Position Arm Pads.
2. 1" thick pad, unlimited fore/aft adjustment.

ELBOW BLOCKS

Left Elbow Block⁽¹⁾

Select Pad Size and Both Extensions:

☐ EBXL	Size XS (3.25" x 4.25")	\$195
☐ EBSL	Size S (3.5" x 5.25")	\$195
☐ EBML	Size M (3.75" x 5.75")	\$195
☐ EBL	Size L (4.25" x 6.5")	\$195

☐ LV3	Left Vertical Extension 3"
☐ LV5	Left Vertical Extension 5"
☐ LH5	Left Horizontal Extension 5"
☐ LH7	Left Horizontal Extension 7"

Right Elbow Block⁽¹⁾

Select Pad Size and Both Extensions:

☐ EBXR	Size XS (3.25" x 4.25")	\$195
☐ EBSR	Size S (3.5" x 5.25")	\$195
☐ EBMR	Size M (3.75" x 5.75")	\$195
☐ EBLR	Size L (4.25" x 6.5")	\$195

☐ RV3	Right Vertical Extension 3"
☐ RV5	Right Vertical Extension 5"
☐ RH5	Right Horizontal Extension 5"
☐ RH7	Right Horizontal Extension 7"

1. Prevents arm from slipping.

ADDITIONAL ARM ACCESSORIES

☐ MACESL	Multi Axis Upper Extremity Support ⁽¹⁾ (each) Left Side	\$225
☐ MACESR	Multi Axis Upper Extremity Support ⁽¹⁾ (each) Right Side	\$225
☐ OAS2L	1" Offset Armrest Spacer Left(each)	\$125
☐ OAS2R	1" Offset Armrest Spacer Right(each)	\$125
☐ OAS1L	1/2" Offset Armrest Spacer Left(each)	\$125
☐ OAS1R	1/2" Offset Armrest Spacer Right(each)	\$125
☐ MPAPA	Adapter Block ⁽²⁾ (each) Quantity ()	\$35

1. Allows rotation of the armpad, only compatible with Gel Arm pad, Flat Arm pad or Ergonomic Arm Troughs.
2. Not available with Standard Arm Pads or Waterfall Arm Pads

MAXX LATERAL TRUNK SUPPORTS⁽¹⁾

Lateral Trunk Support Left HCPCS code E0956

☐ LATXL	Size XS (3.25" x 4.25")	\$100
☐ LATSL	Size S (3.5" x 5.25")	\$100
☐ LATML	Size M (3.75" x 5.75")	\$100
☐ LATLL	Size L (4.25" x 6.5")	\$100

Add:

☐ SALATL	MAXX Style Swing Away Hardware - Left ⁽²⁾ HCPCS code E1028	\$195
☐ ATLL	Additional Telescoping Link - Left ⁽³⁾	\$75

Lateral Trunk Support Right HCPCS code E0956

☐ LATXR	Size XS (3.25" x 4.25")	\$100
☐ LATSR	Size S (3.5" x 5.25")	\$100
☐ LATMR	Size M (3.75" x 5.75")	\$100
☐ LATLR	Size L (4.25" x 6.5")	\$100

Add:

☐ SALATR	MAXX Style Swing Away Hardware - Right ⁽²⁾ HCPCS code E1028	\$195
☐ ATLR	Additional Telescoping Link - Right ⁽³⁾	\$75

1. Comes standard with Fixed Mount Hardware. Not available with Matrix PB, PB Elite Deep, and PB Deep Backs.
2. Comes standard with 1" telescoping link which gives 0-2.75" medial offset.
3. The amount of offset available may be impacted by back thickness.

MATRIX FIXED LATERAL TRUNK SUPPORTS⁽¹⁾

HCPCS code E1028

☐ E2FLL	Matrix E2 Fixed Laterals-Left(6.5x5") ⁽²⁾	\$145
☐ E2FLR	Matrix E2 Fixed Laterals-Right(6.5x5") ⁽²⁾	\$145
☐ PBSSL	Matrix Small Lateral Fixed-Left (3x5.5")	\$100
☐ PBSSR	Matrix Small Lateral Fixed-Right (3x5.5")	\$100
☐ PBLSL	Matrix Standard Lateral Fixed-Left (7.5x5")	\$100
☐ PBLSR	Matrix Standard Lateral Fixed-Right (7.5x5")	\$100
☐ PBOLL	Matrix Offset Lateral Fixed Left (7.5x5")	\$110
☐ PBOLR	Matrix Offset Lateral Fixed Right (7.5x5")	\$110

1. Not available with Matrix PB, PB Elite Deep, and PB Deep Backs.
2. Only available with E2 backs.

MATRIX SWING AWAY LATERAL TRUNK SUPPORTS⁽¹⁾

Matrix Lateral Trunk Support Left HCPCS code E0956

☐ EPADXL	Size XS (3.25" x 4.25")	\$100
☐ EPADSL	Size S (3.5" x 5.25")	\$100
☐ EPADML	Size M (3.75" x 5.75")	\$100
☐ EPADLL	Size L (4.25" x 6.5")	\$100

Must select one option:

☐ SAHEBL	Matrix Elite Swing Away Hardware Left HCPCS code E1028	\$195
☐ OSAHL	Matrix Offset Elite Swing Away Hardware Left HCPCS code E1028	\$195

Matrix Lateral Trunk Support Right HCPCS code E0956

☐ EPADXR	Size XS (3.25" x 4.25")	\$100
☐ EPADSR	Size S (3.5" x 5.25")	\$100
☐ EPADMR	Size M (3.75" x 5.75")	\$100
☐ EPADLR	Size L (4.25" x 6.5")	\$100

Must select one option:

☐ SAHEBR	Matrix Elite Swing Away Hardware Right HCPCS code E1028	\$195
☐ OSAHR	Matrix Offset Elite Swing Away Hardware Right HCPCS code E1028	\$195

1. Not available with Matrix PB, PB Elite Deep, and PB Deep Backs.

Ultra Low Maxx CG Power Positioning System

Storm Torque RX20

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

HIP AND KNEE SUPPORTS

Lateral Hip Support Left HCPCS code E0953

- ☐ HSPXLH Size XS (5.5" W x 4" H) \$120
- ☐ HSPSMLH Size S/M (7" W x 4" H) \$120
- ☐ HSPSLH Size S (7" W x 6" H) \$120
- ☐ HSPMLH Size M (9" W x 4" H) \$120
- ☐ HSPLLH Size L (13" W x 4" H) \$120

Lateral Hip Support Right⁽¹⁾ HCPCS code E0953

- ☐ HSPXRH Size XS (5.5" W x 4" H) \$120
- ☐ HSPSMRH Size S/M (7" W x 4" H) \$120
- ☐ HSPSRH Size S (7" W x 6" H) \$120
- ☐ HSPMRH Size M (9" W x 4" H) \$120
- ☐ HSPLRH Size L (13" W x 4" H) \$120

Lateral Knee Support Left⁽¹⁾ HCPCS code E0953

- ☐ HSPXLK Size XS (5.5" W x 4" H) \$120
- ☐ HSPSMLK Size S/M (7" W x 4" H) \$120
- ☐ HSPSLK Size S (7" W x 6" H) \$120
- ☐ HSPMLK Size M (9" W x 4" H) \$120
- ☐ HSPLLK Size L (13" W x 4" H) \$120

Lateral Knee Support Right⁽¹⁾ HCPCS code E0953

- ☐ HSPXRK Size XS (5.5" W x 4" H) \$120
- ☐ HSPSMRK Size S/M (7" W x 4" H) \$120
- ☐ HSPSRK Size S (7" W x 6" H) \$120
- ☐ HSPMRK Size M (9" W x 4" H) \$120
- ☐ HSPLRK Size L (13" W x 4" H) \$120

Add Removable Hardware:

MAXX Style Quick Release, Multi-Axis Mounting

Hardware for Hip Supports HCPCS code E1028

- ☐ MHSHLH Left mounted at Hip \$215
- ☐ MHSHLK Left mounted at Knee \$215
- ☐ MHSHRH Right mounted at Hip \$215
- ☐ MHSHRK Right mounted at Knee \$215

Swing Away Hardware for Hip Supports HCPCS code E1028

- ☐ SHSHLH Left mounted at Hip \$215
- ☐ SHSHLK Left mounted at Knee \$215
- ☐ SHSHRH Right mounted at Hip \$215
- ☐ SHSHRK Right mounted at Knee \$215

Extended Maxx Style Quick Release, Multi-Axis Mounting

Hardware for Hip Supports⁽²⁾ HCPCS code E1028

- ☐ EMQHLH Left mounted at Hip \$215
- ☐ EMQHLK Left mounted at Knee \$215
- ☐ EMQHRH Right mounted at Hip \$215
- ☐ EMQHRK Right mounted at Knee \$215

Lift Off Hardware for Hip Support

(Hip Position Only) HCPCS code E1028

- ☐ LOHL Left \$195
- ☐ LOHR Right \$195

1. Comes standard with Fixed Mount Hardware.

2. Recommended when finished seat cushion thickness will be greater than 3.5" high.

POWER CENTER MOUNT FOOT PLATFORM

MUST BE COMPLETED	<i>Measurement in Inches</i>
<i>Knee to Heel:</i>	
<i>Cushion Thickness: (Subtract)</i>	
<i>Cushion Compression: (Add)</i>	
<i>Seat Pan to Footplate: (Total)</i>	

LNx Power Center Mount Foot Platform HCPCS code E1012

- ☐ LNX 9" Articulation⁽¹⁾⁽²⁾ \$2,995
 - ☐ LSPTF1 10" - 14" Seat pan to footplate
 - ☐ LSPTF2 13" - 19" Seat pan to footplate
- ☐ LNXA 6" Articulation⁽¹⁾⁽²⁾⁽⁶⁾ \$2,995
 - ☒ LASPTF2 10" - 14" Seat pan to footplate

Options:

- ☐ LIFPN Rubber coated footplate 10.5" W x 10" D STD
- ☐ LICP Individual calf pads STD
- ☐ OPCP One Piece Malleable Calf Panel⁽³⁾ \$250
- ☐ LNXPF Planter Flex Option for LNX Foot platform⁽⁴⁾ \$250

Latitude Style Foot Platform

- ☐ LSFP Latitude Style Foot Platform⁽⁵⁾ \$2,995
For finished knee to heel measurements of 6.5"-12"
- ☐ LSPFP Latitude Style Padded Foot Plate \$250

1. 3" ground clearance required.

2. Only available on adult range seat systems 16" deep or greater.

For shorter seat depths call Customer Service for other options.

3. Not available with Planter Flex Option for LNX Foot platform, or with LNXA.

4. Only available with Foot Platform Options.

5. Not available with TRRO and TRRO-E.

6. For clients under 4' only.

MANUAL ELEVATING CENTER MOUNT LEGREST

- ☐ MECML Manual Elevating Center Mount Legrest
With Planter Flex \$595

Select legrest set up position:

- ☐ 90SU Set at 90° (elevating range 90-40°) STD

Rubber Coated Footplate:

- ☐ MRCPN Rubber Coated Footplate 10.5" W x 10" D \$165

Select seat pan to footplate for One Piece Foot Plate

- ☐ CSPTF1 11.5" - 15" N/C
- ☐ CSPTF2 14" - 19" N/C

Individual Footplates:


- ☐ IFSL Small - (5" x 7.5" deep) - Left \$165
- ☐ IFSR Small - (5" x 7.5" deep) - Right \$165
- ☐ IFML Medium - (5.5" x 9.5" deep) - Left \$165
- ☐ IFMR Medium - (5.5" x 9.5" deep) - Right \$165
- ☐ IFLL Large - (6" x 11.5" deep) - Left \$165
- ☐ IFLR Large - (6" x 11.5" deep) - Right \$165

Select seat pan to footplate for Individual Footplates

- ☐ CSPTF3 9.5" - 14.5" N/C
- ☐ CSPTF4 12" - 19" N/C

Optional:

- ☐ CICP Individual calf pads \$200



Ultra Low Maxx CG Power Positioning System

Storm Torque RX20

Dealer Name: _____

Dealer Account Number: _____

PO #: _____


Tag: _____

STANDARD FRONT RIGGING

- ☐ FCMP Fixed Center Mount Platform⁽¹⁾N/C
☐ FIFPN Comes with a rubber coated footplate
 10.5" W x 10" D.N/C

Select seat pan to footplate:

- ☐ FSPTF1 9" - 13"N/C
☐ FSPTF14 13" -19"N/C



- ☐ I70SAF Invacare 70° Swingaway Footrests (pair)N/C
☐ I70-CF Comes with composite foot plates.N/C

1.3" ground clearance required


MANUAL CENTER MOUNT FRONT RIGGING

- ☐ MSFCMFP Maxx Style Fixed Center Mount
 Foot Platform⁽¹⁾⁽²⁾\$395

- ☐ MIFPN Comes with a rubber coated footplate
 10.5" W x 10" D.N/C
☐ MICP Comes standard with individual calf padsN/C

Select seat pan to footplate:


- ☐ MSPTF1 9" - 13"N/C
☐ MSPTF14 13" -19"N/C



1. Depth adjustable, and easy slide adjustment for seat pan to foot plate extension.

2. 3" ground clearance required.





FOOT PLATFORM OPTIONS FOR FIXED/LNX CMFP



- ☐ SFFP Small Foot Platform (12"w x 8.5"d)\$180
☐ MFFP Medium Platform (12"w x 10.75"d)\$180
☐ LFFP Large Platform (14.5"w x 13"d)\$250
☐ XLFFP Extra Large Platform (17"w x 13"d)\$250
☐ XXLFFP Extra Extra Large Platform (21"w x 13"d)\$250

1. May need to adjust leg position to avoid interference with front casters

INDIVIDUAL FOOT PLATES FOR FIXED/LNX CMFP



- ☐ IFSL Small (5" W x 7.5" D) - Left\$165
☐ IFSR Small (5" W x 7.5" D) - Right\$165
☐ IFML Medium (5.5" W x 9.5" D) - Left\$165
☐ IFMR Medium (5.5" W x 9.5" D) - Right\$165
☐ IFLL Large (6" W x 11.5" D) - Left\$165
☐ IFLR Large (6" W x 11.5" D) - Right\$165

Add:

- ☐ TSSL Left Toe Strap Small 6.5"\$125
☐ TSSR Right Toe Strap Small 6.5"\$125
☐ TSML Left Toe Strap Medium 8.5"\$125
☐ TSMR Right Toe Strap Medium 8.5"\$125
☐ TSLL Left Toe Strap Large 10.5"\$125
☐ TSLR Right Toe Strap Large 10.5"\$125

- ☐ MFOL Multi-Axis Footplate (6"W x 10"D) - Left\$225
☐ MFOR Multi-Axis Footplate (6"W x 10"D) - Right\$225

1. Foot plates may interfere with casters in certain configurations

MAXX INDEPENDENT ADJUSTABLE POSITIONING



- 5.5" - 6" WIDTH ADJUSTMENT
- 9.5" - 11.5" DEPTH ADJUSTMENT
- UP TO 9° ADDUCTION
- UP TO 45° ABDUCTION

- ☐ MPFPL Maxx Independent Adjustable Positioning
 Foot plate⁽¹⁾ - Left\$200
☐ MPFPR Maxx Independent Adjustable Positioning
 Foot plate⁽¹⁾ - Right\$200
☐ MTSL Left Toe Strap (each) HCPCS code E0952\$35
☐ MTSR Right Toe Strap (each) HCPCS code E0952\$35

1. Foot plates may interfere with casters in certain configurations.

Ultra Low Maxx CG Power Positioning System

Storm Torque RX20

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

INVACARE CENTER MOUNT⁽¹⁾

Invacare Center Mount with Footboard⁽³⁾

Select seat pan to footplate:

- ☐ CTMT-FJ-FB Invacare Center Mount 9.25" - 12" \$760
- ☐ CTMT-F-FB Invacare Center Mount 12.25" - 15" \$760

Invacare Center Mount with Individual Footplates

Select seat pan to footplate:

- ☐ CTMT-FJ Invacare Center Mount 9.25" - 12" \$595
- ☐ CTMT-F Invacare Center Mount 12.25" - 15" \$595

Select Footplates:

- ☐ SM Individual Foot Plates (5.25"W x 7.5"D) \$165
- ☐ MDM Individual Foot Plates (5.25"W x 10"D) \$165
- ☐ LRG Individual Foot Plates (5.9"W x 12"D) \$165

Optional:

- ☐ CALF Invacare Calf Pads for use with CTMT⁽²⁾ \$268

1. 3" ground clearance required.

2. Not available with CTMT-FJ, CTMT-FJ-FB

3. Footboard Included

MAXX STYLE POWER ELEVATING LEGRESTS

- ☐ MSPPL MAXX Swing Away
Power Elevating⁽¹⁾ HCPCS code E1010 \$2,490

Select Power Legrest Function:

- ☐ INDA Individual Movement N/C
- ☐ COMA Combined Movement N/C
- ☐ SYNA Synchronized/Recline Movement N/C

Select Receivers:

- ☐ MRS Straight Receivers N/C
- ☐ MRF Flared 1" Outward Receivers N/C
- ☐ MRF2 Flared 2" Outward Receivers \$250

Select Seat Pan to Footplate:

- ☐ PSPTF5 Length: 10"-14" N/C
- ☐ PSPTF2 Length: 13"-17" N/C
- ☐ PSPTF7 Length: 16"-20" N/C

Optional:

- ☐ PMCP Moulded I-Skin Calf Pads STD
- ☐ PFCPMS Flat Calf Pad Option (pair) \$275
- ☐ PCCPMS Curved Calf Pad Option⁽²⁾ (pair) \$275

1. Must choose footplate on page 15.

2. Not available on 14 and 15" seat widths. For 16" and 17" seat widths you must select flared receivers.

MAXX STYLE MANUAL ELEVATING LEGRESTS

- ☐ MMPPSA MAXX Swing Away
Manual Elevating⁽¹⁾ HCPCS code E0990 \$680

Select Receivers:

- ☐ ERS Straight Receivers N/C
- ☐ ERF Flared 1" Outward Receivers N/C
- ☐ ERF2 Flared 2" Outward Receivers \$250

Select Seat Pan to Footplate:

- ☐ ESPTF5 Length: 10"-14" N/C
- ☐ ESPTF2 Length: 13"-17" N/C
- ☐ ESPTF7 Length: 16"-20" N/C

Optional:

- ☐ EMCP Moulded I-Skin Calf Pads STD
- ☐ EFCPMS Flat Calf Pad Option (pair) \$275
- ☐ ECCPMS Curved Calf Pad Option⁽²⁾ (pair) \$275

1. Must choose footplate on page 15.

2. Not available on 14 and 15" seat widths. For 16" and 17" seat widths you must select flared receivers.

SWING AWAY LEGREST

- ☐ HD70 Heavy Duty 70° Swing Away Footrest(pair)⁽¹⁾⁽²⁾ . \$495

Select Receivers:

- ☐ RS Straight Receivers N/C
- ☐ RF Flared 1" Outward Receivers N/C
- ☐ RF2 Flared 2" Outward Receivers \$250

Select Seat Pan to Footplate:

- ☐ HSPTF11 Length: 6"-10" N/C
- ☐ HSPTF8 Length: 10"-14" N/C
- ☐ HSPTF2 Length: 12"-16" N/C
- ☐ HSPTF7 Length: 15"-19" N/C

Options:

- ☐ HD70FCP Flat Calf Pad Option (pair) \$275
- ☐ HD70CCP Curved Calf Pad Option⁽³⁾ (pair) \$275

- ☐ I70SAFT Invacare 70° Swingaway Footrests \$325
- ☐ I70-CFT Comes with composite foot plates N/C

- ☐ RLSL Residual Limb Support - Left (8.5" x 6")⁽⁴⁾
HCPCS code E1020 \$495
- ☐ RLSR Residual Limb Support - Right (8.5" x 6")⁽⁴⁾
HCPCS code E1020 \$495

1. Must choose footplate on page 15.

2. Heavy duty durable construction (4.5lbs ea.).

3. Not available on 14 and 15" seat widths. For 16" and 17" seat widths you must select flared receivers.

4. Swing away and removable, 50° of angle adjustment, 4.5" height adjustment, 3" horizontal and 2" of for/aft adjustment.

DUAL FRONT RIGGING ADAPTOR⁽¹⁾

- ☐ DFRA Dual Front Rigging Receivers \$450

1. Must purchase each set of front rigging separately. Not available on no charge front rigging options.

Ultra Low Maxx CG Power Positioning System

Storm Torque RX20

Dealer Name: _____





Dealer Account Number: _____

PO #: _____

Tag: _____

FOOTPLATES - SWING AWAY/ELEVATING LEGRESTS⁽¹⁾

HCPCS code K0040






Angle Adjustable Footplates

☐ AFSL	Left Small Footplate (4.25"W x 8"D)	\$165
☐ AFSR	Right Small Footplate (4.25"W x 8"D)	\$165
☐ AFML	Left Medium Footplate (5.25"W x 8"D)	\$165
☐ AFMR	Right Medium Footplate (5.25"W x 8"D)	\$165
☐ AFLL	Left Large Footplate (6.25"W x 8"D)	\$165
☐ AFLR	Right Large Footplate (6.25"W x 8"D)	\$165
☐ AFXL	Left Extra Large Footplate (9"W x 11"D)	\$165
☐ AFXR	Right Extra Large Footplate (9"W x 11"D)	\$165

Multi-Axis Adjustable Angle Footplates

☐ MFSL	Left Small Footplate (4.25"W x 8"D)	\$225
☐ MFSR	Right Small Footplate (4.25"W x 8"D)	\$225
☐ MFML	Left Medium Footplate (5.25"W x 8"D)	\$225
☐ MFMR	Right Medium Footplate (5.25"W x 8"D)	\$225
☐ MFL	Left Large Footplate (6.25"W x 8"D)	\$225
☐ MFLR	Right Large Footplate (6.25"W x 8"D)	\$225
☐ MFXL	Left Extra Large Footplate (9"W x 11"D)	\$225
☐ MFXR	Right Extra Large Footplate (9"W x 11"D)	\$225



Cut Out Angle Adjustable Footplates

☐ AFCSL	Left Small Footplate (7.25"W x 9"D)	\$165
☐ AFCSR	Right Small Footplate (7.25"W x 9"D)	\$165
☐ AFCML	Left Medium Footplate (7.5"W x 11"D)	\$165
☐ AFCMR	Right Medium Footplate (7.5"W x 11"D)	\$165
☐ AFCLL	Left Large Footplate (10"W x 14"D)	\$165
☐ AFCLR	Right Large Footplate (10"W x 14"D)	\$165

Cut Out Multi-Axis Adjustable Angle Footplates

☐ MFCSL	Left Small Footplate (7.25"W x 9"D)	\$225
☐ MFCSR	Right Small Footplate (7.25"W x 9"D)	\$225
☐ MFCML	Left Medium Footplate (7.5"W x 11"D)	\$225
☐ MFCMR	Right Medium Footplate (7.5"W x 11"D)	\$225
☐ MFCLL	Left Large Footplate (10"W x 14"D)	\$225
☐ MFCLR	Right Large Footplate (10"W x 14"D)	\$225

1. Certain footplate sizes may not fit all seat widths depending on legrests selected.

HEEL LOOPS, CUPS AND TOE STRAPS

☐ HLL	Left Heel Loop (each) HCPCS code E0951	\$20
☐ HLR	Right Heel Loop (each) HCPCS code E0951	\$20
☐ HCL	Left Cup Loop (each) HCPCS code E0951	\$20
☐ HCR	Right Cup Loop (each) HCPCS code E0951	\$20
☐ TSL	Left Toe Strap (each) HCPCS code E0952	\$20
☐ TSR	Right Toe Strap (each) HCPCS code E0952	\$20

CALF STRAPS AND CALF PANELS⁽¹⁾

☐ CSM	Body Point Airmesh Calf Strap for Chairs 15"-18" wide with 3.5" H Pad ⁽¹⁾	\$150
☐ CSL	Body Point Airmesh Calf Strap for Chairs 18"-22" wide with 3.5" H Pad ⁽¹⁾	\$150
☐ CPM	Body Point Airmesh Calf Panel for Chairs 15"-18" wide with 35.5" W x 9" H Pad ⁽¹⁾	\$220
☐ CPL	Body Point Airmesh Calf Panel for Chairs 18"-22" wide with 38.5" W x 9" H Pad ⁽¹⁾	\$235

1. Not available for Center Mounted Front Rigging.

STANDARD FRONT RECEIVERS

☐ PHDR	Pride HD Drop-in Receivers Only	N/C
☐ PSAR	Pride Style #8 Swing Away Receivers Only	N/C
☐ IAR	Invacare 70° Swing Away Receivers Only	N/C
☐ Q70R	Quickie 70° Swing Away Receivers Only	N/C
☐ Q70HDR	Quickie 70HD Swing Away Receivers Only	N/C

Ultra Low Maxx CG Power Positioning System

Storm Torque RX20

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

SINGLE-DUAL FUNCTION ELECTRONICS

- ☐ SFCB Single/Dual Function Electronics⁽¹⁾ STD
- Select Switch:
- ☐ SPBSO Single Push Button
- ☐ STO Single Toggle
- ☐ DPBO Dual Push Button
- ☐ OSO Omit Switch
- Optional:
- ☐ RDE Reduced Drive for Elevating Seat⁽²⁾ \$400
- ☐ SFCBA Single/Dual Function Drive Control⁽¹⁾
- HCPCS code E2310 \$1,795
- Optional:
- ☐ RDES Reduced Drive for Elevating Seat⁽²⁾ \$400
- ☐ DFCBD Dual Function Through Drive Control⁽¹⁾
- HCPCS code E2311 \$2,495
- Optional:
- ☐ RDED Reduced Drive for Elevating Seat⁽²⁾ \$400

1. These electronics include drive lock-out and tilt limit.

2. If reduced drive is not selected the system will be in drive lock-out when elevated.

MULTI- FUNCTION ELECTRONICS

- ☐ MFCBS Multi-Function through Switch⁽¹⁾ N/C
- Select Switch:
- ☐ QPBO Quad Push Button N/C
- ☐ FWT Four-Way Toggle Switch \$150
- ☐ EWR 8-Way Rocker \$205
- Optional:
- ☐ RDEM Reduced Drive for Elevating Seat⁽²⁾ \$400
- ☐ ACT4S Multi-Function Act4-LiNX Through Switch⁽¹⁾⁽⁴⁾ ... N/C
- Select Switch:
- ☐ AQPBO Quad Push Button N/C
- ☐ AFWTA Four-Way Toggle Switch \$150
- ☐ MFCBD Multi-Function M610 Through Drive Control⁽¹⁾⁽³⁾ HCPCS code E2311 \$2,995
- Up to 6 Seating Functions
- On Board Programming Seating Functions
- Specialty Products (i.e. PSB Power Sliding Back)
- Includes reduced drive for Elevating Seat
- ☐ ACT4 Multi-Function Act4-LiNX Through Drive Control⁽¹⁾⁽⁴⁾⁽⁵⁾ HCPCS code E2311 \$2,995
- Up to 5 Seating Functions
- Programming through LAK key/Cell Phone/Laptop
- Includes reduced drive for Elevating Seat

1. Requires Expandable Controller (EXPC) & Harness for Expandable System (PWH), select on page 2.

2. If reduced drive is not selected system will be in drive lock-out when elevated.

3. These electronics include seating attendant control, drive lock-out, tilt/recline limit and actuator speed control.

4. Only available on Multi-function Seating Systems.

5. The LAK key is not included, please select on page 3.

ELECTRONIC SWITCH MOUNTING

- ☐ ESML Left Side (Armpad)
- ☐ ESMR Right Side (Armpad)
- ☐ ESMJL Joystick Left (Mounts to Joystick Bracket)
- ☐ ESMJR Joystick Right (Mounts to Joystick Bracket)
- ☐ ESMD Do not mount

ELECTRONICS OPTIONS

Additional Switches

- ☐ SPBSA Single Push Button \$150
- ☐ STA Single Toggle Switch⁽⁴⁾ \$150
- ☐ DPBA Dual Push Button⁽⁴⁾ \$150
- ☐ QPBA Quad Push Button⁽⁴⁾ \$150
- ☐ FWT4 Four-Way Toggle Switch⁽⁴⁾ \$150
- ☐ EWRA 8-Way Rocker⁽¹⁾⁽³⁾ \$205
- ☐ TWRA 10-Way Rocker⁽¹⁾⁽³⁾ \$205
- ☐ EWPA 8-Way Push Button⁽¹⁾⁽⁴⁾ \$255
- ☐ 10WA 10-Way Switch⁽¹⁾⁽⁴⁾ \$255
- ☐ ASL320R ASL Pulse Mechanical Switch - Red \$125
- ☐ ASL320G ASL Pulse Mechanical Switch - Green \$125
- ☐ ASL320BK ASL Pulse Mechanical Switch - Black \$125
- ☐ SESG Stealth Egg Switch - Green \$135
- ☐ SESB Stealth Egg Switch - Black \$135
- ☐ FTSQ Feather Touch - Single Quarter Size \$295
- ☐ FTMD Feather Touch - Dual Dime Size⁽⁴⁾ \$295
- ☐ FTMQ Feather Touch - Quad Dime Size⁽⁴⁾ \$295

Specify Use for Additional Switch:

- ☐ STG Seating
- ☐ JM Joystick Mode⁽²⁾⁽³⁾
- ☐ JOO Joystick On/Off⁽²⁾⁽³⁾

Specify Mounting for Additional Switch:

- ☐ ESMLA Left Side (Armpad)
- ☐ ESMRA Right Side (Armpad)
- ☐ ESMJLA Joystick Left (Mounts to Joystick Bracket)
- ☐ ESMJRA Joystick Right (Mounts to Joystick Bracket)
- ☐ ESMDA Do not mount

Miscellaneous

- ☐ MJBS Motion Junction Box for Switch Control⁽¹⁾⁽⁵⁾⁽⁶⁾
- For Mode and/or Emergency Stop \$175
- ☐ S46L Single to Dual Switch Adapter \$265
- Allows for dual use of a single switch
- ☐ M616S M616 Converter⁽¹⁾ HCPCS code E2399 \$125
- Allows up to 4 mono phono switches into Multi-Function box
- ☐ LJS Lazarus Jump Starter HCPCS code A9270 \$150
- ☐ APPS Accessory Port Power Supply⁽¹⁾
- 3 amp capacity \$195
- ☐ TPSMH 12 Piece Switch Mounting Hardware Kit (additional kits) \$39

1. Not available with Single-Dual Function Electronics.

2. IN200 is required with REM210, REM211, REM216 joysticks.

3. Not available with REM110 joystick.

4. Not available for Joystick Mode or Joystick On/Off. Must Select Seating (STG)

5. Not available on Multi Function Act4-LiNX Electronics.

6. Cannot be used for power On/Off.

Ultra Low Maxx CG Power Positioning System

Storm Torque RX20

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

ARTICULATING TRAYS⁽¹⁾

- ☐ USAT Universal Style Articulating Tray \$1,350
14"W x Adjustable 6"-10"D
- Add:
- ☐ TMA Trilogy Mounting Adapter \$50
- ☐ AVT Articulating Vent (LTV) Tray \$1,350
- ☐ AVBH Articulating Vent Bag Holder⁽²⁾ \$1,350
- Add:
- ☐ ABMC Accessory Bag Mounting and Clips \$100

1. The user weight limit is 225 lbs. The maximum weight capacity of the vent tray is 35 lbs. Vents are not standard configuration on Elevating or Precline Systems or with Single Post Backs. Not available with TRRO and TRRO-E.
2. Motion Concepts does not supply the Vent Bag.

GIMBAL VENT TRAYS⁽¹⁾

- ☐ GUT Gimbal Universal Tray \$1,350
14"W x Adjustable 6"-10"D
- Add:
- ☐ GTMA Trilogy Mounting Adapter \$50
- ☐ GVT Gimbal Vent (LTV) Tray \$1,350
- ☐ GVBH Gimbal Vent Bag Holder⁽²⁾ \$1,350
- Add:
- ☐ GABMC Accessory Bag Mounting and Clips \$100

1. The Gimbal Vent Trays are not available with Tilt Only and/or Elevate Only Systems. The user weight limit is 200 lbs. Vents are not standard configuration on Elevating or Precline Systems or with Single Post Backs. The maximum weight capacity of the vent tray is 20 lbs. Not available with TRRO and TRRO-E.
2. Motion Concepts does not supply the Vent Bag.

O2 HOLDERS⁽¹⁾

- ☐ NAO2H Non-Articulating O2 Gas Holder \$395
- ☐ AOHG Articulating O2 Gas Holder Cylinder Style⁽²⁾⁽³⁾ \$1,250
- ☐ AOHL Articulating Liquid O2 Holder (6.5" x 4.5")⁽²⁾⁽³⁾ \$1,250
- ☐ NAIP L Non-Articulating IV Pole User's Left⁽⁴⁾ \$395
- ☐ NAIP R Non-Articulating IV Pole User's Right⁽⁴⁾ \$395
- ☐ AIVP Articulating IV Pole⁽²⁾⁽³⁾⁽⁴⁾ \$1,250

1. The user weight limit is 225 lbs. The maximum weight capacity of the O2 holder is 20 lbs. Not available with TRRO and TRRO-E.
2. Required for use on Recline Systems
3. Articulating O2 Holders are not standard configuration on Elevating or Precline Systems or with Single Post Backs.
4. The maximum weight capacity of the IV pole is 5 lbs.

VENTILATOR INFORMATION

Required when ordering a Vent Tray

- ☐ BIB Built in battery
- ☐ DTSVPS Dealer to supply vent power source

ACCESSORIES

- ☐ MUSB Motion Dual USB Charger N/C
- ☐ UCH Universal Cup Holder⁽⁴⁾ \$35
- ☐ PH X-Grip Phone Holder \$75
- ☐ AMC Accessory Mount with Clip (QTY ____) ⁽¹⁾ \$50
- ☐ CRVMR Rear View Mirror Right Mounted \$45
- ☐ CRVML Rear View Mirror Left Mounted \$45
- ☐ TAG571850 Wheelchair Cover⁽¹⁾ \$105
- ☐ HK Hydration Kit⁽²⁾ \$250
- ☐ MUB Medical Utility Bag w/ Mounting Hardware \$250
- ☐ THS Transfer Handles 4" - 6.5" \$350
- ☐ THL Transfer Handles 9" - 11.5" \$350
- ☐ FFLTL Flip and Fold Away 1/2 Lap Tray
Left Mounting(14"W x 10"D)⁽¹⁾⁽²⁾⁽³⁾ \$795
- ☐ FFLTR Flip and Fold Away 1/2 Lap Tray
Right Mounting(14"W x 10"D)⁽¹⁾⁽²⁾⁽³⁾ \$795
- ☐ CFH Flag Holder⁽¹⁾ \$65
- ☐ Left Mount
- ☐ Right Mount
- ☐ CSCC Cane and Crutch Holder \$395
- ☐ Left Mount
- ☐ Right Mount

1. Not available with TRRO, TRRO-E brackets
2. Not available on joystick side and only available with a flat arm pad or modular arm pad. Not available with Outback Cantilever Arms.
3. The maximum weight capacity of the tray is 10 lbs.
4. Not available on the joystick side.

BATTERIES

- ☐ M50-12 SLD M MK Group 50 Battery (ea) Quantity()
HCPCS code E2361 \$245
- ☐ M24 SLD G FT MK Group 24 Gel Battery (ea) Quantity()
HCPCS code E2363 \$365

INSTALLATION

- ☐ MBI Motion Concepts Installation⁽¹⁾ N/C

1. Batteries may be required if items on Accessories order form are selected at time of installation, please call to confirm.

NOTES

RETAIL \$ _____

SUBMIT

PRINT FORM

SAVE

CLEAR FORM