



Yes, you can.®

# Account Information

Request Type: Quote                      Order

Date: \_\_\_\_\_

Account #: \_\_\_\_\_

Company: \_\_\_\_\_

SHIP TO  
Name: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_

Country: \_\_\_\_\_

Purchase Order #: \_\_\_\_\_

CONTACT

Name: \_\_\_\_\_

Back Up Contact: \_\_\_\_\_

Phone: \_\_\_\_\_

Email: \_\_\_\_\_

Comments: \_\_\_\_\_

## Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

### IVC™ TRACER® SX5 WHEELCHAIRS

TRSX56FBP 16" x 16" Frame with Flip-Back Fixed Height  
Desk Length Arm ..... 1,086.96 K0003

TRSX58FBP 18" x 16" Frame with Flip-Back Fixed Height  
Desk Length Arm ..... 1,086.96 K0003

TRSX58FBFP 18" x 16" Frame with Flip-Back Fixed  
Height Full Length Arm ..... 1,086.96 K0003

TRSX50FBP 20" x 16" Frame with Flip-Back Fixed  
Height Desk Length Arm ..... \$1,086.96 K0006 + E2201

TRSX50FBFP 20" x 16" Frame with Flip-Back Fixed  
Height Full Length Arm ..... \$1,086.96 K0006 + E2201

TRSX52FBP 22" x 16" Frame with Flip-Back Fixed  
Height Desk Length Arm..... \$1,086.96 K0006 + E2201

### STOCK FRONT RIGGING OPTIONS

T93HAP Hemi Sw. Footrests Aluminum Footplates  
w/ Heel Loops ..... \$99.13 E0951

T93HCP Hemi Sw. Footrests Composite Footplates  
w/ Heel Loops ..... \$99.13 E0951

T94HAP Hemi Sw. Elev Legrests Aluminum Footplates &  
Padded Calf Pad..... \$191.86 K0195

T94HCP Hemi Sw. Elev Legrests Composite Footplates  
& Padded Calf Pad ..... \$191.86 K0195

AHL4A/U67 Hemi Smart Leg Articulating Legrests  
Aluminum Footplates .....\$341.09 K0053

T93HEP Hemi Sw. Footrests Composite Footplates  
(no heel loop)..... \$69.28

T94HEP Hemi Sw. Legrest Composite Footplates &  
Non-Padded Calfpad..... \$117.25 K0195



Yes, you can.®

# FRAME AND SEAT SELECTIONS

## CHAIR TYPE

TRSX5 Tracer SX5 Builder Wheelchair ..... \$1,446.00 K0003

## FRAME TYPE

The below selections are based on 24" wheels and 8" casters, see seat-to-floor matrix.

ADULT Adult 19 1/2" ..... No Charge

HEMI Hemi 17 1/2" ..... No Charge

## FRAME WIDTH AND DEPTH

WD46 14" Wide x 16" Deep ..... No Charge E1236  
Wheelchair base code changes to E1236 due to 14" seat width.

WD66 16" Wide x 16" Deep ..... No Charge

WD86 18" Wide x 16" Deep ..... No Charge

WD06 20" Wide x 16" Deep ..... No Charge K0006 + E2201  
20" and 22" widths have 300 lbs. weight limit requiring 1228 casters (upcharge necessary, see "Casters" section), and 24" rear wheels (no upcharge). Wheelchair base code changes to K0006 (instead of K0003).

WD26 22" Wide x 16" Deep ..... \$152.08 K0006 + E2201  
20" and 22" widths have 300 lbs. weight limit requiring 1228 casters (upcharge necessary, see "Casters" section), and 24" rear wheels (no upcharge). Wheelchair base code changes to K0006 (instead of K0003).

# ARM AND BACK OPTIONS

## ARM TYPE

The below arm types are flip-back, non-removable.

18FB Full Length Fixed Height Flip-Back 10" .. No Charge

19FB Full Length Adj. Height Flip-Back 10"-14" .....  
..... \$139.23 E0973

28FB Desk Length Fixed Height Flip-Back 10" .....  
..... No Charge

29FB Desk Length Adj. Height Flip-Back 10"-14" .....  
..... \$139.23 E0973

## BACK TYPE

BH16 Back Height 16" ..... Standard

BH18 Back Height 18" ..... No Charge

# COLOR OPTIONS

## UPHOLSTERY COLORS

U240 Black Nylon Upholstery ..... Standard

U550 Midnight Blue Vinyl Upholstery ..... No Charge

U67 Black Vinyl ..... No Charge

U86 Ruby Red Vinyl ..... No Charge



## POSITIONING STRAP OPTIONS

### SEAT POSITIONING STRAPS

1515 Push Button Style Positioning Strap .....\$53.55 E0978

1321BK Hook and Loop Seat Positioning Strap .....  
..... \$42.84 E0978

## CASTERS

### CASTER OPTIONS

AT902 6" Composite Urethane 3 Spoke ..... No Charge

AT903 8" x 1 1/4" Composite Urethane ..... No Charge

1255 8" x 1" Solid Rubber Dark Grey ..... Standard

1228 8" x 1 3/4" Semi Pneumatic Casters Lt Grey .....  
.....\$59.98

## WHEELS

### WHEEL OPTIONS

U1763BK 20" Composite Urethane Tires ..... \$121.02

U1764BK 22" Composite Urethane Tires ..... \$121.02

2244BK 24" Composite Pneumatic Tires .....  
..... \$148.51 E2211 + E2212

4244BK 24" Comp Pneumatic w/FF Inserts .....  
..... \$240.98 E2211 + E2212

U2222C 24" Composite Urethane Corded Tires .....  
..... Standard

U2244BK 24" Composite Urethane Treaded Tires .....  
.....\$128.52

### WHEEL OPTIONS (ONE ARM DRIVE)

The below selections are not available with the amputee attachment.

9456L One Arm Drive with U2222 LH ..... \$649.03 E0958

9456R One Arm Drive with U2222 RH..... \$649.03 E0958

### HANDRIM OPTIONS

COM Composite Handrims ..... Standard  
Available with 24" rear wheels only.

CHDRM Chrome Plated Handrims .....No Charge

AHDRM Aluminum Handrims ..... No Charge  
Not available with U2222C rear wheels.

### WHEEL LOCKS

9690B Pull-To-Lock Wheel Locks (Std on Hemi) .....  
..... No Charge

9691B Push-To-Lock Wheel Locks (Std on Adult) .....  
..... No Charge



Yes, you can.®

# ACCESSORIES

## BACK ACCESSORIES

STC Stencil Back Upholstery .....\$31.06

1160 Carrying Pocket ..... \$57.83

## SEAT ACCESSORIES

L8651 Omit Seat Upholstery .....Less \$10.00

## MISCELLANEOUS ACCESSORIES

1360 Adjustable Rear Anti-Tipper For 8" Casters .....  
.....\$87.82 E0971

1514 Overhead Anti-Theft Device ..... \$49.27  
Not available on 14" and 22" widths.

1360A Adjustable Rear Anti-Tipper For 6" Casters .....  
.....\$87.82 E0971

1497 Oxygen Holder Package Assembly ..... \$161.72 E2208  
Not available with One Arm Drive.

9357L Amputee Attachment ..... \$100.78 E0959

1496 IV/Oxygen Holder Package Assembly .....  
..... \$169.22 K0105 + E2208  
Not available with One Arm Drive.

1329BK Toggle Wheel Lock Extension (8") .....  
..... \$41.77 E0961

EXT2 Seat Extension Kit 2" ..... No Charge

1335 Non-Folding Device .....\$42.85  
Not available on 20" and 22" widths.  
Not available with One Arm Drive.

# FRONT RIGGING OPTIONS

## FRONT RIGGINGS

T93HE Hemi Sw. Footrests Composite Footplates  
(no heel loop).....\$69.28  
Not available with 14" wide.

AHL4A Hemi Smart Leg Articulating Legrests  
Aluminum Footplates ..... \$341.09 K0053  
Not available with 14" wide.  
All calf pads are black.

T93HA Hemi Sw. Footrests Aluminum Footplates  
w/ Heel Loops ..... \$99.13 E0951  
Not available with 14" wide.

AHL4C Hemi Smart Leg Articulating Legrests  
Composite Footplates ..... \$341.09 K0053  
Not available with 14" wide.  
All calf pads are black.

T93HC Hemi Sw. Footrests Composite Footplates  
w/ Heel Loops ..... \$99.13 E0951  
Not available with 14" wide.

AHS4A Hemi Smart Leg Articulating Legrestst-14"  
W Chairs ..... \$341.09 K0053  
All calf pads are black.

T94HE Hemi Sw. Legrest Composite Footplates &  
Non-Padded Calfpad ..... \$117.25 K0195

T93HS Hemi Footrests Aluminum Footplates w/ Heel ....  
Loops- 14" W Chairs ..... \$95.93 E0951

T94HA Hemi Sw. Elev Legrests Aluminum Footplates  
& Padded Calf Pad ..... \$191.86 K0195  
Not available with 14" wide.  
All calf pads are black.

T94HS Hemi Elevating Legrests Aluminum Footplates &  
Padded Calf Pad- 14" W Chairs .....\$186.53 K0195  
All calf pads are black.

T94HC Hemi Sw. Elev Legrests Composite Footplates &  
Padded Calf Pad ..... \$191.86 K0195  
Not available with 14" wide.  
All calf pads are black.



Yes, you can.®

### FOOTPLATES TYPE

1350 Extra Large Footplates ..... \$84.24 K0041  
Not available with 14" wide.

AT5543 Adjustable Angle Flip-Up Footplates .....  
.....\$71.42 K0040

### FOOTREST/LEGREST ACCESSORIES

1820 Non - Removable Device ..... \$12.79

1346 Heel Loops ..... \$41.57 E0951  
Already included on swingaway footrests T93HA, T93HC.

SPECIFICATIONS	
Seat Widths	14", 16", 18", 20" or 22"
Overall Width	Seat Width plus 8 1/2"
Seat Depth	16"
Overall Depth	3 1/8" (18 x 16 without Front Riggings)
Back Heights	16" or 18"
Seat-to-Floor	17 1/2" or 19 1/2"
Weight Limit	250 lbs. (300 lbs. on 20" or 22" widths)
Chair Weight	34 3/4" (18 x 16 without Front Riggings)
Shipping Weight	44 1/4 lbs.

20-22 Wide (300 lb. weight limitation) Wheelchairs only: 8 x 1 3/4 inch semi-pneumatic front casters MUST be installed. DO NOT replace with any other size/style front caster. Otherwise, injury or damage may occur.


### 6-Inch Front Casters

SEAT-TO-FLOOR HEIGHT (IN INCHES)	FRONT CASTER MOUNTING POSITION	REAR WHEEL SIZE	REAR WHEEL MOUNTING POSITION
15 1/2	Top	20 Inch	Top
16 1/2	Middle	22 Inch	Top
17 1/2	Bottom	24 Inch	Top
17 1/2	Bottom	20 Inch	Bottom

### 8-Inch Front Casters

SEAT-TO-FLOOR HEIGHT (IN INCHES)	FRONT CASTER MOUNTING POSITION	REAR WHEEL SIZE	REAR WHEEL MOUNTING POSITION
17 1/2	Top	20 Inch	Bottom
17 1/2	Top	24 Inch	Top
18 1/2	Middle	22 Inch	Bottom
19 1/2	Bottom	24 Inch	Bottom

Rear Wheel Mounting Positions - Adult/Hemi Frame Style



6 or 8-Inch Fork Mounting Positions

