



Yes, you can.®

Account Information

Request Type: Quote Order

Date: _____

Purchase Order #: _____

Account #: _____

CONTACT Name: _____

Company: _____

Back Up Contact: _____

SHIP TO Name: _____

Phone: _____

Address: _____

Email: _____

City: _____ State: _____ Zip: _____

Comments: _____

Country: _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

PRE-CONFIGURED STOCK WHEELCHAIRS

STOCK CHAIRS

TRSX5RC6P 16" x 16" Recliner Frame with Fixed Height Conventional Desk Arm \$1,270.29 K0003 + E1226

TRSX5RC8P 18" x 16" Recliner Frame with Fixed Height Conventional Desk Arm \$1,270.29 K0003 + E1226

STOCK FRONT RIGGINGS

T93HAP Hemi Sw. Footrests Aluminum Footplates w/ Heel Loops \$99.13 E0951

AHL4A/U67 Hemi Smart Leg Articulating Legrests Aluminum Footplates \$341.09 K0053

T93HCP Hemi Sw. Footrests Composite Footplates w/ Heel Loops \$99.13 E0951

T93HEP Hemi Sw. Footrests Composite Footplates (no heel loop) \$69.29

T94HAP Hemi Sw. Elev Legrests Aluminum Footplates & Padded Calf Pad \$191.86 K0195

T94HEP Hemi Sw. Legrest Composite Footplates & Non-Padded Calfpad \$117.25 K0195

T94HCP Hemi Sw. Elev Legrests Composite Footplates & Padded Calf Pad \$191.86 K0195

FRAMES

CHAIR TYPE

TRSX5RC Tracer SX5 Recliner Builder Wheelchair \$1,410.00

Note: See Frame Width & Depth for Base HCPC Codes



Yes, you can.®

FRAME WIDTH AND DEPTH

WD46 14" Wide x 16" Deep	\$96.39 E1236 + E1226
WD48 14" Wide x 18" Deep	\$126.38 E1236
WD66 16" Wide x 16" Deep	No Charge
WD68 16" Wide x 18" Deep	\$126.38
WD86 18" Wide x 16" Deep	No Charge
WD88 18" Wide x 18" Deep	\$126.38
WD06 20" Wide x 16" Deep	No Charge E2201
20" and 22" widths have 300 lbs. weight limit requiring 1228 casters (upcharge necessary, see "Casters" section), and 24" rear wheels (no upcharge).	

WD08 20" Wide x 18" Deep \$272.03 E2201
20" and 22" widths have 300 lbs. weight limit requiring 1228 casters (upcharge necessary, see "Casters" section), and 24" rear wheels (no upcharge).

WD26 22" Wide x 16" Deep No Charge E2201
20" and 22" widths have 300 lbs. weight limit requiring 1228 casters (upcharge necessary, see "Casters" section), and 24" rear wheels (no upcharge).

WD28 22" Wide x 18" Deep \$272.03 E2201
20" and 22" widths have 300 lbs. weight limit requiring 1228 casters (upcharge necessary, see "Casters" section), and 24" rear wheels (no upcharge).

ARMS

ARM TYPE

All arm types below are removable, conventional and nonflip-back.

18 Full Length Fixed Height Conv. 10"	No Charge
19 Full Length Adj. Height Conv. 10" - 14"	
.....	\$139.23 E0973

28 Desk Length Fixed Height Conv. 10" No Charge

29 Desk Length Adj. Height Conv. 10" - 14" \$139.23 E0973

FRAME TYPE

R Recliner Frame	\$439.11 E1226
90 degree to 180 degree of recline. 24" back height with 10" headrest extension.	

COLOR OPTIONS

UPHOLSTERY COLORS

U240 Black Nylon Upholstery	No Charge
U550 Midnight Blue Vinyl Upholstery	Standard
U67 Black Vinyl	No Charge
U86 Ruby Red Vinyl	No Charge



Yes, you can.®

CASTERS

CASTER OPTIONS

AT903 8" x 1.25" Composite Urethane Standard
1228 8" x 1.75" Semi Pneumatic Casters Lt Grey ...\$59.98

4252 8" x 1.25" Pneumatic Casters w/FF Inserts
.....\$191.71 E2217

WHEELS

WHEEL OPTIONS

U1763BK 20" Composite Urethane Tires \$121.02
U1764BK 22" Composite Urethane Tires \$121.02
2244BK 24" Composite Pneumatic Tires
..... \$149.94 E2211+E2213
4244BK 24" Composite Pneumatic with FF Inserts
..... \$240.98 E2211+E2213

U2222C 24" Composite Urethane Corded Tires
..... Standard E2211 + E2212
U2244BK 24" Composite Urethane Treaded Tires
.....\$128.52

WHEEL LOCKS

9690B Pull To Lock Wheel Locks (Std on Hemi)
..... No Charge

9691B Push To Lock Wheel Locks (Std on Adult).....
..... No Charge

HANDRIM OPTIONS

COM Composite Handrims Standard
Available with 24" rear wheels only.

AHDRM Aluminum Handrims No Charge
Not available with U2222C rear wheels.

POSITIONING STRAP OPTIONS

SEAT POSITIONING STRAPS

1515 Push Button Style Positioning Strap\$53.35 E0978

1321BK Hook and Loop Seat Positioning Strap
..... \$42.84 E0978

ACCESSORIES

BACK ACCESSORIES

L8601 Omit Back UpholsteryLess \$0.00

STC Stencil Back Upholstery \$31.06
3/4" Block Letters.

SEAT ACCESSORIES

L8651 Omit Back UpholsteryLess \$0.00



Yes, you can.®

MISCELLANEOUS ACCESSORIES

1322SXR IV Rod and Bracket, Telescopic, Recliner \$117.81

MUPH Matching Upholstery No Charge
Applies to arm pads and calf pads.

1329BK Toggle Wheel Lock Extension (8") ... \$41.77 E0961

FRONT RIGGING OPTIONS

FRONT RIGGINGS

T93HE Hemi Sw. Footrests Composite Footplates
(no heel loop) \$69.28

T93HA Hemi Sw. Footrests Aluminum Footplates w/
Heel Loops \$99.13 E0951
Not available with 14" wide.

T93HC Hemi Sw. Footrests Composite Footplates w/
Heel Loops \$99.13 E0951
Not available with 14" wide.

T94HE Hemi Sw. Legrest Composite Footplates &
Non-Padded Calfpad \$117.25 K0195

T94HA Hemi Sw. Elev Legrests Aluminum Footplates &
Padded Calf Pad \$191.86 K0195
Not available with 14" wide.
All calf pads are black.

T94HC Hemi Sw. Elev Legrests Composite Footplates
& Padded Calf Pad \$191.86 K0195
Not available with 14" wide.
All calf pads are black.

AHL4A Hemi Smart Leg Articulating Legrests
Aluminum Footplates \$341.09 K0053
Not available with 14" wide.
All calf pads are black.

AHL4C Hemi Smart Leg Articulating Legrests
Composite Footplates \$341.09 K0053
Not available with 14" wide.
All calf pads are black.

AHS4A Hemi Smart Leg Articulating Legrestst-14"
W Chairs \$341.09 K0053
All calf pads are black.

T93HS Hemi Sw. Footrests Aluminum Footplates
w/ Heel Loops- 14" W Chairs \$95.93 E0951

T94HS Hemi Sw. Elevating Legrests Aluminum
Footplates & Padded Calf Pad- 14" W Chairs \$186.53 K0195
All calf pads are black.

FOOTPLATE TYPE

1350 Extra Large Footplates \$84.21 K0041
Not available with 14" wide.

AT5543 Adjustable Angle Flip-Up Footplates \$71.42 K0040

FOOTREST/LEGREST ACCESSORIES

1820 Non - Removable Device \$12.79

1346 Heel Loops \$41.57 E0951
Standard on Swingaway Footrests T93HA and T93HC.



Yes, you can.®

Seat-to-Floor 8" Casters

SEAT-TO-FLOOR HEIGHT	FRONT CASTER MOUNTING POSITION	REAR WHEEL SIZE	REAR WHEEL MOUNTING POSITION
17 1/2	Top	20"	Recliner Axle Mounting Bracket
18 1/2	Middle	22"	Recliner Axle Mounting Bracket
19 1/2	Bottom	24"	Recliner Axle Mounting Bracket

Standard on SX5 Recliner

- Rear anti-tippers
- Patient only operated push-to-lock wheel lock

Weight Limitation 250 lbs., 300 lbs. 20" & 22" width

Invacare recommends all wheelchairs be fitted to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.

